009 1561J776R

ADHEISIVE DENTAL RESTORATIVE

CLEARFIL FII NEW BOND

2°C /
36°F8°C /
46°FCE
0197

ENGLISH INSTRUCTIONS FOR USE

GENERAL INFORMATION

CLEARFIL FII is a chemically hardening hybrid composite, specifically developed for the construction of anterior composite restorations. If you combine this composite with CLEARFIL NEW BOND, you will ensure a good adhesion to the enamel and dentin, which will in turn provide good marginal sealing. It is classified as a Type 2 and Class 1 material by ISO 4049. The general clinical benefit of this product is to restore tooth function for the following INDICATIONS FOR USE.

I. INDICATIONS FOR USE

- Class III, IV, V restorations of anterior teeth.

II. Contraindication

- Patients with a history of hypersensitivity to methacrylate monomers

III. Precautions

[Safety precautions]

- This product contains substances that may cause allergic reactions. Avoid use of the product in patients with known allergies to methacrylate monomers or any other components.
- Wear protective gloves or take other appropriate measures to prevent the occurrence of hypersensitivity due to contact with this product.
- If any hypersensitivity, such as a rash or dermatitis occurs, discontinue the use of the product and consult a doctor.
- Use caution to prevent the product from coming in contact with soft tissue or the skin, or getting into eye.
- If the product comes in contact with the human body, take the following measures:
 - If the product gets in the eye, wash the eye with copious amounts of running water and consult a doctor.
 - If the product comes in contact with soft tissue or the skin, wipe it away with a cotton pledge or gauze moistened with alcohol and wash with copious amounts of water.
- Dispose of this product as a medical waste to prevent infection.

[Storage precautions]

- Do not use after expiration date. Note expiration date on outside of package.
- The product must be refrigerated when not in use: bring it to room temperature before using.
- Keep away from extreme heat or direct sunlight.

IV. Components

1. Components

Please see the outside of the package for contents and quantity.

CLEARFIL FII NEW BOND

- Catalyst paste
- Universal paste

2. Ingredients

- Catalyst Paste**
Silanated silica filler (60-90%), Bisphenol A diglycidylmethacrylate (Bis-GMA) (5-25%), Triethylene glycol dimethacrylate (TEGDMA) (< 7%), Colloidal silica (< 3%), Benzoyl peroxide (< 1%)
- Universal Paste**
Silanated silica filler (60-90%), Bisphenol A diglycidylmethacrylate (Bis-GMA) (10-30%), Colloidal silica (< 3%), Neopentylglycol dimethacrylate (1-5%), 1,2-Bis(3-methacryloxy-2-hydroxypropoxy)ethane (1-5%)

Units in parentheses are mass %.

The total amount of inorganic filler is approx. 59 vol%. The particle size of inorganic fillers ranges from 0.01 µm to 54 µm.

V. Clinical procedure

1. Isolation

Rubber dam is recommended for isolation.

2. Cavity preparation

Determine the presence of caries and remove any infected dentin if present. CARIES DETECTOR is recommended to serve as a guideline for safe caries removal.

3. Pulpal protection

Generally there is no need for cement lining or basing. However, any actual or near-pulpal exposure area should be covered with a hard setting calcium hydroxide material.

4. Etching

It is recommended to etch the cavity for 30 to 60 seconds with the aid of K-ETCHANT GEL. Then thoroughly rinse the cavity with a spray for at least 20 second after which it must be carefully blow-dried. The enamel should now have a mat, white surface.

From this moment onwards it is important for the cavity to remain absolutely dry and clean, in order to make an optimum adhesion process possible.

Should contamination of the etched surface occur unexpectedly (for instance caused by blood or saliva), it is necessary for the entire etching procedure to be repeated.

5. Applying CLEARFIL NEW BOND

Mix two components: Universal Liquid and Catalyst Liquid in the proportion of 1:1. Apply a thin coating of the mixture with a sponge on the entire wall of the cavity, the enamel and dentine. Evaporate the solvent gently and carefully with a flow of air.

Place the mixed composite in the cavity within one minute.

6. Mixing the pastes

Take the relevant quantity of Universal paste from the jar with a plastic spatula and put in the mixing dish. Take a similar quantity of Catalyst paste with the other end of the spatula. Mix both pastes vigorously and thoroughly for 15 seconds in the mixing dish.

7. Applying the composite into the cavity

You can apply the composite into the cavity by means of hand instruments as well as the syringe technique. It is recommended to use synthetic matrix strips for this.

Mixing time (at 25°C) 15 secs.
↓

Working time (at 25°C) 105 secs.
↓

Setting time (at 25°C) 30 secs.

8. Finishing and polishing

Contour the restoration and adjust occlusion and articulation using a fine diamond point. Polish with silicon rubber points or polishing discs.

WARNING

Kuraray Noritake Dental Inc. will replace any product that is proven to be defective. Kuraray Noritake Dental Inc. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

[NOTE]

If a serious incident attributable to this product occurs, report it to the manufacturer's authorized representative shown below and the regulatory authorities of the country in which the user/patient resides.

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"CLEARFIL" is a registered trademark or trademark of KURARAY CO., LTD.
La Kuraray Noritake Dental Inc. sostituisce tutti i prodotti che dovessero risultare difettosi.
"CLEARFIL" è un marchio registrato o un marchio di fabbrica di KURARAY CO., LTD.
Kuraray Europe GmbH (EU Importer)