CLEARFIL™ LINER BOND 2V

EN	ENGLISH	English	RESIN-BASED DENTAL ADHESIVE SYSTEM
FR	FRANÇAIS	French	SYSTÈME DE COLLAGE DENTAIRE À BASE DE RÉSINE
ES	ESPAÑOL	Spanish	SISTEMA ADHESIVO DENTAL CON BASE DE RESINA
IT	ITALIANO	Italian	SISTEMA ADESIVO DENTALE A BASE DI RESINA
NL	NEDERLANDS	Dutch	DENTAAL ADHESIEF SYSTEEM OP BASIS VAN KUNSTHARS
DE	DEUTSCH	German	DENTALES ADHÄSIVSYSTEM AUF KUNSTHARZBASIS
SV	SVENSKA	Swedish	RESINBASERAT DENTALT ADHESIVSYSTEM
NO	NORSK	Norwegian	DENTALT ADHESIVSYSTEM BASERT PÅ RESIN
FI	SUOMI	Finnish	YHDISTELMÄMUOVIPOHJAINEN HAMPAAN KIINNITYSAINE
DA	DANSK	Danish	RESINBASERET DENTALT ADHÆSIVSYSTEM
PT	PORTUGUÊS	Portuguese	SISTEMA ADESIVO DENTÁRIO À BASE DE RESINA
EL	ΕΛΛΗΝΙΚΑ	Greek	ΟΔΟΝΤΙΚΟ ΣΥΓΚΟΛΛΗΤΙΚΟ ΣΥΣΤΗΜΑ ΜΕ ΒΑΣΗ ΤΗ ΡΗΤΙΝΗ
TR	TÜRKÇE	Turkish	REZİN BAZLI DENTAL ADEZİV SİSTEM

Kuraray Noritake Dental Inc.

1621 Sakazu, Kurashiki, Okayama 710-0801, Japan

EMERGO EUROPE

Westervoortsedijk 60,
6827 AT Arnhem, The Netherlands


Kuraray Europe GmbH (EU Importer)

Philipp-Reis-Str. 4,
65795 Hattersheim am Main, Germany
Phone: +49 (0)69 305 35835 Fax: +49 (0)69 305 98 35835
URL: <https://www.kuraraynoritake.eu>

1561R519R

CLEARFIL™ LINER BOND 2V

I. INTRODUCTION

CLEARFIL LINER BOND 2V is a light- and dual- cure bonding system, and consists of a self-etching primer and a bonding agent. The primer offers simultaneous treatment of both dentine and enamel. The general clinical benefit of this product is to restore tooth function for the following INDICATIONS FOR USE.

II. INDICATIONS FOR USE

CLEARFIL LINER BOND 2V is indicated for the following cases:

- Direct filling restorations using light-cure or chemical-cure composite resin
- Bonded amalgam restorations
- Treatment of hypersensitive and/or exposed root surfaces
- Cavity sealing as a pretreatment for indirect restorations
- Intraoral repairs of fractured porcelain facing crowns using light-cure composite
- Cementing laminated veneer, inlays and onlays made of porcelain (or composite) using composite resin cement

III. CONTRAINDICATIONS

Patients with a history of hypersensitivity to methacrylate monomers

IV. POSSIBLE SIDE EFFECT

The mucous membrane may turn whitish when contacted with primer due to the coagulation of protein. This is a temporary phenomenon that will disappear in a few days.

V. INCOMPATIBILITY

- 1) Do not use eugenol-containing materials for pulp protection or temporary sealing, since the eugenol could retard the bonding system curing process.
- 2) Do not use hemostatic agents, especially those containing ferric compounds, since these materials may impair adhesion and may cause discoloration at the tooth margin or surrounding gingiva due to ferric ions that may remain.

VI. PRECAUTIONS

1. Safety precautions

1. Avoid use of the product for patients with a history of hypersensitivity to methacrylate monomers.
2. If any hypersensitivity such as dermatitis occurs, discontinue use of the product and consult a physician.
3. Use caution to prevent the product from coming in contact with the skin or getting into the eye. Before using the product, cover the patient's eyes with a towel to protect the patient eyes from splashing material.
4. If the product comes in contact with the human body, take the following measures:
 - <If the product gets in the eyes>
Immediately wash the eye with copious amounts of water and consult a physician.
 - <If the product comes in contact with the skin>
Immediately wipe it off with a cotton pledget moistened with alcohol or gauze and wash with copious amounts of water.
5. Use caution to prevent the patient from accidentally swallowing the product.
6. Avoid looking directly at the curing light when curing the product.
7. Avoid using the same disposable brush tip for different patients to prevent cross infection. Discard the tip after use and sterilize the brush tip handle.
8. Dispose of this product as a medical waste to prevent infection.

2. Handling and manipulation precautions

1. The emitting tip of the visible light-cure activator should be held as near and vertical to the resin surface as possible. If a large resin surface is to be light cured, it is advisable to divide the area into several sections and light-cure each section separately.
2. Low intensity of light causes poor adhesion. Check the lamp for service life and the visible light guide tip for contamination. It is advisable to check the light intensity of light-cure activator using an appropriate lightchecker at periodic intervals.
3. Use caution not to drop the system tray. Do not disassemble it.
4. The use of this product is restricted to dental professionals.

3. Storage precautions

1. The product must be used by the expiration date indicated on the package.
2. The product must be stored in a refrigerator (2 - 8°C / 36 - 46°F) when not in use ; bring it to room temperature before using.
3. Keep away from extreme heat or direct sunlight.
4. Replace the bottle cap as soon as possible after resin has been dispensed from the bottle.
5. The product must be carefully stored so that only a licensed dental professional can access it.

VII. COMPONENTS

Please see the outside of the package for contents and quantity.

1) PRIMER Liquid A

Principal ingredients:

- 2-Hydroxyethyl methacrylate (HEMA) (10-25%)
- 10-Methacryloyloxydecyl dihydrogen phosphate (MDP) (20-40%)
- Hydrophilic aliphatic dimethacrylate (1-15%)
- dl-Camphorquinone (< 1%)
- Water (30-50%)
- Accelerators (< 15%)
- Dyes (< 0.1%)

2) PRIMER Liquid B

Principal ingredients:

- 2-Hydroxyethyl methacrylate (HEMA) (20-40%)
- Hydrophilic aliphatic dimethacrylate (5-15%)
- Water (45-65%)
- Catalysts (< 9%)

3) BOND Liquid A

Principal ingredients:

- Bisphenol A diglycidylmethacrylate (Bis-GMA) (25-45%)
- 2-Hydroxyethyl methacrylate (HEMA) (20-40%)
- 10-Methacryloyloxydecyl dihydrogen phosphate (MDP) (1-10%)
- Hydrophobic aliphatic dimethacrylate (10-30%)
- Colloidal silica (5-15%)
- dl-Camphorquinone (< 2%)
- Accelerators (< 5%)

4) BOND Liquid B

Principal ingredients:

- Bisphenol A diglycidylmethacrylate (Bis-GMA) (30-50%)
- 2-Hydroxyethyl methacrylate (HEMA) (20-40%)
- Benzoyl peroxide (< 2%)
- Hydrophobic aliphatic dimethacrylate (10-30%)
- Colloidal silica (3-10%)

Units in parentheses are mass %.

5) Accessories

- Sponge pledget
- Mixing dish*
- Light blocking plate*
- System tray*
- Brush tip handle*
- Disposable brush tips
- *Consumables

VIII. CLINICAL PROCEDURES

A. Direct filling restorations, cavity sealing and restoration of root surfaces using light-cure composite

A-1. Cleaning of tooth structure

An insufficiently cleaned cavity may impair the adhesive performance. Make sure the cavity is sufficiently cleaned.

A-2. Moisture control

In order to produce optimal result, avoid contamination of the treatment area (saliva, blood). A rubber dam is recommended to keep the tooth dry.

A-3. Cavity preparations

Remove any infected dentine using CRIES DETECTOR and prepare the cavity in the usual manner.

A-4. Pulp protection

Any actual or near pulp exposure area should be covered with a hard setting calcium hydroxide material. There is no need for cement lining or basing. Do not use eugenol materials for pulp protection.

A-5. Acid etching on uncut enamel

If there is a possibility of resin spreading over uncut enamel, apply K-ETCHANT GEL to the enamel, let it remain for 10 seconds, wash with water, and then dry.

[CAUTION]

The use of PRIMER alone does not sufficiently condition uncut enamel. Overfilling of resin onto unetched uncut enamel could cause marginal discoloration.

A-6. Tooth surface treatment

1. Dispense equal amounts of PRIMER Liquids A and B into a well of the mixing dish and mix immediately before application.
2. Apply the mixture to the entire cavity wall with a sponge or a disposable brush tip. Leave in place for 30 seconds. Use caution not to allow saliva or exudate to contact the treated surfaces for at least 30 seconds.
3. After conditioning the tooth surface for 30 seconds, evaporate the volatile ingredients with a mild oil-free air stream.

[CAUTION]

Avoid pooling of the primer. Also do not wash. Observe the specified mixing amount, drying method, and treatment time to ensure optimum adhesion. Avoid touching treated surface. If the treated surface is contaminated, wash it with water, dry, or clean with alcohol, and treat with PRIMER again.

[NOTE]

PRIMER must be mixed immediately before use. After mixing, the mixture must be stored in dark area and used within 3 minutes.

A-7. Bonding

1. Dispense the necessary amount of BOND Liquid A into a well of the mixing dish.

[NOTE]

Use it as soon as possible after dispensing; it will harden if exposed for any extended period to the operating light.

2. Apply BOND Liquid A to the entire surface of the cavity with a sponge or a disposable brush tip.
3. After application, make the bond film as uniform as possible using a gentle oil-free air stream.

[CAUTION]

A strong air blast will scatter the bonding agent, resulting in poor adhesion.

4. Light-cure the BOND for 20 seconds with a dental curing light (irradiation wave length : 400-515nm, light intensity : >300mW/cm²).

[NOTE]

BOND Liquid A is light-cure type. Use it as soon as possible after dispensing; it will harden if exposed for any extended period to the operating light.

A-8a. Light-cure composite filling

Apply composite (e.g. CLEARFIL MAJESTY Esthetic, CLEARFIL MAJESTY Posterior) into the cavity, light-cure, finish and polish according to the Manufacturer's instructions for use.

A-8b. Cavity sealing and restoration of root surface

Apply a thin coat of flowable composite resin (e.g. CLEARFIL MAJESTY Flow) to the tooth, and light-cure the composite according to its' Instructions for Use.

Remove unpolymerized resin with a cotton pledget moistened with alcohol.

B. Direct filling restorations using chemical-cure composite or bonded amalgam restorations

B-1. Cleaning of tooth structure

B-2. Moisture control

B-3. Cavity preparation

B-4. Pulp protection

B-5. Acid etching on uncut enamel

B-6. Tooth surface treatment

Use the same procedures described in A-1 - A-6 for the above steps.

B-7. Bonding

1. Dispense equal amounts of BOND Liquids A and B into a well of the mixing dish and mix immediately before application.
2. Apply the mixture to the entire surface of the cavity with a sponge or a disposable brush tip.
3. After application, make the bond film as uniform as possible using a gentle oil-free air stream.

[CAUTION]

A strong air blast will scatter the bonding agent, resulting in poor adhesion.

[CAUTION]

Observe the specified mixing amount, applying method, and curing time; otherwise poor adhesion will result.

[NOTE]

BOND Liquids A and B should be mixed immediately before use. The mixture must be stored in a dark area and used within 5 minutes after mixing.

B-8a. Chemical-cure composite filling

Mix the composite pastes, apply the mixture into the cavity, finish and polish according to the manufacturer's instructions for use.

B-8b. Amalgam filling

Mix amalgam powder and mercury, condense the mixture into the cavity, finish and polish according to the manufacturer's instructions for use.

[NOTE]

Unpolymerized layer of BOND is necessary in order to ensure bonded amalgam restoration.

C. Intraoral repairs of fractured porcelain facing crowns using light-cure composite

C-1. Preparation of fractured surfaces

- 1) Porcelain surface
Using a diamond point, remove a layer of the porcelain, and place a bevel at the marginal area.
- 2) Metal surface
Roughen the metal surface with a diamond point.

C-2. Acid etching of porcelain surface

Apply K-ETCHANT GEL to the porcelain surface and leave it in place for 5 seconds before washing and drying.

If the adherent surface extends to uncut enamel, apply K-ETCHANT GEL and leave it in place for 10 seconds before washing and drying.

C-3. Silane treatment

1. Dispense one drop each of PRIMER Liquids A, B and CLEARFIL PORCELAIN BOND ACTIVATOR into a well of the mixing dish and mix them.

[NOTE]

PRIMER and CLEARFIL PORCELAIN BOND ACTIVATOR should be mixed immediately before application.

2. Apply the mixture to the porcelain and metal surfaces with a sponge or a disposable brush tip.
3. Leave it in place for 5 seconds. Use caution not to allow saliva or exudate to contact the treated surfaces.
4. If the adherent surface extends to the tooth, leave the mixture in place for 30 seconds.
5. After applying the primer, evaporate the volatile ingredients with a mild oil-free air stream. The treated metal surface must be sufficiently dried; if not, adhesion might be seriously impaired.

[CAUTION]

Do not wash treated surface. If the treated surface is contaminated with saliva, wash it with water, dry, or clean with alcohol, and treat with the mixture again.

C-4. Bonding

1. Dispense the necessary amount of BOND Liquid A into a well of the mixing dish.

[NOTE]

Use it as soon as possible after dispensing; it will harden if exposed for any extended period to the operating light.

2. Apply the liquid to the entire cavity wall with a sponge or a disposable brush tip.
3. Use a light air stream to make the bond film as uniform as possible.
4. Light-cure the BOND for 20 seconds with a dental curing light (irradiation wave length : 400-515nm, light intensity : >300mW/cm²).

C-5. Light-cure composite filling

Apply composite (e.g. CLEARFIL MAJESTY Esthetic, CLEARFIL MAJESTY Posterior) into the cavity, light-cure, finish and polish according to the Manufacturer's instructions for use.

[NOTE]

Use opaque resin for metal surface.

D. Cementing laminated veneer, inlays and onlays made of porcelain (or composite) using composite resin cement

D-1. Treatment of prosthetic restorations

Perform the following steps according to the procedures C-2 and C-3.

1. Acid etching of porcelain surface
2. Silane treatment

D-2. Tooth surface treatment and cementing

Treat the tooth surface and cement the prosthetic restorations according to the manufacture's instructions of the composite resin cement.

[WARRANTY]

Kuraray Noritake Dental Inc. will replace any product that is proved to be defective. Kuraray Noritake Dental Inc. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

[NOTE]


If a serious incident attributable to this product occurs, report it to the manufacturer's authorized representative shown below and the regulatory authorities of the country in which the user/patient resides.

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"CLEARFIL" and "CLEARFIL MAJESTY" are registered trademarks or trademarks of KURARAY CO., LTD.
"CLEARFIL" ve "CLEARFIL MAJESTY", KURARAY CO., LTD.'nin tescilli ticari markaları veya ticari markalarıdır.