







# ENGLISH



POLYMER-BASED CROWN AND BRIDGE MATERIAL

***ESTENIR™ C&B*****C O N T E N T S**

<b>I</b>	<b>Introduction</b>	3
<b>II</b>	<b>Contraindication</b>	3
<b>III</b>	<b>Precautions</b>	3
<b>IV</b>	<b>Indications for use</b>	6
<b>V</b>	<b>Kit components</b>	7
<b>VI</b>	<b>Basic laboratory procedures</b>	11
1.	Flowchart	11
2.	Curing unit, polymerization time and polymerization depth	12
3.	Color matching	13
4.	Bridge frameworks	20
5.	Fabrication of jacket crowns	21
6.	Fabrication of inlays and onlays	25
7.	Fabrication of bridges with frameworks	27
8.	Fabrication of facing cast crowns (superstructure of implants)	31
9.	Additional fabrication procedures	34
<b>VII</b>	<b>Clinical procedures</b>	35
1.	Preparing abutments and cavities	35
2.	Loading crown restorations	37
3.	Repairing crown restorations	38

## I Introduction

ESTENIA™ C&B is a polymer-based crown and bridge material for creation of facing cast crowns, facing cast bridges, jacket crowns, inlays, onlays and bridges with frameworks.

## II Contraindication

This product should not be used for patients with a history of hypersensitivity, such as rash or dermatitis, when exposed to methacrylate monomer.

## III Precautions

\*Please review this Instructions for Use carefully before using the product. After reviewing the Instructions for Use, keep them in a convenient place for quick reference if necessary. If this Instructions for use are lost, contact the distributor or supplier shown on the package.

\*This product should only be used by dentists or dental technicians.

\*To prevent the occurrence of trouble with operation or storing, the product should be used according to the instructions specified in this Instructions for Use. Do not use this product for contraindications given in this Instructions for Use.

### Safety precautions

①The use of this product should be determined by a dentist individually for each patient.

②This product should NOT be used when the following conditions are present:

- Malocclusion, clenching or bruxism
- Closed bite or cases where there is no occlusal support for biting except the tooth being treated (If used in these situations, the occlusal contact area should be covered with metal.)
- The need to cover the occlusal surfaces of molars (The product can be used if there is some occlusal support other than the tooth under treatment.)
- Clasped teeth for partial dentures (The tooth area which comes in contact with the clasp should be covered with metal.)
- Patients requiring bridges replacing more than two missing teeth, or with an edentulous span of more than 15 mm.

③If any hypersensitivity, such as rash or dermatitis occurs, discontinue the use of the product and get medical attention.

④Avoid any direct contact of material components with the human body to prevent hypersensitivity.


Wear gloves or take other appropriate measures when using the product. In particular, avoid touching unhardened resin.


⑤Use caution to prevent the product from coming in contact with the oral soft tissues or the skin or getting into the eyes. If the product comes in contact with oral tissue or the skin, wipe it off with a cotton pledge or a piece of gauze moistened with alcohol and then wash immediately with copious amounts of water. If the products get into an eye, immediately wash the eye with copious amounts of water and consult an ophthalmologist.

⑥Use caution to prevent anyone from accidentally swallowing the product.

⑦Use safety equipment, such as a local vacuum unit, protective shields, and an appropriately approved dustproof mask, to prevent the operator from inhaling dust when shaping, finishing or polishing the product.

⑧When any of the related materials covered in this manual are used with this product, read the Instructions for Use for that material carefully.

## Handling and manipulation precautions

### ① Avoid cracked or broken crowns or inlays

- Assure that the crown or inlay is sufficiently thick, using the methods for preparing abutments and cavities described in this manual.
- Assure correct curing by using the light curing units and heat curing units as specified in this manual. Observe the specified curing conditions.
- Review the part of additional fabrication procedures in this manual for additional fabrication.
- Do not mix different types of Body resin together, or mix Body resin with other materials, in order to prevent impairing the physical properties of the materials.

### ② Avoid broken or cracked bridges with frameworks

- Assure that the bridge is sufficiently thick, using the methods for preparing abutments and cavities described in this manual. Take particular care not to use this product in cases where sufficient thickness cannot obtain.
- Assure correct curing by using the light curing units and heat curing units as specified in this manual. Observe the specified curing conditions.
- To prevent the impairment of product physical properties, do not use framework or crown materials supplied by other companies.

### ③ Avoid the separation of resin from metal frames

- Use spherical retention beads with a particle size of 100 to 200 µm on the facing surface of the metal frame wax mold.
- Sandblast the facing surface of the metal frame with alumina particles, 50 µm diameter.
- If the metal frame is a precious metal alloy, apply a metal adhesive primer (eg. ALLOY PRIMER) or use a dental tin plating unit to coat it.
- To securely harden the Opaque resin and assure strong adhesion of the Opaque resin to the metal, first apply Opaque Primer to the adhering surface of the metal. Apply the Opaque resin after being sure the Opaque Primer has dried completely.
- The standard procedure is to apply the Opaque twice. Be sure the first layer of the Opaque has hardened completely before applying Opaque Primer on the first layer of the Opaque. Dry the second Primer coat, and then apply a second layer of the Opaque.
- Use Opaque Modifier for color adjustment of the Body Opaque.

### ④ Avoid separation of the Body resin from the EG Fiber


- Sandblast the EG Fiber with alumina particles, 50 µm diameter.
- Do not touch the EG Fiber with your bare hands or fingers when sandblasting or applying the primer; the fiber surface will be contaminated, preventing the EG Fiber from adhering to the Body.

### ⑤ Avoid roughness in the oral cavity

- When heat curing the pastes coats, be sure to use the heating temperature and length of time specified in these Instructions for Use.
- To minimize the formation of unpolymerized resin, which leads to a rough surface, heat curing should be performed in the 100°C - 110°C/212°F - 230°F range for 15 minutes.
- Before final light curing, apply the Air Barrier Paste to the restoration surface and keep it on until the heat curing is over.
- Shaping, finishing and polishing should be performed after heat curing. If heat curing is performed after shaping, finishing or polishing, insufficient surface hardness could result, thus leading to premature roughening of the restoration while in the oral cavity.
- Poor initial smoothness will cause premature surface roughness after the restoration is cemented in the oral cavity. Observe the following conditions:
  1. Remove rough scores and scratches with a silicone point before polishing.
  2. When polishing the restoration, use the polishing agent, brush and felt wheel supplied with the product.
  3. The restoration should be, and can be, polished until the surface is as smooth as glazed porcelain.

### ⑥ Avoid dislodged jacket crowns, bridges with frameworks, inlays and onlays

- Use a dental adhesive resin cement to lute jacket crowns, esthetic facing cast crowns, inlays, onlays or bridges with frameworks. Do not use zinc phosphate cement for restorations other than metal-core facing crowns.
- If a silicone-based material is to be used for preliminary test seating, use a polycondensation type to avoid hindering the bond between the inner surface of the crown and the abutment.  
After testing, clean the restoration and tooth surface with a cotton pledge moistened with alcohol or a tooth surface cleaner.
- To seat temporary restorations, use a temporary sealing or cementing agent that does not contain eugenol.
- Sandblast the inner surface of the jacket crown, bridge with a framework, inlay or onlay using low pressure (1 to 2 atmospheres) to prevent poor adhesion that could be caused by the separating agent, and as a pretreatment process.

## ⑦ Avoid poor polymerization

- The Opacious Dentin should be built up in thin layers and each layer light cured for the specified length of time; otherwise poor polymerization will result.
- After curing the EG Fiber using EG Core instrument, remove the EG Core instrument and then light cure the EG Fiber again to ensure optimum polymerization.
- In very deep cavity preparations, the inner surface at the floor of the cavity may not cure sufficiently. If the cavity is deeper than the polymerization depth indicated, build up the paste layers incrementally, and light cure each installation.

## ⑧ Avoid poor separation from plaster models

- Use the separating agent supplied with the product kit; otherwise the restoration may not release properly or the plaster model may break.
- Be careful not to let the cured resin extend over the margins; otherwise the restoration may not release properly or the plaster model may break during release. Scrape the cured resin from the margins before releasing the restoration from the model.

## ⑨ Avoid polymerization of pastes in ambient light

- The Body resins and EG Fiber can harden if they are used under strong light, such as by a window or under a laboratory light. Use the materials in more moderate light conditions.

## ⑩ Avoid entrapped air bubbles

- Apply Modeling Liquid on layered surfaces whenever adding layering Body resin.
- Do not mix different Body resins.
- Be careful to prevent air bubbles from being trapped between the EG Fiber and the Body resins when placing the Body resins on the EG Fiber.

## ⑪ Other points to avoid

- Do not use the same brush to apply Opaque Primer and Opaque resin; the Opaque resin will cure on the brush before application when the same brush is used. Brushes used to apply Opaque resin and Opaque Primer should be washed in ethyl alcohol or a self-curing resin monomer liquid after use.
- Be careful not to burn yourself when softening EG Core in hot water.

## Storage precautions

### ① Storage


- Do not leave the caps off the Opaque Primer, Add-On Primer and CR Sep III. After use, the cap should be replaced as soon as possible to prevent premature setting or contamination.
- The EG Fiber should be put into an aluminum bag very soon after use to prevent it from hardening.
- The Opaque Primer, Add-On Primer and CR Sep III are flammable. Do not use or store them near an open flame.
- The product should not be stored in direct sunlight or near an open flame. It should be stored at room temperature (2°C-25°C/39°F - 77°F).


### ② Expiration date

The product should be used by the expiration date\* indicated on the package.

[The expiration date is established by Kuraray Medical's own certification (from Kuraray Medical research data).]

\*(Example: Exp. 2005-11 means that the expiration date is Nov 2005.)

## IV Indications for use

### Indications

This material is indicated for use of the restoring crowns and defects.

### Restoration Applications

#### 1) Facing cast crowns and facing cast bridges

(|5 : Superstructure of implant)

#### 2) Jacket crowns

(|6 : Jacket crown)

#### 3) Inlays and onlays


(|56 : Inlays)

#### 4) Bridges with frameworks

(|567 : Crown bridge)

## V Kit components

### Body resin

The Body resin is a light/heat-cured, crown and bridge restorative material that has surface-treated glass powder and surface-treated alumina-based microfiller densely mixed in a matrix of multifunctional methacrylate monomers (polyurethane methacrylate monomer and other methacrylate monomers). The Body resin, when cured, has the mechanical strength and wears resistance to withstand occlusal pressures in the posterior area of the mouth.

### Opaque resin

Opaque resin is a single-paste, light-cured composite resin containing methacrylate monomers including Bis-GMA, surface-treated quartz powder, and surface-treated organic composite filler. The opaque nature of this resin makes it especially useful for screening out the unwanted glare of metal through the product.

### Opaque Primer

This primer contains the metal-adhesive phosphate monomer (MDP), methacrylate monomer, and solvent. It contributes to the hardening of Opaque resin.


### Modeling Liquid

Modeling Liquid contains methacrylate monomer, polyurethane methacrylate monomer, and MDP. The liquid is applied to instruments or a resin surface on which another layer of paste is to be applied, to make it easier to shape the paste.

### Add-On Primer


This primer contains a silane coupling agent and is used when adding on an additional amount of resin after shaping.

## V Kit components


### EG Fiber

This is a fiber framework material for use in bridge spans, consisting of methacrylate monomers including polyurethane methacrylate monomer, TEGDMA, surface-treated glass fiber, and surface-treated microfiller.

### EG Flow

EG Flow is a single-paste, light-cured composite resin containing methacrylate monomers such as polyurethane methacrylate monomer and TEGDMA. It is used for repairing frameworks during fabrication.

### Accessories

#### ● Air Barrier Paste

This paste contains a polymerization accelerator. It is applied to the surface of Body resin before polymerization, to enhance the polymerization of the resin.

#### ● Jacket Separate Kit

A plaster separating agent for jacket crowns. It consists of Jacket Spacer, which creates a small gap for the cement, and makes it easier to release the crown from the mold, and Margin Sep, which helps separate the plaster from the jacket crown.

#### ● CR Sep III

A plaster separating agent for inlays.


#### ● Polishing set

This polishing set is supplied with the ESTENIA™ C&B system. It consists of a polishing agent containing diamond powder, a polishing brush, and a felt wheel.

#### ● EG Core


EG Core is a core material with light transmission properties. It is used to fabricate bridge frameworks with EG Fiber.

#### ● Instruments

These instruments are used to work with the Body resins. The instruments have shapes that minimize paste sticking which helps during the crown shaping procedure.

## System components of ESTENIA™ C&B

The system components of ESTENIA™ C&B are shown on each package.

### Single items


The following products can be purchased separately.

#### ● Body resin

Shade group	Content	Shade
<b>Transparent</b>	(6.4 g/2.6 ml each)	T0,T1,T2,TLV
<b>Enamel</b>	(6.4 g/2.6 ml each)	E0,E1,E2,E3,E4
<b>Dentin</b>	(6.4 g/2.6 ml each)	DA1,DA2,DA3,DA3.5,DA4,DB1,DB2,DB3,DB4,DC1,DC2,DC3,DC4,DD2,DD3,DD4, DNW0,DNW0.5,DNP1.5,DNP2.5
<b>Cervical</b>	(6.4 g/2.6 ml each)	CE1,CE2,CE3,CE4,CE5,CE6,CE7,CE8
<b>Opacious Dentin</b>	(6.4 g/2.6 ml each)	ODA1,ODA2,ODA3,ODA3.5,ODA4,ODB1,ODB2,ODB3,ODB4,ODC1,ODC2,ODC3, ODC4,ODD2,ODD3,ODD4,ODN0,ODN0.5,ODNP1.5,ODNP2.5
<b>Cervical Transparent</b>	(6.4 g/2.6 ml each)	CT1,CT2,CT3,CT4,CT5,CT6
<b>Transparent Effect</b>	(6.4 g/2.6 ml each)	TG,AM,AM-Y,ME,TB
<b>Enamel Effect</b>	(6.4 g/2.6 ml each)	CE,CE-O,CE-Y,WE,MA-1,MA-2,IE
<b>Cervical Dentin Effect</b>	(6.4 g/2.6 ml each)	CDE1,CDE2,CDE3,CDE4,CDE5
<b>Gingival</b>	(6.4 g/2.6 ml each)	P1,P2,P3,P4,P5


#### ● Opaque resin

Shade group	Content	Shade
<b>Body Opaque</b>	(2.0 ml each)	OA1,OA2,OA3,OA3.5,OA4,OB1,OB2,OB3,OB4,OC1,OC2,OC3,OC4,OD2, OD3,OD4,ONW0,ONW0.5,ONP1.5,ONP2.5
<b>Cervical Opaque</b>	(2.0 ml each)	CO1,CO2,CO3,CO4,CO5
<b>Opaque Modifier</b>	(1.0 ml each)	W,PRP,DPG,YBR,O,BR,PUR

**V Kit components**



- |  |                                |
|--|--------------------------------|
| ● EG Fiber for anterior teeth (12 cm×2)  | ● Polishing set                |
| ● EG Fiber for posterior teeth (12 cm×2) | • Polishing compound (10 g)    |
| ● EG Flow (1.0 ml)                       | • Polishing brush (12 pcs.)    |
| ● Opaque Primer (9 ml)                   | • Felt wheel (12 pcs.)         |
| ● Modeling Liquid (6 ml)                 | ● Polishing compound (10 g)    |
| ● Add-On Primer (6 ml)                   | ● Polishing brush (12 pcs.)    |
| ● Jacket Separate Kit                    | ● Felt wheel (12 pcs.)         |
| •Jacket Spacer (5 ml)                    | ● EG Core (6 pcs.)             |
| •Margin Sep (5 ml)                       | ● Instrument No.1              |
| ● Jacket Spacer (5 ml)                   | ● Instrument No.2              |
| ● Margin Sep (5 ml)                      | ● Instrument No.3              |
| ● CR Sep III (20 ml)                     | ● Instrument set               |
| ● Air Barrier Paste (10 ml)              | ● Mixing plate (#954) (8 pcs.) |
|  | ● Small brush (#901)           |
|  | ● Small brush (#932)           |
|  | ● Small brush (#933)           |

## VI Basic laboratory procedures

### 1. Flowchart

The fabrication of restorations using ESTENIA™ C&B is outlined in the following chart.

## 2. Curing unit, polymerization time and polymerization depth

This system is designed to be cured by light and heat.

The following are the conditions for light- and heat-curing the ESTENIA™ C&B system: The recommended effective wavelength of light curing unit is 400 - 515nm and the showed unit is an example.

### 1) Curing unit and polymerization time

- **Light curing condition** [The figures in parentheses are given for inlays and onlays.]

Light curing unit	EG Fiber	EG Flow	Opaque	Body		(sec)
				Preliminary polymerization	Final polymerization	
Light Curing-300 (TOESCO)	270	90	180	30(270)	270(270)	

- **Heat curing condition**

100 - 110°C / 212 - 230°F for 15 minutes [Heat curing unit: Heat Curing-110 (TOESCO)]

### 2) Polymerization depth


- **Body resin**

Light curing unit	Polymerization time (sec)	Polymerization depth (mm)									
		T0	E1	DA3	CE1	ODA3	CT2	ME	WE	CDE4	P1
Light Curing-300 (TOESCO)	30	3.9	2.5	1.5	1.4	1.1	2.3	3.0	3.1	1.2	1.9
	270	7.9	4.9	3.6	2.6	1.6	3.9	6.1	5.1	2.0	3.5

- **Opaque resin**

Light curing unit	Polymerization time (sec)	Polymerization depth (mm)				
		OA1	OA4	CO1	W	YBR
Light Curing-300 (TOESCO)	180	0.3	0.2	0.1	1.0	0.1

(They are determined based on ISO10477)

## 3. Color matching

### 1) List of shades

The ESTENIA™ C&B system contains 16 shades that correspond to the Vita shades, plus 4 new shades. It is advisable to use shades as outlined in the following general list of shades to assure the creation of optimal color match for the patient.

	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4	NW0	NW0.5	NP1.5	NP2.5
<b>Body Opaque</b>	OA1	OA2	OA3	OA3.5	OA4	OB1	OB2	OB3	OB4	OC1	OC2	OC3	OC4	OD2	OD3	OD4	ONW0	ONW0.5	ONP1.5	ONP2.5
<b>Cervical Opaque</b>	—	CO1	—	CO2	—	CO3	—	CO4	CO5	—	—	—	—	—	—	—	—	—	—	
<b>Opaque Modifier</b>	W	P	RP	DP	G	YBR	O	BR	PUR											
<b>Opacious Dentin</b>	ODA1	ODA2	ODA3	ODA3.5	ODA4	ODB1	ODB2	ODB3	ODB4	ODC1	ODC2	ODC3	ODC4	ODD2	ODD3	ODD4	ODNW0	ODNW0.5	ODNP1.5	ODNP2.5
<b>Dentin</b>	DA1	DA2	DA3	DA3.5	DA4	DB1	DB2	DB3	DB4	DC1	DC2	DC3	DC4	DD2	DD3	DD4	DNW0	DNW0.5	DNP1.5	DNP2.5
<b>Cervical</b>	—	CE1	CE2	—	CE3	CE4	—	CE5	CE6	CE7	CE8	—	—	—	—	—	—	—	—	
<b>Enamel</b>	E1	E2	E4	E1	E3	E1	E3	E1	E3	E1	E3	E0	E1							
<b>Transparent</b>	T0	T1	T2	TLV																
<b>Cervical Transparent</b>	CT1	CT2	CT3	CT4	CT1	CT5	CT6	CT1	CT2	CT4	CT2	CT4	—	—	—	—	—	—	—	
<b>Cervical Dentin Effect</b>	CDE1	CDE2	CDE3	CDE4	CDE5															
<b>Enamel Effect</b>	CE	CE-O	CE-Y	WE	MA-1	MA-2	IE													
<b>Transparent Effect</b>	TG	AM	AM-Y	ME	TB															
<b>Gingival</b>	P1	P2	P3	P4	P5															

#### Using the Opaque resin

##### ● Body Opaque

Used for masking the color of the substrate and produce the basic shade of the crown.

A-line shades Body Opaque

##### ● Cervical Opaque

Used for the cervical area. Especially effective when the cervical area is narrow.

A and B-line shades Cervical Opaque


##### ● Opaque Modifier

Used for adjusting Opaque resin shades or used independently.

Available shade: W(white), P(pink), RP(regular pink), DP(dark pink), G(gray), YBR(yellow-brown), O(orange), BR(brown), PUR(purple)


## VI Basic laboratory procedures

### Using the Body resin

#### ● Dentin

Used to produce the basic color of the crown.

A-line shades Dentin

#### ● Opacious Dentin

Paste used to produce the basic color of the crown and to maintain the shade of the substrate.


A-line shades Opacious Dentin

#### ● Cervical


Used to produce cervical color.

A-line shades Cervical

#### ● Enamel


Paste used to produce enamel color.


#### ● Transparent

Paste used to produce transparency in cervical areas.


Shade	Used for producing:
T0	Enamel of high transparency
T1	Enamel of ordinary transparency
T2	Cloudy enamel
TLV	Slightly dark enamel

### ● Cervical Transparent


Used to produce the basic shade of the inlay. Used also to produce transparency in cervical areas.

### ● Cervical Dentin Effect

Used to characterize the cervical region. Usually layered on Cervical paste.


Shade	Used for producing:
CDE1	A-line shades
CDE2	B-line shades
CDE3	Metal frame
CDE4	Bright shades of A and B lines
CDE5	Further brighter shades than CDE4 such as whitening


### ● Enamel Effect

Paste used to characterize enamel

Shade	Used for producing:
CE	Crests on the occlusal surface of molars or enamel color of low transparency
CE-O	Orange enamel color of low transparency or incisal halo. It is also effective for producing colors in contact areas.
CE-Y	Yellow enamel color of low transparency or incisal halo. It is also effective for producing colors in contact areas.
WE	Cloudy enamel color seen at the angular areas of anterior teeth, white band, and the cusps of molars
MA-1	Bright mamelon color
MA-2	Mamelon color of slightly high chroma
IE	Transparency for a long abutment or incisal edge with metal backing


## VI Basic laboratory procedures

### ● Transparent Effect

Paste used for producing a variety type of crown transparency

Shade	Used for producing:
<b>ME</b>	Inlay margins
<b>AM</b>	Producing amber color enamel
<b>AM-Y</b>	Producing amber color enamel in which yellow is slightly stronger than AM
<b>TG</b>	Producing a gray transparent layer
<b>TB</b>	Producing a blue transparent layer


### ● Gingival

Paste used to produce gingival color

Shade	Used to produce:
<b>P1</b>	Pink
<b>P2</b>	Dark pink
<b>P3</b>	Orangish pink
<b>P4</b>	Reddish orange
<b>P5</b>	Violet

## 2) Color matching procedure

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When fabricating a jacket crown, you will be able to reproduce subtle shading if the core is made using composite resin. In such a case, Opacious Dentin is used to form the substrate. If the abutment is metal, it is advisable to use Body Opaque. You can reproduce the Vita shade you want by using the basic quadruple-layer method that layers the substrate shade (Body Opaque or Opacious Dentin), Dentin, Transparent, and Enamel. When enough space for forming is available, Cervical Transparent may be used to reproduce subtle shading (Applied method).

You can also use Transparent or Effect for the surface layer.

#### Jacket crown

##### ● Opacious Dentin (or Opaque resin)

Form the paste, about 0.2 mm thick, into a scaled-down contour of the crown. Use Opaque resin if the abutment is metal.

##### ● Dentin

Form the paste on the area from the margin to 1/3 of the crown length to the occlusal surface into the same as the final contour and make the paste thickness at the remaining 2/3 area gradually thinner to the occlusal surface.

##### ● Transparent



Form the paste thinly on an area of 1/4 - 1/5 of the crown length to the occlusal surface.

##### ● Cervical Transparent

If enough space for forming is available, this paste may be used thinly between Dentin and Enamel to reproduce subtle shading.

##### ● Enamel

Form the paste thinly on an area of 1/3 of the crown length to the occlusal surface.

#### Facing crown

##### ● Opaque resin

Apply the paste on the facing surface to mask out the metal color of the substrate.

##### ● Dentin

Form the paste on an area from the margin to 1/3 of the crown length to the occlusal surface into the same as the final contour. Make the paste thickness at the remaining 2/3 area gradually thinner to the occlusal surface.

##### ● Transparent




Form the paste on an area of 1/4 - 1/5 of the crown length to the occlusal surface.

##### ● Cervical Transparent

If enough space for forming is available, this paste may be used thinly between Dentin and Enamel to reproduce subtle shading.

##### ● Enamel

Form the paste thinly on an area of 1/3 of the crown length to the occlusal surface.

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When fabricating a jacket crown for the anterior region, you will be able to reproduce subtle shading if the core is made using composite resin, as in the case of the posterior restoration.

In such a case, Opacious Dentin is used for the substrate instead of Opaque resin. If the abutment is metal, it is advisable to use Opaque resin. Basically, you can reproduce the Vita shade you want by using the triple-layer method which layers the substrate shade (Opaque or Opacious Dentin), Dentin, and Enamel. You may use Cervical for the cervical area, Transparent for incisal edges, and Effect pastes when you think it necessary.

\*Enamel should not be used for more than 1/2 of the crown length at the marginal ridge or more than 1/3 of the crown length to the incisal edge at the center of the crown.


\*Cervical should be formed in such a way that the paste is made gradually thinner from the margin to 1/4 - 1/5 of the crown length to the cervical region.

The paste should be formed sparingly; otherwise you may not obtain the correct Vita shade.

### Jacket crown

- **Opacious Dentin (or Opaque resin)**

Form the paste, about 0.2 mm thick, into a scaled-down contour of the crown. Use Opaque resin if the abutment is metal.

- **Dentin**

Form the paste on the area from the margin to 1/3 of the crown length to the incisal edge into the same as the final contour and make the paste thickness for the remaining 2/3 area gradually thinner to the incisal edge.


- **Enamel**

Form the paste thinly on an area of 1/3 of the crown length to the incisal edge.

### Facing crown

- **Opaque resin**

Apply the paste on the facing surface to mask out the metal color of the substrate.


- **Dentin**

Form the paste on the area from the margin to 1/3 of the crown length to the incisal edge into the same as the final contour and make the paste thickness at the remaining 2/3 area gradually thinner to the incisal edge.

- **Enamel**

Form the paste thinly on an area of 1/3 of the crown length to the incisal edge.


## Inlay and onlay

### Basic method

#### ● Dentin

Form the paste in the proximal area if the cavity preparation contains a proximal surface. If the cavity preparation has a lining material at the bottom, make a thin layer of Dentin or Opacious Dentin on it.

#### ● Cervical Transparent

Form the paste from 0.2 mm inside the marginal line.

#### ● Enamel

Form the paste on the occlusal surface. Cervical Transparent

### Applied method (Subtle shade reproduction)

#### ● Cervical

Form the paste, about 0.2 mm thick, at the bottom of the cavity preparation. If the cavity preparation has a lining material in the bottom, make a thin layer of Dentin or Opacious Dentin on it.

#### ● Dentin

Form the paste at the proximal area if the cavity preparation contains a proximal surface.

#### ● Cervical Transparent

Form the paste from 0.2 mm inside the marginal line.

#### ● Transparent Effect (ME)


Form the paste thinly in a belt-like pattern, starting at the near-end of the margin.

#### ● Enamel

Form the paste on the occlusal surface.



#### ● Enamel Effect (CE)

Form the paste to enhance the brightness at the bulges of the occlusal surface.

## Bridges with frameworks

As in making a jacket crown, use Opaque resin for the substrate color if the abutment is metal or use Opacious Dentin if the abutment is composite resin. The subsequent procedure is the same as that for making a crown.

## 4. Bridge frameworks

Fabricate bridges with frameworks using EG Fiber. EG Fiber comes in two types: one type for anterior restorations and one type for posterior restoration frameworks. Use each type as is appropriate to the site you are going to restore.


### Posterior region framework

The framework should be designed so that the fiber curves from the occlusal surface of the abutment to the mucosal surface of the pontic. This curve strengthens the bridge. However, it is advisable to fabricate the framework so that there is about 1 mm of space between it and the mucosal surface, to prevent the EG Fiber from being exposed on the surface of the bridge. The EG Fiber should be far enough from the occlusal surface to allow a sufficient thickness of Body pastes to be applied.

- **Crown bridge**

- **Inlay bridge**

### Frameworks in the anterior region

The framework should be designed so that the fiber falls in the area from the lingual side of the abutment to the center of the pontic. It is preferable to locate it on the lingual side to the extent that the occlusal relationship is permits, in consideration of the esthetics of the labial side.


- **Crown bridge**

- **Inlay bridge**

\*Be careful not to expose EG Fiber on the restoration's surface when fabricating a framework.

## 5. Fabrication of jacket crowns

The product should be cured according to the section of Curing units, polymerization time and polymerization depth in this manual.

### 1 Fabrication of a working model

Fabricate a working model in the usual manner, and then trim the model.  
Do not trim the dies too much; too strong an undercut under the margin could cause the plaster model to break when the crown is released from the mold.  
It is advisable to block out with wax any area that has been trimmed too much.

### 2 Application of the separating agent

#### ① Applying and drying Jacket Spacer

To make it easy to release the crown from the mold, apply Jacket Spacer uniformly to all the area, except the margins, with a brush, and then dry the applied Jacket Spacer by leaving it or by blowing air gently over the coating.


#### ② Applying and drying Margin Sep

Apply Margin Sep on the margins of the abutment and dry by leaving it or by blowing air gently over it.

### 3 Forming and preliminary polymerization of Opacious Dentin (when abutment is composite resin)

When the abutment is composite resin, use Opacious Dentin for the substrate color.


### 4 Application and light curing of Opaque resin (when abutment is metal)

When the abutment is metal, use Opaque resin.

#### ① Applying and drying Opaque Primer

Apply Opaque Primer on the abutment tooth and evaporate the volatile content of the primer by mildly blowing with air or leaving it for about 30 seconds.


#### ② Applying and light curing Opaque resin

Apply Body Opaque paste of the shade specified on the abutment surface thinly and light cure the resin for the specified length of time using a laboratory-type light curing unit.

When you want to enhance the color in the cervical area, apply Cervical Opaque in the cervical area of the abutment.

Opaque Modifier may be used to adjust the Body Opaque color by mixing with Body Opaque paste or by applying it alone.

**VI Basic laboratory procedures****③Applying and light curing the 2nd layer of Opaque Primer and Opaque resin**

Apply Opaque Primer on the cured 1st layer of Body Opaque and dry it. Then, apply another layer of Body Opaque and light cure it in the same manner as for the 1st layer.

Repeat this process until the spacer color is masked out.

**5****Forming and preliminary polymerization of Body resin****①Forming and preliminarily light curing Cervical**

Form Cervical to characterize the cervical area as necessary.

Form the paste to an area from the margin to 1/4 - 1/5 of the crown length to the cervical area, making the paste thickness gradually thinner. Then, light cure the resin preliminarily for the specified length of time.

**②Applying Modeling Liquid**

Apply a thin coat of Modeling Liquid on the surface of the cured resin for better wetting and improved adaptation when layering another layer of paste.

Do not apply Modeling Liquid too much; otherwise Dentin paste might stick to the cured resin. Wipe away any excess liquid with tissues.

**③Forming and preliminarily light curing Dentin**

To produce the basic color of the crown, form a generous amount of Dentin paste and light cure it preliminarily for the specified length of time.

**④Forming and preliminarily light curing Transparent**


Apply Modeling Liquid to the surface of cured Dentin resin using caution not to entrap air bubbles as necessary. Then, apply a thin coat of Transparent paste to 1/4 to 1/5 of the crown length to the occlusal surface, and then light cure the paste preliminarily for the specified length of time.


**⑤Forming and preliminarily light curing Enamel**

Form Enamel paste to 1/3 of the crown length to the occlusal surface and then light cure the paste preliminarily for the specified length of time.

**6****Application of Air Barrier paste**


After forming Enamel paste in the contact area, apply Air Barrier paste on the surface layer of the crown to reduce the generation of unpolymerized resin. There may be a situation where some ingredients separate from others in the paste, but it does not have any adverse effect on the performance of the paste.

## 7 Final light curing

Light cure the pastes for the specified length of time

## 8 Heat curing

After releasing the crown from the mold, heat cure the restoration for 15 minutes at a temperature of 100°C - 110°C / 212°F - 230°F.

## 9 Shaping

### ① Adjustment of the proximal surface contact area

Adjust the proximal surface contact area using a silicone point or other instruments.

### ② Adjustment of the occlusion surface

Adjust the occlusion using a carborundum point or other instruments.

### ③ Adjustment of the crown contour

Adjust the crown contour using a carborundum point and a heatless stone.


### ④ Carving of pits and fissures

Carve pits and fissures in the occlusal surface using a carbide inverted cone and a white point.

## 10 Finishing

Remove scratches and grooves thoroughly from the crown surface using a silicone point.

## ⑪ Glazing

Glaze the restoration using the polishing set supplied with the kit.

### ① Polishing with brush

Apply the polishing compound to the brush and polish the occlusal surface and areas around it with the brush.

### ② Finishing


Apply the polishing compound to the felt wheel and glaze the entire crown using the felt wheel.

## ⑫ Completion

Sandblast the inner surface of the crown using alumina particles, 30 - 50 µm diameter, at a pressure of 1 to 2 kg/cm<sup>2</sup>


## 6. Fabrication of inlays and onlays

The product should be cured according to the section of Curing units, polymerization time and polymerization depth in this manual.

### 1 Fabrication of a working model

Fabricate a working model in the usual manner. If there is an undercut in the cavity, block it out using plaster or wax.

### 2 Application and drying of the resin separating agent


Apply CR Sep III to areas around the cavity, adjacent teeth or opposing tooth.

### 3 Forming and preliminary polymerization of Body resin


#### ①Application and preliminary polymerization of Dentin

If the cavity preparation has a proximal area, form Dentin paste in that proximal area and light cure the paste preliminarily for the specified length of time. Then, form Enamel paste and adjust the contour.

#### ②Application and preliminary polymerization of Cervical Transparent

Form Cervical Transparent paste on the occlusal surface of the cavity and then light cure the paste preliminarily for the specified length of time.


#### ③Forming and preliminary polymerization of Enamel

Form Enamel paste on the occlusal surface and light cure the paste preliminarily for the specified length of time.

### 4 Application of Air Barrier paste

After forming Enamel paste in the contact area, apply Air Barrier paste on the surface layer of the inlay to reduce the generation of unpolymerized resin. There may be a situation where some ingredients separate from others in the paste, but it does not have any adverse effect on the performance of the paste.

## 5 Final light curing

Light cure the pastes for the specified length of time.


## 6 Heat curing

After releasing the inlay from the mold, heat cure the restoration for 15 minutes at a temperature of 100°C - 110°C / 212°F - 230°F.

## 7 Shaping, finishing and glazing


As in making a jacket crown, shape, finish and glaze the restoration.

## 8 Completion

Remove the separating agent and sandblast the inner surface of the inlay using alumina particles, 30 - 50 µm diameter at a pressure of 1 to 2 kg/cm<sup>2</sup>, for pretreatment for adhesion.

## 7. Fabrication of bridges with frameworks

The product should be cured according to the section of Curing units, polymerization time and polymerization depth in this manual.

### ① Fabrication of a working model

Fabricate a working model in the usual manner. If there is an undercut, block it out with plaster or wax as necessary.

To provide a uniform space between the framework and the mucosal surface of the pontic, keep the pontic about 1 mm above the mucosal level.

### ② Making of the frame fabrication core

#### ①Waxing up

Using the sprue wax (about 2 mm diameter for the anterior region, about 2.5 - 3 mm diameter for the posterior region), reproduce the framework contour. The framework contour should be designed so that it is located under the pontic as much as possible. Check the occlusal relationship to be sure sufficient gap is attained.

#### ②Blocking out with silicone putty

To prevent the framework from displacing due to placement of EG Core, block out the area under the wax frame and the axial surface of the abutment with silicone putty.


#### ③Softening EG Core

Cut EG Core into a proper length and immerse it in a hot water bath (about 80°C/ 176°F) for about 3 minutes to soften the Core material.


#### ④Placing EG Core

Place softened EG Core onto the wax frame of the abutment.

#### ⑤Completion of the frame fabrication core

After curing, remove EG Core from the model.

### ③ Fabrication of the framework

#### ① Applying the separating agent

Apply CR Sep III on the plaster model.

#### ② Preparing EG Fiber

Cut EG Fiber into a size long enough to cover the entire occlusal surface of the bridge.

#### ③ Placing EG Fiber into the inner surface of the core

Place EG Fiber into the inner surface of EG Core or on the model.

#### ④ Placing the core under pressure

Return EG Core onto the model. After that, light cure the core preliminarily for 60 seconds using the laboratory-type light curing unit.


#### ⑤ Light curing the framework


Remove the core from the model and light cure it for the specified length of time to cure EG Fiber.

#### ⑥ Adjusting the framework

Remove EG Fiber from the model and trim using a carborundum point. Return EG Fiber onto the model and check the occlusal surface for proper clearance.

## ④ Application of the resin separating agent

### ① Applying and drying Jacket Spacer

To make it easy to release the bridge from the mold, apply Jacket Spacer uniformly to all the area except the margins with a brush, and then dry the applied Jacket Spacer by leaving it or by blowing air gently over the coating.

### ② Apply and drying Margin Sep

Apply Margin Sep on the margins of the abutment and the mucosal surface of the pontic, and dry it by leaving or by blowing air gently.

## ⑤ Designing of the framework

### ① Sandblasting

Sandblast the entire framework with alumina particles.

### ② Bonding treatment

Apply Add-On Primer on the framework surface, After drying it, apply Modeling Liquid.

### ③ Forming Opacious Dentin

Apply a thin coat of Opacious Dentin onto the basal surface of the pontic and light cure the paste for the specified length of time. Form Opacious Dentin thinly on the abutment. If the abutment is metal, it is advisable to apply Opaque Primer and Opaque resin, and light cure the pastes for the specified length of time in advance.

### ④ Placing the framework

Form Dentin or Opacious Dentin on the pontic area and place the framework on the abutment area on the model.



### ⑤ Applying EG Flow

Apply EG Flow on the connection areas of the frame and the joint part of the framework and paste.

### ⑥ Light curing

Light cure the paste for the specified length of time to secure the framework to the abutment.

## VI Basic laboratory procedures

## ⑥ Forming Body resin

## ① Forming and preliminarily light curing Dentin

Form a generous amount of Dentin paste to produce the basic color of the crown and light cure the paste preliminarily for the specified length of time.

## ② Characterization

Using CHROMA ZONE™ Color Stain, characterize the resin as necessary.

## ③ Forming and preliminarily light curing Transparent

Form the paste thinly to 1/4 to 1/5 of the crown length to the occlusal surface.

## ④ Forming and preliminarily light curing Enamel

Form the paste thinly to 1/3 of the crown length to the occlusal surface and light cure the paste preliminarily for the specified length of time.

## ⑦ Applying Air Barrier paste

After forming Enamel paste at the contact area, apply the Air Barrier paste on the surface layer of the bridge to reduce the generation of unpolymerized resin. There may be a situation where some ingredients separate from others in the paste, but it does not have any adverse effect on the performance of the paste.

## ⑧ Final light curing and heat curing

Light cure the restoration for the specified length of time. After releasing the bridge from the mold, heat cure the restoration for 15 minutes at a temperature of 100°C - 110°C / 212°F - 230°F.

## ⑨ Shaping, finishing and glazing


As in making a jacket crown, shape, finish and glaze the restoration using the polishing compound supplied with the kit.

## ⑩ Completion

As in making a jacket crown, sandblast the inner surface of the crown to complete the restoration.

## 8. Fabrication of facing cast crown (superstructure of implants)

The product should be cured according to the section of Curing units, polymerization time and polymerization depth in this manual.

### 1 Fabrication of a working model

Fabricate a working model in the usual manner.

### 2 Fabrication of metal frame

Fabricate a metal frame by waxing up and casting in the usual manner.

Apply retention beads with a diameter of 100 - 200 µm to the facing surface

### 3 Pretreatment of metal frame

#### ① Sandblasting

Sandblast the metal frame with alumina particles, 50 µm diameter. After that, clean the surface with ultrasound for 2 minutes and dry.

#### ② Pretreatment for metal adhesion

If a precious metal alloy is used, apply the metal adhesive primer (eg. Alloy Primer) to the facing surface of the metal frame before drying.

#### ③ Apply and light curing Opaque resin

##### 1) Apply and dry Opaque Primer

Apply Opaque Primer to the facing surface and dry the primer by mildly blowing with air or leaving it for about 30 seconds to evaporate the volatile content of the primer.


##### 2) Apply and light cure Opaque resin

Apply Body Opaque paste of the shade specified on the abutment surface thinly and light cure the resin for the specified length of time using a laboratory-type light curing unit.

When you want to enhance the color in the cervical area, apply Cervical Opaque in the cervical area of the abutment.

Opaque Modifier may be used to adjust the Body Opaque shade by mixing with a Body Opaque paste or by applying it alone.

## VI Basic laboratory procedures

**3) Apply and light cure the 2nd layer of Opaque Primer and Opaque resin**


Apply Opaque Primer on the cured 1st layer of Body Opaque and dry it. Then, apply another layer of Body Opaque and light cure it in the same manner as for the 1st layer. Repeat this process until the spacer color is masked out.

**4 Forming and polymerization of Body resin****①Forming and preliminarily light curing Cervical**

Form the paste to an area from the margin to 1/4 - 1/5 of the crown length to the cervical area, making the paste thickness gradually thinner. Then, light cure the paste for the specified length of time.

**②Applying Modeling Liquid**

Apply a thin coat of Modeling Liquid on the surface of the cured resin for better wetting and improved adaptation when layering another layer of paste. Do not apply Modeling Liquid too much; otherwise the Dentin paste might stick to the cured resin. Wipe away any excess liquid with tissues.

**③Forming and preliminarily light curing Dentin**

To produce the basic color of the crown, form a generous amount of Dentin paste and light cure it preliminarily for the specified length of time.

**④Forming and preliminarily light curing Transparent**

As necessary, apply Modeling Liquid to the surface of the cured Dentin resin using caution not to entrap air bubbles. Then, apply a thin coat of the Transparent paste to 1/4 to 1/5 of the crown length to the occlusal surface, and then light cure the paste preliminarily for the specified length of time.


**⑤Forming and preliminarily light curing Enamel**

Using CHROMA ZONE™ Color Stain, characterize the resin as necessary. Form Enamel paste to 1/2 to 1/3 of the crown length to the occlusal surface and then light cure the paste preliminarily for the specified length of time.

**5 Application of Air Barrier paste**


After forming Enamel paste at the contact area, apply Air Barrier paste on the surface layer of the crown to reduce the generation of unpolymerized resin. There may be a situation where some ingredients separate from others in the paste, but it does not have any adverse effect on the performance of the paste.


## 6 Final light curing and heat curing

After forming Transparent to the contact areas, light cure the paste for the specified length of time. After that, heat cure the restoration for 15 minutes at a temperature of 100°C - 110°C/ 212°F - 230°F.

## 7 Shaping, finishing, glazing and completion

As in making a jacket crown, shape, finish and glaze the facing crown using the polishing compound supplied with the kit to complete the facing crown.

## 9. Additional fabrication procedures

If an additional fabrication is needed after shaping and glazing, use the following procedures:

### 1 Exposure of a fresh surface where an additional forming is to be made

Using a carborundum point, cut the area where you want to add paste to expose a fresh surface. Remove cuttings by blowing with air.

### 2 Applying and drying Add-On Primer

Apply Add-On Primer to the area where you want to add paste.

### 3 Applying Modeling Liquid

Apply a thin coat of Modeling Liquid.

### 4 Forming, light curing and heat curing Body resin


As in the fabrication of a jacket crown, form and light cure Body resin and then heat cure it.

### 5 Completion

Shape and glaze the restoration.


## VII Clinical procedures



### 1. Preparing abutments and cavities

Use the following procedures when preparing abutments and cavities with ESTENIA™ C&B system.  
An insufficient thickness or a too sharp margin will lead to breakage.

#### Jacket crown



#### Inlay, onlay



## VII Clinical procedures

## Bridges with frameworks (Crown bridges)



### Preparation at the anterior region

### Preparation at the posterior region

## Inlay bridge

## 2. Loading crown restorations

When loading restorations fabricated using ESTENIA™ C&B system, use adhesive resin cements which contain fillers and feature excellent hardness and wear resistance, such as PANAVIA™ F 2.0. If a soft cement or non-adhesive cement is used, the restoration might break.

Restorations fabricated using ESTENIA™ C&B system need to be treated with a silane coupling agent when bonding. They do not bond well if an adhesive resin cement is used alone.

### Luting with PANAVIA™ F 2.0

#### Treating the inner surface of ESTENIA™ C&B restorations

##### ① Sandblasting


- 30 - 50 µm alumina powder, pressure of 1 - 2 kg/cm<sup>2</sup>

##### ② Acid etching

- Apply phosphoric acid treatment agent (eg. K-Etchant Gel).
- Leave it for 5 seconds.
- Wash with water and dry.

##### ③ Silane treatment

- Mix one drop each of CLEARFIL SE BOND™ Primer and Porcelain Bond Activator.
- Apply the mixture and leave it for 5 seconds.
- Dry by mildly blowing with air.

#### Preparing the abutment

##### ① Treating metal surface (when metal core is used)


- Apply metal-adhesive primer (eg. Alloy Primer).

- Dry.

##### ② Treating the tooth surface

- Mix one drop each of the ED Primer II components.

- Apply the mixture (to tooth, metal core, resin core).

- Leave it for 30 seconds.

- Dry completely by blowing with air.

#### Light cure

- Light cure each point for 20 seconds.

#### Chemical cure

- Apply OXYGUARD™ II and leave it for 3 minutes before washing.

### 3. Repairing crown restorations

If an ESTENIA™ C&B restoration breaks in the mouth and needs to be repaired, use the following procedures. Before repairing, check the cause of the breakage and give considerations to the occlusal relationship.

#### Repairing a crown restoration

##### ①Preparing a fractured surface

- Cut one layer of tooth structure from the fractured surface.
- Roughen the metal surface. (in the case of facing cast crown)


##### ⑤Applying CLEARFIL SE BOND™

- Apply the bond to the adherent surface.
- Even out the bond by mildly blowing with air.
- Light cure the bond for 10 seconds.

##### ②Acid etching the surface of the ESTENIA™ C&B restoration


- Apply phosphoric acid treatment agent (eg. K-Etchant Gel).
- Leave it for 5 seconds.
- Wash with water and dry.

##### ⑥Masking metal color (in the case of facing cast crown)

- Apply CLEARFIL ST™ Opaquer.
- Light cure the paste for 40 seconds.

##### ③Treating metal surface (when facing crown is made of precious metal)

- Apply the metal-adhesive primer (eg. Alloy Primer).
- Dry.


##### ⑦Filling, light curing and polishing


- Fill CLEARFIL AP-X™ or CLEARFIL ST™.
- Light cure and polish.



##### ④Silane treatment of adherent surface

- Mix one drop each of CLEARFIL SE BOND™ Primer and Porcelain Bond Activator.

- Apply the mixture and leave it for 5 seconds.
- Dry by mildly blowing with air.

## [WARRANTY]

KURARAY MEDICAL INC. will replace any product that is proved to be defective. KURARAY MEDICAL INC. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

## [NOTE]

ESTENIA, CLEARFIL, CLEARFIL SE BOND, CLEARFIL AP-X, CLEARFIL ST, PANAVIA, OXYGUARD and CHROMA ZONE are trademarks of KURARAY CO., LTD.

Manufactured by:

**KURARAY MEDICAL INC.**

1621 Sakazu, Kurashiki, Okayama 710-8622, Japan

European Union Representative

**KURARAY EUROPE GmbH**









Schiess-Strasse 68, 40549 Düsseldorf, Germany

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