CLEARFIL™ S³ BOND Universal

[**CLEARFIL TRI-S BOND Universal**]

CE
0197

EN	ENGLISH	English	DENTAL BONDING AGENT
FR	FRANÇAIS	French	AGENT DE LIAISON DENTAIRE
ES	ESPAÑOL	Spanish	AGENTE DE ADHESIÓN DE USO DENTAL
IT	ITALIANO	Italian	AGENTE DI ADESIONE DENTALE
NL	NEDERLANDS	Dutch	DENTALE BONDING
DE	DEUTSCH	German	Dentaler Haftvermittler
SV	SVENSKA	Swedish	DENTALT BONDINGSMATERIAL
NO	NORSK	Norwegian	DENTALT BINDINGSMIDDEL
FI	SUOMI	Finnish	HAMPAAN SIDOSAINE
DA	DANSK	Danish	DENTALT BONDINGMATERIALE
PT	PORTUGUÊS	Portuguese	AGENTE DE ADESÃO DENTÁRIO
EL	ΕΛΛΗΝΙΚΑ	Greek	ΟΔΟΝΤΙΚΟΣ ΣΥΓΚΟΛΛΗΤΙΚΟΣ ΠΑΡΑΓΟΝΤΑΣ
TR	TÜRKÇE	Turkish	DENTAL BONDING AJAN
PL	POLSKI	Polish	STOMATOLOGICZNY SYSTEM WIAŻĄCY
RO	ROMÂNĂ	Romanian	AGENT DENTAR ADEZIV
HR	HRVATSKI	Croatian	DENTALNI ADHEZIV
HU	MAGYAR	Hungarian	FOGÁSZATI KÖTŐANYAG
SL	SLOVENŠČINA	Slovenian	DENTALNO VEZIVNO SREDSTVO
CS	ČESKY	Czech	DENTÁLNÍ POJIVO
SK	SLOVENSKY	Slovak	DENTÁLNE TMELIACE ČINIDLO
BG	БЪЛГАРСКИ	Bulgarian	СТОМАТОЛОГИЧЕН БОНДИНГ АГЕНТ
UK	УКРАЇНСЬКА	Ukrainian	СТОМАТОЛОГІЧНИЙ АДГЕЗИВ
ET	EESTI	Estonian	DENTAALNE ADHESIIV
LV	LATVIEŠU	Latvian	ZOBĀRSTNIECĪBAS SAISTVIELA
LT	LIETUVIŠKAI	Lithuanian	ODONTOLOGINIS RIŠIKLIS

Kuraray Noritake Dental Inc.

1621 Sakazu, Kurashiki, Okayama 710-0801, Japan

EC REP

EMERGO EUROPE

Westervoortsedijk 60,
6827 AT Arnhem, The Netherlands

Kuraray Europe GmbH (EU Importer)

Philipp-Reis-Str. 4,
65795 Hattersheim am Main, Germany

Phone: +49 (0)69 305 35835 Fax: +49 (0)69 305 98 35835


URL: <https://www.kuraraynoritake.eu>

1562R675R

CLEARFIL™ S³ BOND Universal

[CLEARFIL TRI-S BOND Universal]

I. INTRODUCTION

"CLEARFIL TRI-S BOND Universal" is a single-component, light-cured bonding agent that allows simultaneous treatment of both dentin and enamel. Depending on the indication, this product is used as self-etching or with phosphoric acid for selective enamel etching or total-etching procedures. This product is intended to be used for both direct and indirect restorations. "CLEARFIL DC Activator" activates the dual-curing mechanism of this product; however, the addition of "CLEARFIL DC Activator" to this product is not required when using with "CLEARFIL DC CORE PLUS" or self-adhesive resin cements manufactured by Kuraray Noritake Dental Inc. (Kuraray's self-adhesive cements). The general clinical benefit of this product is to restore tooth function for the following INDICATIONS FOR USE.

II. INDICATIONS FOR USE

"CLEARFIL TRI-S BOND Universal" is indicated for the following uses:

- [1] Direct restorations using light-cured composite resin
- [2] Cavity sealing as a pretreatment for indirect restorations
- [3] Treatment of exposed root surfaces
- [4] Treatment of hypersensitive teeth
- [5] Intraoral repairs of fractured restorations
- [6] Post cementation and core build-ups
- [7] Cementation of inlays, onlays, crowns, bridges and veneers

III. CONTRAINDICATIONS

Patients with a history of hypersensitivity to methacrylate monomers

IV. POSSIBLE SIDE EFFECTS

The oral mucosal membrane may turn whitish when contacted by the product due to the coagulation of protein. This is usually a temporary phenomenon that will disappear in a few days. Instruct patients to avoid irritating the affected area while brushing.

V. INCOMPATIBILITIES

- [1] Do not use eugenol-containing materials for pulp protection or temporary sealing, since the eugenol can retard the curing process.
- [2] Do not use hemostatics containing ferric compounds, since these materials may impair adhesion and may cause discoloration of the tooth margin or surrounding gingiva, due to remaining ferric ions.
- [3] When using hemostatics containing aluminum chloride, minimize the quantity; take care to prevent contact with the adherend surface. Failure to do so might weaken the bond strength to the tooth structure.

VI. PRECAUTIONS

1. Safety precautions

1. This product contains substances that may cause allergic reactions. Avoid use of the product in patients with known allergies to methacrylate monomers or any other components.
2. If the patient demonstrates a hypersensitivity reaction such as rash, eczema, features of inflammation, ulcer, swelling, itching or numbness, discontinue use of the product and seek medical attention.
3. Take care to prevent the product from coming in contact with the skin or getting into the eye. Before using the product, cover the patient's eyes with a towel to protect them in the event of splashing material.
4. If the product comes in contact with human body tissues, take the following actions:
 - <If the product gets in the eye>
Immediately wash the eye with copious amounts of water and consult a physician.
 - <If the product comes in contact with the skin or the oral mucosa>
Immediately wipe the area with a cotton pellet or a gauze pad moistened with alcohol; rinse with copious amounts of water.
5. Take care to prevent the patient from accidentally swallowing the product.
6. Avoid looking directly at the dental curing light when curing the product.
7. Avoid using the same product dispensed into a well of the dispensing dish and the same applicator brush for different patients to prevent cross contamination. The applicator brush is single-use only. Discard the applicator brush after use.
8. Wear gloves or take other appropriate protective measures to prevent the occurrence of hypersensitivity that may result from contact with methacrylate monomers or any other components.
9. Dispose of this product as a medical waste to prevent infection.

2. Handling and manipulation precautions

[CLEARFIL TRI-S BOND Universal]

1. The product must not be used for purposes other than specified in [II.INDICATIONS FOR USE].
2. The use of this product is restricted to dental professionals.
3. The product contains ethanol, a flammable substance. Do not use it near an open flame.
4. Use a pulp capping agent in a cavity close to the pulp or in the event of accidental pulp exposure.
5. To prevent poor performance or poor handling characteristics, observe the specified light-curing times and other handling requirements.
6. Clean the cavity sufficiently to prevent poor bonding. If the adherend surface is contaminated with saliva or blood, wash it thoroughly and dry before bonding.
7. Use the light blocking plate to avoid exposing the material to an operating light or natural light, and use within the working time indicated below after dispensing or mixing:

Material	Working time
BOND	7 minutes
BOND + CLEARFIL DC Activator	90 seconds

As the volatile ethanol which is contained in BOND evaporates, the viscosity increases, thereby making it difficult to apply.

8. During the application of BOND to the adherend surface and the duration of rubbing for 10 seconds, move the lighting-spot out of the mouth or turn off the light to prevent the applied BOND from being exposed to the operating light. In addition, after treatment for 10 seconds, dry sufficiently as soon as possible by blowing mild air.
9. BOND contains ethanol and water. "CLEARFIL DC Activator" contains ethanol. Dry the entire adherend surface sufficiently by blowing mild air for more than 5 seconds until BOND or the mixture of BOND and "CLEARFIL DC Activator" does not move; otherwise, the adhesion effect will be impaired. In order to dry sufficiently, adjust the air pressure according to the shape and size of the cavity and the prosthetic appliance. Use a vacuum aspirator to prevent BOND or the mixture from scattering.
10. If the treated surface is contaminated, wash it with water, dry, or clean with alcohol, and treat again with the product.
11. Do not mix BOND with other bonding agents except "CLEARFIL DC Activator".
12. The container should be tightly capped immediately after use to reduce the evaporation of the volatile solvent (ethanol contained in the product). If the liquid does not flow easily from the nozzle, the nozzle may have plugged up. Do not dispense the liquid forcibly from the plugged container.
13. If the liquid may not flow easily, since it has not been used for a long term, shake the container before use.

[Dental light-curing unit]

1. Low light intensity causes poor adhesion. Check the lamp for service life and the dental curing light guide tip for contamination. It is advisable to check the dental curing light intensity using an appropriate light evaluating device at appropriate intervals.
2. The emitting tip of the dental curing unit should be held as near and vertical to the resin surface as possible. If a large resin surface has to be light-cured, it is advisable to divide the area into several sections and light-cure each section separately.

3. Storage precautions

1. The product must be used by the expiration date indicated on the package.
2. The product must be refrigerated (2-8°C/ 36-46°F) when not in use, and should be brought to room temperature for more than 15 minutes before using. After taking out of the refrigerator, the product must be left standing until it comes to room temperature; otherwise, the excessive amount of liquid might be dispensed or the liquid might ooze after use.
3. Keep away from extreme heat, direct sunlight or a flame.
4. The product must be stored in a proper place where only dental practitioners can access.

VII. COMPONENTS

Please see the outside of the package for contents and quantity.

1) BOND

Principal ingredients:

- 10-Methacryloyloxydecyl dihydrogen phosphate (MDP) (5-15%)
- Bisphenol A diglycidylmethacrylate (Bis-GMA) (25-50%)
- 2-Hydroxyethyl methacrylate (HEMA) (10-25%)
- Hydrophilic aliphatic dimethacrylate (1-10%)
- Colloidal silica (1-10%)
- Silane coupling agent (< 2%)
- dl-Camphorquinone (< 3%)
- Ethanol (10-25%)
- Water (10-25%)
- Phenyl bis(2,4,6-trimethylbenzoyl)-phosphine oxide (< 1%)

Units in parentheses are mass %.

2) Accessories

- Applicator brush (fine <silver>)
- Dispensing dish*
- Light blocking plate*
- *Consumables

VIII. CLINICAL PROCEDURES

A. Standard procedure I

- [1] Direct restorations using light-cured composite resin
- [2] Cavity sealing as a pretreatment for indirect restorations
- [3] Treatment of exposed root surfaces
- [4] Treatment of hypersensitive teeth

A-1. Isolation and Moisture control

Avoid contamination of the treatment area from saliva or blood to produce optimal results. A rubber dam is recommended to keep the tooth clean and dry.

A-2. Cavity or root surface preparations

Remove any infected dentin and prepare the cavity in the usual manner. When treating hypersensitive teeth, clean the root surface as usual. Then, clean thoroughly by spraying with water, and dry with air or with cotton pellets.

A-3. Pulp protection

Any actual or near pulp exposure could be covered with a hard setting calcium hydroxide material. There is no need for cement lining or basing. Do not use eugenol materials for pulp protection.

A-4. Pretreatment of tooth

Choose either etching procedure before applying BOND.

[NOTE]

When treating hypersensitive teeth, Choose the following A-4a before applying BOND.

A-4a. Self-etching procedure

Move to the section A-5 without etching by a phosphoric acid.

A-4b. Selective enamel etching procedure

Apply a phosphoric acid etching gel (e.g. K-ETCHANT Syringe) to the uncut and/or cut enamel. Leave it in place for 10 seconds, then rinse and dry.

A-4c. Total etching procedure

Apply a phosphoric acid etching gel (e.g. K-ETCHANT Syringe) to the entire cavity (enamel and dentin), leave it in place for 10 seconds, then rinse and dry.

A-5. Application of BOND

1. Dispense the necessary amount of BOND into a well of the dispensing dish immediately before application.

[CAUTION]

Use the light blocking plate to avoid exposing the material to an operating light or natural light, and use within 7 minutes after dispensing.

2. Apply BOND to the entire cavity wall with the applicator brush and rub it in for 10 seconds. Use caution not to allow saliva or exudate to contact the treated surfaces.
3. Dry the entire cavity wall sufficiently by blowing mild air for more than 5 seconds until BOND does not move. Use a vacuum aspirator to prevent BOND from scattering.
4. Light-cure BOND with a dental curing unit (see table "Dental curing unit and curing time").

Table: Dental curing unit and curing time

Type	Light source	Light Intensity	Light-curing time
Halogen	Halogen lamp	More than 400 mW/cm ²	10 seconds
LED	Blue LED*	800 – 1400 mW/cm ²	10 seconds
		More than 1500 mW/cm ²	5 seconds

The effective wavelength range of each dental curing unit must be 400 - 515 nm.

*Peak of emission spectrum: 450 - 480 nm

A-6. Placement of composite resin restorative or treatment of hypersensitive teeth

A-6a. Direct restorations using light-cured composite resin

Place composite resins (e.g. CLEARFIL Performance Pro, CLEARFIL MAJESTY ES Flow) into the cavity, light-cure, finish and polish according to the manufacturer's instructions.

A-6b. Cavity sealing and treatment of exposed root surfaces

Place a thin coat of composite resins (e.g. CLEARFIL MAJESTY ES Flow) onto the tooth, and light-cure according to the manufacturer's instructions. Remove unpolymerized resin with a cotton pellet or a gauze moistened with alcohol.

A-6c. Treatment of hypersensitive teeth

Remove the unpolymerized layer of BOND with a cotton pellet or a gauze moistened with alcohol.

B. Standard procedure II

[5] Intraoral repair of fractured restorations

B-1. Preparation of adherend surfaces

Roughen the adherend surfaces using a diamond point or blasting with 30 to 50 µm alumina powder at air pressure of 0.1-0.4 MPa (14-58 PSI/ 1-4 kgf/cm²). The air pressure should be properly adjusted to suit the material and/or shape of the prosthetic restoration, using caution to prevent it from chipping. And place a bevel at the marginal area.

B-2. Phosphoric acid treatment of adherend surfaces

Apply a phosphoric acid (e.g. K-ETCHANT Syringe) to the adherend surface (including tooth structure). Leave in place for 5 seconds, then rinse and dry.

B-3. Application of BOND

Apply BOND to the entire adherend surfaces. Refer to section A-5.

[NOTE]

For optimal performance, a silane coupling agent (e.g. CLEARFIL CERAMIC PRIMER) to the surface of the silica-based ceramics before applying BOND can be used according to the manufacturer's instructions.

B-4. Place composite resin restorative

Place composite resins (e.g. CLEARFIL Performance Pro) into the cavity, light-cure, finish and polish according to the manufacturer's instructions.

[NOTE]

Use an opaque resin (e.g. CLEARFIL ST OPAQUER) to mask metal color.

C. Standard procedure III

[6] Post cementation and core build-ups

When using "CLEARFIL DC CORE PLUS", the use of "CLEARFIL DC Activator" is not necessary.

C-1. Isolation and Moisture control

Avoid contamination of the treatment area from saliva or blood to produce optimal results. A rubber dam is recommended to keep the tooth clean and dry.

C-2. Preparing the root canal

Prepare and clean the root canal opening in the usual manner.

C-3. Post preparation

Choose either C-3a or C-3b based on the post you use. Please follow the Instructions for Use of the restoration material. In the absence of specific instructions, we recommend the following procedure:

C-3a. For Glass Fiber Posts

Apply a phosphoric acid (e.g. K-ETCHANT Syringe) to the post surface. Leave it in place for 5 seconds, then rinse and dry.

[CAUTION]

- Do not blast glass fiber posts with alumina powder, or the posts could be damaged.

- Avoid all contamination of the surfaces to be treated during pre-treatment and until the final core build-up.

C-3b. For Metal Posts

Roughen the adherend surface by blasting with 30 to 50 µm alumina powder at an air pressure of 0.1-0.4 MPa (14-58 PSI/ 1-4 kgf/cm²). The air pressure should be properly adjusted to suit the material. After blasting, clean the metal post by using ultrasound for 2 minutes followed by drying it with an air stream.

C-4. Post surface treatment

Choose either procedure based on the material you use.

C-4a. When using with CLEARFIL DC CORE PLUS

1. Prepare BOND according to section A-5.1.

2. Apply BOND to the entire post surface with the applicator brush and leave it in place for 5 seconds.

3. Dry the entire adherend surface sufficiently by blowing mild air for more than 5 seconds until BOND does not move. Use a vacuum aspirator to prevent BOND from scattering.

[NOTE]

To dry thoroughly, adjust the air pressure according to the shape and size of the adherend surface.

C-4b. When using with other dual- /self-cured resin core material or a dual- /self-cured resin cement

1. Dispense one drop each of BOND and "CLEARFIL DC Activator" into a well of the dispensing dish and mix them together with the applicator brush.

[CAUTION]

Use the light blocking plate to avoid exposing the material to an operating light or natural light, and use within 90 seconds after mixing.

2. Apply the mixture to the post surface, and leave it in place for 5 seconds.

3. Dry the entire adherend surface sufficiently by blowing mild air for more than 5 seconds until the mixture does not move. Use a vacuum aspirator to prevent the mixture from scattering.

[NOTE]

To dry thoroughly, adjust the air pressure according to the shape and size of the adherend surface.

4. Light-cure the mixture with a dental curing unit (see table "Dental curing unit and curing time" in A-5).

[CAUTION]

Working time will be shortened dramatically when not light-curing the mixture on the post.

C-5. Pretreatment of tooth

As necessary, apply a phosphoric acid. Refer to the section A-4.

C-6. Bonding

Choose either procedure based on the material you use.

C-6a. When using with CLEARFIL DC CORE PLUS

1. Apply BOND to the entire root canal with the applicator brush and rub it in for 10 seconds. Use caution not to allow saliva or exudate to contact the treated surfaces.

2. Dry the entire adherend surface sufficiently by blowing mild air for more than 5 seconds until BOND does not move. Use a vacuum aspirator to prevent BOND from scattering. Remove the excess mixture with a paper point. After wiping away the excess mixture, dry the adherend surface again, if necessary.

[NOTE]

To dry thoroughly, adjust the air pressure according to the shape and size of the adherend surface.

3. Light-cure BOND with a dental curing unit (see table "Dental curing unit and curing time" in A-5).

C-6b. When using with other dual- /self-cured resin core material or a dual- /self-cured resin cement

1. Apply the mixture to the entire root canal with the applicator brush and rub it in for 10 seconds. Use caution not to allow saliva or exudate to contact the treated surfaces.

2. Dry the entire adherend surface sufficiently by blowing mild air for more than 5 seconds until the mixture does not move. Use a vacuum aspirator to prevent the mixture from scattering. Remove the excess mixture with a paper point. After wiping away the excess mixture, dry the adherend surface again, if necessary.

[NOTE]

To dry thoroughly, adjust the air pressure according to the shape and size of the adherend surface.

3. Light-cure the mixture with a dental curing unit (see table "Dental curing unit and curing time" in A-5).

[CAUTION]

Working time will be shortened dramatically when not light-curing the mixture on the root canal.

C-7. Placing the post and core build-up

Place the post and core build-up using "CLEARFIL DC CORE PLUS" or other resin material according to the manufacturer's instructions.

D. Standard procedure IV

[7] Cementation of inlays, onlays, crowns, bridges and veneers

When using with Kuraray's self-adhesive cements, the use of "CLEARFIL DC Activator" is not necessary.

D-1. Conditioning the cavity and core (tooth, metal, composite) surfaces

1. Remove the temporary sealing material and temporary cement in the usual manner, and clean the cavity with moisture control.

2. Trial fit the prosthetic restoration to check its fit on the cavity or core (tooth, metal, composite). When using a try-in paste to check the color, follow the manufacturer's instructions.

D-2. Surface preparation of prosthetic restorations

Choose either D-2a or D-2b based on the restoration you use. Follow the Instructions for Use of the restoration material. In the absence of specific instructions, we recommend the following procedure:

D-2a. For silica-based glass ceramics (e.g. Lithium disilicate)

Etch the glass ceramic surfaces with hydrofluoric acid in accordance with manufacturer's instructions, and thoroughly wash and dry the surface.

D-2b. For Metal oxides (e.g. Zirconia), metals or composite resins

Roughen the adherend surface by blasting with 30 to 50µm alumina powder at an air pressure of 0.1-0.4 MPa (14-58 PSI/ 1-4 kgf/cm²). The air pressure should be properly adjusted to suit the material and/ or shape of the prosthetic restoration, using caution to prevent chipping. After blasting, clean the prosthetic restoration by using ultrasound for 2 minutes followed by drying it with an air stream.

D-3. Pretreatment of prosthetic restorations

Choose either procedure based on the material you use.

D-3a. When using with Kuraray's self-adhesive cements

1. Prepare BOND according to section A-5.1.
2. Apply BOND to the entire adherend surface with the applicator brush and leave it in place for 5 seconds.
3. Dry the entire adherend surface sufficiently by blowing mild air for more than 5 seconds until BOND does not move. Use a vacuum aspirator to prevent BOND from scattering.

[NOTE]

- To dry thoroughly, adjust the air pressure according to the shape and size of the adherend surface.
- Application of BOND to metal oxides or metals is optional since Kuraray's self-adhesive cements bonds strongly to these surfaces.

D-3b. When using with other dual- or self-cured resin cement

1. Dispense one drop each of BOND and "CLEARFIL DC Activator" into a well of the dispensing dish and mix them together with the applicator brush.

[CAUTION]

Use the light blocking plate to avoid exposing the material to an operating light or natural light, and use within 90 seconds after mixing.

2. Apply the mixture to the adherend surface, and leave it in place for 5 seconds.
3. Dry the entire adherend surface sufficiently by blowing mild air for more than 5 seconds until the mixture does not move. Use a vacuum aspirator to prevent the mixture from scattering.

[NOTE]

To dry thoroughly, adjust the air pressure according to the shape and size of the adherend surface.

4. Light-cure the mixture with a dental curing unit (see table "Dental curing unit and curing time" in A-5).

[CAUTION]

Working time will be shortened dramatically when not light-curing the mixture on the adherend surface.

D-4. Pretreatment of tooth

As necessary, apply a phosphoric acid. Refer to the section A-4.

D-5. Bonding

Choose either procedure based on the material you use.

D-5a. When using with Kuraray's self-adhesive cements

1. Apply BOND to the entire cavity wall with the applicator brush and rub it in for 10 seconds. Use caution not to allow saliva or exudate to contact the treated surfaces.
2. Dry the entire adherend surface sufficiently by blowing mild air for more than 5 seconds until BOND does not move. Use a vacuum aspirator to prevent BOND from scattering.

[NOTE]

- To dry thoroughly, adjust the air pressure according to the shape and size of the adherend surface.
- Light-curing of BOND is an option for optimal performance. Refer to table "Dental curing unit and curing time" in A-5.

D-5b. When using with other dual- or self-cured resin cement

1. Apply the mixture to the entire cavity wall with the applicator brush and rub it in for 10 seconds. Use caution not to allow saliva or exudate to contact the treated surfaces.
2. Dry the entire adherend surface sufficiently by blowing mild air for more than 5 seconds until the mixture does not move. Use a vacuum aspirator to prevent the mixture from scattering.

[NOTE]

To dry thoroughly, adjust the air pressure according to the shape and size of the adherend surface.

3. Light-cure the mixture with a dental curing unit (see table "Dental curing unit and curing time" in A-5).

[CAUTION]

Working time will be shortened dramatically when not light-curing the mixture on the adherend surface.

D-6. Cementing

Cement the prosthetic restoration using Kuraray's self-adhesive cements or other resin cement according to the manufacturer's instructions.

[WARRANTY]

Kuraray Noritake Dental Inc. will replace any product that is proven to be defective. Kuraray Noritake Dental Inc. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

[NOTE]


If a serious incident attributable to this product occurs, report it to the manufacturer's authorized representative shown below and the regulatory authorities of the country in which the user/patient resides.

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"CLEARFIL", "CLEARFIL MAJESTY", "CLEARFIL ST" and "CLEARFIL DC CORE PLUS" are registered trademarks or trademarks of KURARAY CO., LTD.