

TREATMENT OF A YOUNG PATIENT WITH ZIRCONIA VENEERS



MDT DANIELE
RONDONI



DR. ENZO
ATTANASIO

Veneers made of zirconia? In some cases, like the one presented below, monolithic zirconia veneers may be an option. Reasons for selecting a latest-generation zirconia such as “KATANA™ Zirconia” YML include its very high translucency and a wall thickness of only 0.3 to 0.4 mm supporting minimally invasive tooth preparation. Due to a highly automated production procedure, the manual effort involved may be reduced, while highly aesthetic outcomes are possible.



Fig. 1: Initial situation: Young female patient with misshaped and misaligned maxillary incisors.



Fig. 2: Digital smile design revealing the ideal proportions and positions of the anterior teeth.



Fig. 3: Ideal tooth proportions and positions displayed over a picture of the teeth after orthodontic treatment and the creation of a mock-up. The positions are ideal and the tooth shapes obtained with the mock-up only need some minor adjustments.



Fig. 4: Detailed view of the maxillary teeth with the ideal shapes blended in. The decision is made to produce six veneers made of "KATANA™ Zirconia" YML for definitive treatment.



Fig. 5a: Digital impression of the maxillary teeth with mock-up imported into the design software.




Fig. 5b: Digital model of the maxilla in the CAD software with full-contour veneers designed for the incisors and canines according to the ideal shapes planned with digital smile design.




Fig. 6a: Occlusal view of the virtual model.




Fig. 6b: Occlusal view of the model with the designed veneers.



Fig. 7a: Facial view of the patient with the mock-up in place.



Fig. 7b: Facial view of the patient with the planned veneers blended in.



Fig. 8: Guided tooth structure removal with the aid of a silicone index. The minimum wall thickness of the selected material – “KATANA™ Zirconia” YML – is 0.4 mm.



Fig. 9: Maxillary teeth immediately after minimally-invasive tooth preparation.




Fig. 10: Matched digital impressions of the maxilla and mandible taken after the tooth preparation.




Fig. 11: Digital impression: occlusal view of the maxillary teeth.




Fig. 12: Virtual veneers ready to be milled.



Fig. 13: Pre-sintered veneers made of “KATANA™ Zirconia” YML.



Fig. 14: Zirconia veneers on the resin model after sintering.



Fig. 15: Monolithic zirconia restorations on the resin model after the 7-hour final sintering.



Fig. 16: Frontal view of the master cast with the six veneers individualized with the liquid ceramic system CERABIEN™ FC Paste Stain.



Fig. 17: Lateral view of the veneers on the master cast.



Fig. 18: Tooth-like translucency of the veneers on the model.



Fig. 19: The minimal wall thickness of 0.4 to 0.6 mm and the positioning of the veneers in the upper (enamel) part of the blank are responsible for this effect.



Fig. 20: Veneers made of "KATANA™ Zirconia" YML ready for try-in.



Fig. 21: Intra-oral try-in with two different shades of the PANAVIA™ V5 Try-in Paste: A2 is used in the right and Clear in the left quadrant.




Fig. 22: Cleaning of the tooth structure with pumice paste after try-in.



Fig. 23: Thorough rinsing with water to remove any paste residues.



Fig. 24: Etching of the prepared tooth structure with phosphoric acid etchant.



Fig. 25: Clean, isolated central incisors ready for adhesive cementation.



Fig. 26: Sandblasting of the veneer's intaglio with alumina oxide (particle size: 50 µm, pressure: 2 bar). This kind of surface conditioning treatment is recommended for zirconia, as it creates the desired micro-retentive surface structure favourable for adhesive bonding¹.



Fig. 27: Application of CLEARFIL™ CERAMIC PRIMER PLUS to the bonding surface of the veneers (followed by mild air-drying).



Fig. 28: Treatment of the etched tooth structure with PANA VIA™ V5 Tooth Primer (followed by mild air-drying).



Fig. 29: Application of PANA VIA™ V5 Paste (shade A2) to the bonding surface of the veneers.




Fig.30: Close-up view of the restoration margin during cementation. Working field isolation is a key success factor of adhesive procedures.



Fig. 31: Lateral view of the cemented veneers. The result is a natural surface texture, which contributes to a natural appearance of the restorations.



Fig. 32: Frontal view of the veneers in place.



Fig. 33: Treatment outcome immediately after rubber dam removal.

FINAL SITUATION



Fig. 34: Treatment outcome with healthy soft tissues two weeks after treatment.



Fig. 35: Gums are healthy and the restorations show a great optical integration with the adjacent posterior teeth.

References

¹ Nishigawa G, Maruo Y, Irie M, Maeda N, Yoshihara K, Nagaoka N, Matsumoto T, Minagi S. Various Effects of Sandblasting of Dental Restorative Materials. PLoS One. 2016 Jan 14;11(1):e0147077.

“KATANA™ ZIRCONIA” YML

VISIT THE PRODUCT PAGE ONLINE FOR MORE DETAILS

WHAT YOU CAN EXPECT TO FIND ONLINE

- ✓ General information about “KATANA™ Zirconia” YML
- ✓ Full Product Assortment
- ✓ Safety Data Sheet



KURARAYNORITAKE.EU/KATANA-ZIRCONIA-YML