



Kuraray Noritake Dental Inc.
1621 Sakazu, Kurashiki,
Okayama 710-0801, Japan

INDICATIONS

- | | | |
|--|---|---|
| [1] Direct restorations using light-cured composite resin | → | Direct restoration |
| [2] Sealing of a prepared cavity or abutment tooth as a pretreatment for indirect restorations | → | Sealing of a prepared abutment tooth |
| [3] Treatment of exposed root surfaces* | → | |
| [4] Treatment of hypersensitive teeth* | → | |
| [5] Intraoral repairs of fractured restorations | → | Intraoral repair |
| [6] Post cementation and core build-ups | → | Post cementation / Core build-up |
| [7] Cementation of indirect restorations | → | Cementation |

* Please refer to the Instructions for Use for [3] and [4] of indications.

Table 1: Dental curing unit and curing time

| Type | Light source | Light Intensity | Light-curing time |
|---------|--------------|-----------------------------------|-------------------|
| Halogen | Halogen lamp | More than 400 mW/cm ² | 10 seconds |
| | | 800-1400 mW/cm ² | 10 seconds |
| LED | Blue LED* | More than 1500 mW/cm ² | 5 seconds |

The effective wavelength range of each dental curing unit must be 400-515nm. * Peak of emission spectrum: 450-480nm

If the treated surface is contaminated, KATANA™ Cleaner can be selected to clean the adherent surfaces. When using KATANA™ Cleaner, follow the Instructions for Use.



Direct restoration using light-cured composite resin

Follow the standard procedures for isolation, moisture control, cavity preparation and pulp protection

- Tooth pretreatment**
Choose either etching procedure
 a. Self-etching (Move to section 2)
 b. Selective-etching
Apply K-ETCHANT Syringe to the uncut and/or cut enamel, then rinse and dry 10sec.
 c. Total-etching
Apply K-ETCHANT Syringe to the entire cavity (enamel and dentin), then rinse and dry 10sec.
- Apply BOND with a rubbing motion
No waiting time
- Dry by blowing mild air until BOND does not move*1
5sec. +
- Light-cure*2
- Place composite resin, light-cure and finish

*1 Use a vacuum aspirator to prevent BOND from scattering.
*2 Refer to Table 1 for light-curing time.

Sealing of a prepared abutment tooth as a pretreatment for indirect restorations

Follow the standard procedure for isolation, moisture control, preparation of abutment tooth

- Tooth pretreatment**
Choose either etching procedure
 a. Self-etching (Move to section 2)
 b. Selective-etching*1
 c. Total-etching*1
 10sec.
- Apply BOND with a rubbing motion
No waiting time
- Dry by blowing mild air until BOND does not move*2
5sec. +
- Light-cure*3
- If necessary, place a thin coat of composite resin (e.g. CLEARFIL MAJESTY ES Flow) onto the tooth, light-cure according to the manufacturer's instructions.
- Wipe the surface to remove the un-polymerized layer (oxygen inhibited layer)*4

*1 Refer to **Direct Restoration**.
*2 Use a vacuum aspirator to prevent BOND from scattering.
*3 Refer to Table 1 for light-curing time.
*4 Use a cotton pellet or a gauze moistened with alcohol.

Intraoral repair of fractured restorations

- Roughen, rinse and air dry
- Apply K-ETCHANT Syringe, then rinse and dry*1
5sec.
- Apply BOND with a rubbing motion
No waiting time
- Dry by blowing mild air until BOND does not move*2
5sec. +
- Light-cure*3
- Place composite resin*4, light-cure and finish

*1 This acid etching is not necessary for non-precious metal and metal oxide ceramic.
*2 Use a vacuum aspirator to prevent BOND from scattering.
*3 Refer to Table 1 for light-curing time.
*4 Use an opaque resin (e.g. CLEARFIL ST OPAQUER) to mask metal color.

Post cementation / Core build-ups with CLEARFIL DC CORE PLUS

Follow the standard procedures for isolation, moisture control and preparation of root canal and cavity

- Post pretreatment**
 For glass fiber post
 [1] Apply K-ETCHANT Syringe, rinse and dry 5sec.
 [2] Apply BOND, then dry by blowing mild air 5sec. Dry
 For metal post
 [1] Blast with alumina powder, then ultrasonic clean and dry
- Tooth pretreatment**
Choose either etching procedure
 a. Self-etching (Move to section 3)
 b. Selective-etching*1
 c. Total-etching*1
 10sec.
- Apply BOND with a rubbing motion
No waiting time
- Dry by blowing mild air and paper point until BOND does not move*2
5sec. +
- Light-cure*3
- Post cementation and core build-up using CLEARFIL DC CORE PLUS according to the manufacturer's instructions

*1 Refer to **Direct Restoration**.
*2 Use a vacuum aspirator to prevent BOND from scattering.
*3 Refer to Table 1 for light-curing time.

Post cementation / Core build-ups with other core material (except for CLEARFIL DC CORE PLUS)

Follow the standard procedures for isolation, moisture control and preparation of root canal

1 Post pretreatment

For glass fiber post

[1] Apply K-ETCHANT Syringe, rinse and dry **5sec.**

[2] Apply the mixture of BOND and CLEARFIL DC Activator*1, then dry by blowing mild air **5sec. Dry**

[3] Light-cure*2 **Note**

For metal post

[1] Blast with alumina powder, then ultrasonic clean and dry **5sec.**

2 Tooth pretreatment

Choose either etching procedure

a. Self-etching (Move to section 3)

b. Selective-etching*3 **10sec.**

c. Total-etching*3 **10sec.**

3 Apply the mixture*1 with a rubbing motion

No waiting time

4 Dry by blowing mild air and paper point until the mixture does not move*4

5sec. +

Note

5 Light-cure*2

Note

Note: Working time will be dramatically shortened when not light-curing

*1 Dispense one drop each of BOND and CLEARFIL DC Activator and mix them.

*2 Refer to Table 1 for light-curing time.

*3 Refer to Direct Restoration

*4 Use a vacuum aspirator to prevent the mixture from scattering.

6 Post cementation and build-up by core material according to the manufacturer's instructions

Cementation of indirect restorations with Kuraray's self-adhesive cements

Clean and dry the tooth surface, and then trial fit the prosthetic restoration

1 Surface preparation of prosthetic restorations

Silica-based glass ceramic (e.g. lithium disilicate)

Apply a hydrofluoric acid, then wash and dry

Metal oxide ceramics (e.g. zirconia), metals or composite resins

Blast with alumina powder (30-50µm/ 0.1-0.4MPa/ 14-58 PSI/ 1-4 bar), then ultrasonic clean and dry

2 Tooth Pretreatment

Choose either etching procedure

a. Self-etching (Move to section 3)

b. Selective-etching*1 **10sec.**

c. Total-etching*1 **10sec.**

3 Apply BOND with a rubbing motion

No waiting time

4 Dry by blowing mild air until BOND does not move*2

5sec. +

*2 Use a vacuum aspirator to prevent BOND from scattering.

5 Cementation using Kuraray's self-adhesive cements according to the manufacturer's instructions

Note: When using a partial light-curing (or "Tack-Cure") technique, the setting time of the excess cement will be shorter.

*1 Refer to Direct Restoration

Cementation of indirect restorations with self-adhesive resin cement without any specific instructions to pretreatment the adherent surface

Clean and dry the tooth surface, and then trial fit the prosthetic restoration

1 Surface preparation of prosthetic restorations

Silica-based glass ceramic (e.g. lithium disilicate)

Apply a hydrofluoric acid, then wash and dry

Metal oxide ceramics (e.g. zirconia), metals or composite resins

Blast with alumina powder (30-50µm/ 0.1-0.4MPa/ 14-58 PSI/ 1-4 bar), then ultrasonic clean and dry

2 Apply the mixture of BOND and CLEARFIL DC Activator*1, then dry until the mixture does not move*2

5sec. + Dry

Note

Note: Working time will be dramatically shortened when not light-curing

*1 Dispense one drop each of BOND and CLEARFIL DC Activator and mix them.

*2 Use a vacuum aspirator to prevent the mixture from scattering.

*3 Refer to Table 1 for light-curing time.

3 Light-cure*3

Note

4 Tooth Pretreatment

Choose either etching procedure

a. Self-etching (Move to section 5)

b. Selective-etching*4 **10sec.**

c. Total-etching*4 **10sec.**

5 Apply BOND with a rubbing motion

No waiting time

6 Dry by blowing mild air until BOND does not move*5

5sec. +

Note

Note: Working time will be dramatically shortened when not light-curing

*5 Use a vacuum aspirator to prevent BOND from scattering.

7 Light-cure*3

Note

8 Cementation using resin cement according to the manufacturer's instructions

Note: When using a partial light-curing (or "Tack-Cure") technique, the setting time of the excess cement will be shorter.

*4 Refer to Direct Restoration