



INDICATIONS

- | | CASE |
|--|--------------------------------------|
| [1] Direct restorations using light-cured composite resin | Direct restoration |
| [2] Sealing of a prepared cavity or abutment tooth as a pretreatment for indirect restorations | Sealing of a prepared abutment tooth |
| [3] Treatment of exposed root surfaces* | |
| [4] Treatment of hypersensitive teeth* | |
| [5] Intraoral repairs of fractured restorations | Intraoral repair |
| [6] Post cementation and core build-ups | Post cementation / Core build-up |
| [7] Cementation of indirect restorations | Cementation |

* Please refer to the Instructions for Use for [3] and [4] of indications.

Table 1: Dental curing unit and curing time


Type	Light source	Light Intensity	Light-curing time
Halogen	Halogen lamp	More than 400 mW/cm ²	10 seconds
		800-1400 mW/cm ²	10 seconds
LED	Blue LED*	More than 1500 mW/cm ²	5 seconds

The effective wavelength range of each dental curing unit must be 400-515 nm. * Peak of emission spectrum: 450-480 nm

If the treated surface is contaminated, KATANA™ Cleaner can be selected to clean the adherent surfaces. When using KATANA™ Cleaner, follow the Instructions for Use.



Preparation



* When snapping off the container cap, do not tilt it to avoid spilling BOND.

Direct restoration using light-cured composite resin

Follow the standard procedures for isolation, moisture control, cavity preparation and pulp protection

<p>1 Tooth pretreatment Choose either etching procedure</p> <p>a. Self-etching (Move to section 2)</p> <p>b. Selective-etching Apply K-ETCHANT Syringe to the uncut and/or cut enamel, then rinse and dry 10sec.</p> <p>c. Total-etching Apply K-ETCHANT Syringe to the entire cavity (enamel and dentin), then rinse and dry 10sec.</p>	<p>2 Apply BOND with a rubbing motion No waiting time</p>	<p>3 Dry by blowing mild air until BOND does not move*1 5sec. +</p>	<p>4 Light-cure*2</p>	<p>5 Place composite resin, light-cure and finish</p>
--	---	---	------------------------------	--

*1 Use a vacuum aspirator to prevent BOND from scattering. *2 Refer to Table 1 for light-curing time.

Sealing of a prepared abutment tooth as a pretreatment for indirect restorations

Follow the standard procedure for isolation, moisture control, preparation of abutment tooth

<p>1 Tooth pretreatment Choose either etching procedure</p> <p>a. Self-etching (Move to section 2)</p> <p>b. Selective-etching*1 10sec.</p> <p>c. Total-etching*1 10sec.</p>	<p>2 Apply BOND with a rubbing motion No waiting time</p>	<p>3 Dry by blowing mild air until BOND does not move*2 5sec. +</p>	<p>4 Light-cure*3</p>	<p>6 Wipe the surface to remove the un-polymerized layer (oxygen inhibited layer)*4</p>
--	---	---	------------------------------	--

*1 Refer to [Direct Restoration](#). *2 Use a vacuum aspirator to prevent BOND from scattering. *3 Refer to Table 1 for light-curing time. *4 If necessary, place a thin coat of composite resin (e.g. CLEARFIL MAJESTY ES Flow) onto the tooth, light-cure according to the manufacturer's instructions. *4 Use a cotton pellet or a gauze moistened with alcohol.

Intraoral repair of fractured restorations

<p>1 Roughen, rinse and air dry</p>	<p>2 Apply K-ETCHANT Syringe, then rinse and dry*1 5sec.</p>	<p>3 Apply BOND with a rubbing motion No waiting time</p>	<p>4 Dry by blowing mild air until BOND does not move*2 5sec. +</p>	<p>5 Light-cure*3</p>	<p>6 Place composite resin*4, light-cure and finish</p>
--	--	---	---	------------------------------	--

*1 This acid etching is not necessary for non-precious metal and metal oxide ceramic. *2 Use a vacuum aspirator to prevent BOND from scattering. *3 Refer to Table 1 for light-curing time. *4 Use an opaque resin (e.g. CLEARFIL ST OPAQUER) to mask metal color.

Post cementation / Core build-ups with CLEARFIL DC CORE PLUS

Follow the standard procedures for isolation, moisture control and preparation of root canal and cavity

1 Post pretreatment

For glass fiber post

[1] Apply K-ETCHANT Syringe, rinse and dry **5sec.**

[2] Apply BOND, then dry by blowing mild air **5sec. Dry**

For metal post

[1] Blast with alumina powder, then ultrasonic clean and dry

2 Tooth pretreatment

Choose either etching procedure

a. Self-etching (Move to section 3)

b. Selective-etching*¹ **10sec.**

c. Total-etching*¹ **10sec.**

*¹ Refer to [Direct Restoration](#)

3 Apply BOND with a rubbing motion

No waiting time

4 Dry by blowing mild air and paper point until BOND does not move*²

5sec.+

5 Light-cure*³

6 Post cementation and core build-up using CLEARFIL DC CORE PLUS according to the manufacturer's instructions

*² Use a vacuum aspirator to prevent BOND from scattering. *³ Refer to Table 1 for light-curing time.

Cementation of indirect restorations with Kuraray's self-adhesive cements

Clean and dry the tooth surface, and then trial fit the prosthetic restoration

1 Surface preparation of prosthetic restorations

Silica-based glass ceramic (e.g. lithium disilicate)

Apply a hydrofluoric acid, then wash and dry

Metal oxide ceramics (e.g. zirconia), metals or composite resins

Blast with alumina powder (30-50µm/ 0.1-0.4MPa/ 14-58 PSI/ 1-4 bar), then ultrasonic clean and dry

2 Tooth Pretreatment

Choose either etching procedure

a. Self-etching (Move to section 3)

b. Selective-etching*¹ **10sec.**

c. Total-etching*¹ **10sec.**

*¹ Refer to [Direct Restoration](#)

3 Apply BOND with a rubbing motion

No waiting time

4 Dry by blowing mild air until BOND does not move*²

5sec.+

5 Cementation using Kuraray's self-adhesive cements according to the manufacturer's instructions

*² Use a vacuum aspirator to prevent BOND from scattering.

Note: When using a partial light-curing (or "Tack-Cure") technique, the setting time of the excess cement will be shorter.