



• Read the IFU (Instructions for Use) and Technical Instructions before the procedure.

### ► Line-Ups

<b>Shede Base</b>	10g, 50g, 200g	SBA1 SBC3	SBA2 SBC4	SBA3 SBD2	SBA3.5 SBD3	SBA4 SBD4	SBB1 SBNW0	SBB2 SBNW0.5	SBB3 SBNP1.5	SBB4 SBNP2.5	SBC1	SBC2
	10g, 50g	SB White										
<b>Margin</b>	10g	MA1 MD4	MA2 MNW0	MA3 MNW0.5	MA3.5 MNP1.5	MA4 MNP2.5	MB2 M Clear	MB3 M Orange	MB4 M Peach	MC2 MDL	MC4 MRP	MD3
<b>Opacios Body</b>	10g, 50g, 200g	OBA1 OBC3 OB1110 OB3110 OB4225	OBA2 OBC4 OB1120 OB3120 OB5110	OBA3 OBD2 OB2015 OB3130 OB5120	OBA3.5 OBD3 OB2025 OB3215 OB5130	OBA4 OBD4 OB2110 OB3225	OBB1 OBNP1.5 OB2120 OB4015	OBB2 OBNP2.5 OB2130 OB4025	OBB3 OB White OB2215 OB4110	OBB4 OB Orange OB2225 OB4120	OBC1 OB Pale Pink OB3015 OB4130	OBC2 OBE OB3025 OB4215
<b>Cervical</b>	10g, 50g, 200g	CV-1	CV-2	CV-3	CV-4							
<b>Body</b>	10g, 50g, 200g	A1B C3B EWB 1110B 3110B 4225B	A2B C4B EWYB 1120B 3120B 5110B	A3B D2B 2015B 3130B 5120B	A3.5B D3B 2025B 3215B 5130B	A4B D4B 2110B 3225B	B1B NW0B 2120B 4015B	B2B NW0.5B 2130B 4025B	B3B NP1.5B 2215B 4110B	B4B NP2.5B 2225B 4120B	C1B EW00B 3015B 4130B	C2B EW0B 3025B 4215B
<b>Enamel</b>	10g, 50g, 200g	E1	E2	E3	Silky E1	Silky E2						
<b>Translucent</b>	10g, 50g, 200g	Tx	T0	T1	T2							
<b>Luster</b>	10g, 50g, 200g	LT0 Creamy White	LT1 ELT1	LT Natural ELT2	LT Super Gray ELT3	T Blue	Aqua Blue 1	Aqua Blue 2	LT Yellow	Incisal Aureola	Sun Bright	Creamy Enamel
<b>Clear Cervical</b>	10g, 50g, 200g	CCV-1	CCV-2	CCV-3	CCV-4							
<b>Mamelon</b>	10g, 50g	Mamelon 1	Mamelon 2									
<b>Tissue</b>	10g, 50g	Tissue 1	Tissue 2	Tissue 3	Tissue 4	Tissue 5	Tissue 6	Tissue 7				
<b>Modifier</b>	10g, 50g	White	Gray	Blue	Green	Yellow	Light Orange	Orange	Brown	Pink	Dark Pink	Coral Pink
<b>Add-on</b>	10g, 50g	AD-T	AD-B									
<b>Internal Stain</b>	3g	White Salmon Pink	Gray Red	Incisal Blue 1 A+	Incisal Blue 2 B+	Mamelon Orange 1 C+	Mamelon Orange 2 D+	Cervical 1 Bright	Cervical 2 Fluoro	Cervical 3	Earth Brown	Reddish Brown
<b>External Stain</b>	3g	Pure White Cervical 3	Gray Earth Brown	Black Reddish Brown	Blue Salmon Pink	Green 1 Pink	Green 2 Red	Yellow A+	Orange 1 B+	Orange 2 C+	Cervical 1 D+	Cervical 2
	10g, 30g	Glaze										
	10g	FL Glaze	VC Glaze									
<b>FC Paste Stain</b>	3g	White Mamelon Orange 2 C+	Value Cervical 1 D+	Light Gray Cervical 2 Fluoro	Dark Gray Cervical 3	Black Earth Brown	Grayish Blue Reddish Brown	Blue Salmon Pink	Green 1 Pink	Green 2 Red	Yellow A+	Mamelon Orange 1 B+
	5g	Glaze	Clear Glaze									

<b>Liquid</b>	CERABIEN™ ZR FORMING LIQUID	100ml	MEISTER LIQUID	100ml	IS LIQUID	10ml / 100ml	ES LIQUID	10ml / 100ml
---------------	-----------------------------	-------	----------------	-------	-----------	--------------	-----------	--------------

### ► Color Combination Table

Target shade	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4	NP1.5	NP2.5	EW00	EW0	EW	EWY	NW0	NW0.5
<b>Shede Base</b>	SBA1	SBA2	SBA3	SBA3.5	SBA4	SBB1	SBB2	SBB3	SBB4	SBC1	SBC2	SBC3	SBC4	SBD2	SBD3	SBD4	SBNP1.5	SBNP2.5	SBWhite			SBB1	SBNW0	SBNW0.5
<b>Margin</b>	MA1	MA2	MA3	MA3.5	MA4	MB2+MDL (1:1)	MB2	MB3	MB4	MC2+MDL (1:1)	MC2	MC4+MDL (1:1)	MC4	MD3+MDL (1:1)	MD3	MD4	MNP1.5	MNP2.5	MDL	MNW0+MDL (1:1)	MNW0+MDL (2:1)	MNW0	MNW0.5	
<b>Opacios Body</b>	OBA1	OBA2	OBA3	OBA3.5	OBA4	OBB1	OBB2	OBB3	OBB4	OBC1	OBC2	OBC3	OBC4	OBD2	OBD3	OBD4	OBNP1.5	OBNP2.5	-	-	-	-	-	-
<b>Cervical (+Body)</b>	-	CV-1+A2B (1:2)	CV-1+A3B (1:1)	CV-1+A3.5B (1:1)	CV-1	-	CV-2+B2B (1:2)	CV-2+B3B (1:1)	CV-2	-	CV-3+C2B (1:2)	CV-3+C3B (1:1)	CV-3	CV-4+D2B (1:2)	CV-4+D3B (1:1)	CV-4	-	CV-1+NP2.5B (1:2)	-	-	-	-	-	-
<b>Body</b>	A1B	A2B	A3B	A3.5B	A4B	B1B	B2B	B3B	B4B	C1B	C2B	C3B	C4B	D2B	D3B	D4B	NP1.5B	NP2.5B	EW00B	EW0B	EWB	EWYB	NW0B	NW0.5B
<b>Enamel</b>	E2		E3			E1	E2	E3		E2	E3		E2		E3		E2		Silky E1			Silky E2		
<b>Translucent / Luster</b>	T1/LT1																ELT2		ELT1			T1/LT1		

**► Color Combination Table** • Converting VITA™ 3D-Master™ Shades to Noritake Value Shades

VITA™ 3D-Master™ Shade	OM1	OM2	OM3	1M1	1M2	2L1.5	2L2.5	2M1	2M2	2M3	2R1.5	2R2.5	3L1.5	3L2.5	3M1	3M2	3M3	3R1.5	3R2.5	4L1.5	4L2.5	4M1	4M2	4M3	4R1.5	4R2.5	5M1	5M2	5M3	
Noritake Value Shade	-	-	-	NV1110	NV1120	NV2015	NV2025	NV2110	NV2120	NV2130	NV2215	NV2225	NV3015	NV3025	NV3110	NV3120	NV3130	NV3215	NV3225	NV4015	NV4025	NV4110	NV4120	NV4130	NV4215	NV4225	NV5110	NV5120	NV5130	
The Ideal Frama Shades	Non Colored (White)			A1, B1, B2									A2, A3						C1, D2						C1,D2					
Shade Base *1	SBNW0	SBNW0	SBNW0+SBNW0.5 (1:1)	SBC1	SBA1	SBB2	SBB2	SBC1	SBB2	SBB3	SBA1	SBA2	SBB3+SBC2 (1:1)	SBB3	SBC1	SBB2	SBB3	SBA1	SBA3	SBC2	SBB3	SBC2	SBB3	SBB4	SBA2	SBB4	SBC2	SBA4	SBA4	
Margin	MNW0	MNW0.5	MNW0.5	MA1	MA1	MB2	MB2+MD4 (1:1)	MC1	MB2	MB3	MA2	MA2	MC2	MB3+MD4 (1:1)	MC2	MB3	MB4	MA3+MC2 (1:1)	MA3+MC2 (1:1)	MA4+MC2 (1:1)	MA4+MB3 (1:1)	MC2	MA4	MA4	MC2+MD3 (1:1)	MA4	MA4	MA4	MA4	
Body	NW0B	NW0.5B	NW0.5B	1110B	1120B	2015B	2025B	2110B	2120B	2130B	2215B	2225B	3015B	3025B	3110B	3120B	3130B	3215B	3225B	4015B	4025B	4110B	4120B	4130B	4215B	4225B	5110B	5120B	5130B	
Opacios Body	-	-	-	OB1110	OB1120	OB2015	OB2025	OB2110	OB2120	OB2130	OB2215	OB2225	OB3015	OB3025	OB3110	OB3120	OB3130	OB3215	OB3225	OB4015	OB4025	OB4110	OB4120	OB4130	OB4215	OB4225	OB5110	OB5120	OB5130	
Enamel	Silky E2			E1									E2																	
Translucent / Luster	T1/LT1																													

VITA™, 3D-Master™ and OM1-5M3 are trademarks of VITA Zahnfabrik H. Rauter GmbH &amp; Co. \*1 In case of masking Zirconia frameworks

**► Baking Schedule**

	Baking Schedule	Porcelain type	Dry-Out Time	Predrying Temperature		Start Vacuum		Heat Rate		Vacuum Level	Release Vacuum		High Temperature		Hold Time		Cool Time		
				min.	°C	°F	°C	°F	°C/min.		°F/min.	kPa*1	°C	°F	°C	°F		min.	min.
Restorations using a zirconia framework or Anatomical zirconia coping	1	Margin	5 *2	600 *2	1112 *2	600	1112	50	90	96	1000	1832	1000	1832	1	1	4		
	2	Wash bake	5	600	1112	600	1112	45	81	96	930	1706	930	1706	-	1	4		
	3	Shade Base	5	600	1112	600	1112	45	81	96	930	1706	930	1706	-	1	4		
	4	Opacious Body *3	7-10	600	1112	600	1112	45	81	96	930-940	1706-1724	930-940	1706-1724	-	1	4		
	5	Cervical / Body / Enamel *3	7-10	600	1112	600	1112	45	81	96	930-940	1706-1724	930-940	1706-1724	-	1	4		
	6	Intenal Stain	5	600	1112	-	-	50	90	0	-	-	900	1652	-	-	4		
	7	Translucent / Luster / Clear Cervical / Tissue *3	7-10	600	1112	600	1112	45	81	96	930-940	1706-1724	930-940	1706-1724	-	1	4		
	8	Minor Adjustment *4	7	600	1112	600	1112	45	81	96	930	1706	930	1706	-	0.5	4		
	9	Self Glaze	5	600	1112	-	-	50	90	0	-	-	930	1706	-	0.5	4		
	10	External Stain (including Glaze) / CZR PRESS Glaze *5	5	600	1112	-	-	50	90	0	-	-	930	1706	-	-	4		
	11	FC Paste Stain (including Glaze, Clear Glaze) *6	5	600	1112	-	-	45	81	0	-	-	910	1670	-	-	4		
	12	MRP & Add-on	5	600	1112	-	-	45	81	0	-	-	880	1616	-	-	4		
	13	CZR PRESS Shade Base Stain 1st *7	5	700	1292	700	1292	65	117	96	1090	1994	1090	1994	-	1	4		
	14	CZR PRESS Shade Base Stain 2nd *7	5	700	1292	700	1292	65	117	96	1080	1976	1080	1976	-	1	4		
Full-zirconia-restorations	15	CZR PRESS Glaze *5	5	600	1112	600	1112	65	117	96	850	1562	850	1562	-	1	4		
	16	External Stain	5	600	1112	-	-	50	90	0	-	-	850	1562	-	-	4		
	17	External Stain FL Glaze, VC Glaze	5	600	1112	600	1112	65	117	96	850	1562	850	1562	-	1	4		
	18	FC Paste Stain	5	500	932	600	1112	45	81	96	750	1382	750	1382	-	1	4		

**Attention** The above program is only a guideline. Baking temperature may vary due to different furnace.

\*1 96 kPa = 72 cmHg (29 inchesHg)

\*2 When Magic Former is used, the dry-out time should be 7 min, predrying temperature 700°C(1292°F).

\*3 The baking conditions are established on the assumption that a single crown is baked. Adjust the dry-out time and baking temperature when baking bridges.

\*4 This is baking temperature for cases when Opacious Body, Cervical, Body, Enamel, Translucent, Luster, Clear Cervical, Mamelon, Modifier and Tissue are applied additionally.

\*5 CZR PRESS Glaze is a CZR PRESS LF item.

\*6 On the occasion of use to CZR Press or CZR Press LF surface, please refer to CERABIEN™ ZR FC Paste Stain IFU (Instructions for Use).

\*7 CZR PRESS Shade Base Stain is a CZR PRESS LF item.


**Kuraray Noritake Dental Inc.**

 300 Higashiyama, Miyoshi-cho  
 Miyoshi, Aichi 470-0293, Japan  
<http://www.kuraraynoritake.com>