CLEARFIL™ AP-X PLT

EN	ENGLISH	English	DENTAL LIGHT- CURED COMPOSITE RESTORATIVE MATERIAL
FR	FRANÇAIS	French	MATÉRIAU DE RESTAURATION DENTAIRE COMPOSITE PHOTOPOLYMERISABLE
ES	ESPAÑOL	Spanish	MATERIAL DE RESTAURACIÓN COMPOSITE FOTOPOLIMERIZABLE DE USO DENTAL
IT	ITALIANO	Italian	MATERIALE COMPOSITO FOTOPOLIMERIZZABILE PER RESTAURI DENTALI
NL	NEDERLANDS	Dutch	DENTAAL LICHTUITHARDEND RESTAURATIEF COMPOSITET MATERIAAL
DE	DEUTSCH	German	DENTALES lichthärtendes Komposit-Restaurationsmaterial
SV	SVENSKA	Swedish	DENTALT LJUSHÄRDANDE KOMPOSIT RESTAURERINGSMATERIAL
NO	NORSK	Norwegian	DENTALT LYSHERDET KOMPOSIT-RESTAURERINGSMATERIALE
FI	SUOMI	Finnish	HAMPAAN VALOKOVETTEINEN YHDISTELMÄMUOVI
DA	DANSK	Danish	DENTALT LYSHÆRDENDE KOMPOSIT-RESTAURERINGSMATERIALE
PT	PORTUGUÊS	Portuguese	MATERIAL DE RESTAURAÇÃO DE COMPÓSITO FOTOPOLIMERIZAVEL PARA USO DENTARIO
EL	ΕΛΛΗΝΙΚΑ	Greek	ΟΔΟΝΤΙΚΟ ΦΩΤΟΠΟΛΥΜΕΡΙΖΟΜΕΝΟ ΥΛΙΚΟ ΑΠΟΚΑΤΑΣΤΑΣΗΣ ΣΥΝΘΕΤΗΣ ΡΗΤΙΝΗΣ
TR	TÜRKÇE	Turkish	IŞIKLA SERTLEŞEN DENTAL KOMPOZİT RESTORATİF MATERYAL
ET	EESTI	Estonian	DENTAALNE VALGUSKÕVASTUV KOMPOSIIT-RESTAURATSIOONIMATERJAL
LV	LATVIEŠU	Latvian	ZOBĀRSTNIECĪBĀ LIETOJAMS AR GAISMAS PALĪDZĪBU CIETĒJOŠS RESTAURĀCIJAS MATERIĀLS
LT	LIETUVIŠKAI	Lithuanian	ODONTOLOGINIO ŠVIESA KIETINAMO KOMPOZITO RESTAURACINĖ MEDŽIAGA

Kuraray Noritake Dental Inc.

1621 Sakazu, Kurashiki, Okayama 710-0801, Japan

EMERGO EUROPE

Westervoortsedijk 60,
6827 AT Arnhem, The Netherlands

Kuraray Europe GmbH (EU Importer)


Philipp-Reis-Str. 4,
65795 Hattersheim am Main, Germany

Phone:+49 (0)69 305 35835 Fax:+49 (0)69 305 98 35835

URL:<https://www.kuraraynoritake.eu>

CLEARFIL™ AP-X PLT

I. INTRODUCTIONS

CLEARFIL AP-X PLT is a light-cured, radiopaque restorative composite resin which provides accurate color-matching, high polishability and excellent physical properties, making it ideal for both anterior and posterior restorations. It is formulated with optimal viscosity assuring easy handling and placement. CLEARFIL AP-X PLT, with the special dispensing system, can be quickly placed directly into the cavity.

It has radiopacity equal to or greater than 1mm aluminum, and is classified as a Type 1 and Class 2 (Group 1) material by ISO 4049.

The general clinical benefit of this product is to restore tooth function for the following INDICATIONS FOR USE.

II. INDICATIONS FOR USE

CLEARFIL AP-X PLT is indicated for the following restorative applications

1. Class I, II, V restorations of posterior teeth
2. Class III, IV, V restorations of anterior teeth
3. Cervical cavities or root surfaces defects

III. CONTRAINDICATION

Patients with a history of hypersensitivity to methacrylate monomers

IV. INCOMPATIBILITY

Do not use eugenol-containing materials for pulp protection or temporary sealing since the eugenol could retard the curing process.

V. PRECAUTIONS

1. Safety precautions

1. This product contains substances that may cause allergic reactions. Avoid use of the product in patients with known allergies to methacrylate monomers or any other components.
2. If any hypersensitivity, such as a rash or dermatitis, occurs from contact with the product, discontinue the use of the product and consult a physician.
3. Wear gloves or take other appropriate protective measures to prevent the occurrence of hypersensitivity that may result from contact with methacrylate monomers.
4. Use caution to prevent the product from coming in contact with the skin or getting into the eyes.
Before using the product, cover the patient's eyes with a towel to protect the patient's eyes from possible splashing material.
5. If the product comes in contact with human body tissues, take the following actions:
<If the product gets in the eyes>
Immediately wash the eye with copious amounts of water and consult a physician.
<If the product comes in contact with the skin>
Immediately wipe the material off using a cotton pledget (or gauze) moistened with alcohol and then wash with copious amounts of water.
6. Use caution to prevent the patient from accidentally swallowing the product.
7. Avoid looking directly at the curing light when curing the product.
8. For infection control reasons, PLT tips are for single use only.
9. Dispose of this product as a medical waste to prevent infection.

2. Handling and manipulation precautions

1. If necessary, use a rubberdam in order to prevent moisture contamination
2. Do not use the product in conjunction with the other composite resin. Mixing materials may cause a change in physical properties, possibly a decrease, from the properties expected.
3. If the adherent surface extends to uncut enamel, apply a phosphoric acid etching agent (e.g. K-ETCHANT GEL) and leave it in place for 10 seconds before washing and drying.
4. Amalgam or other lining material remaining in the cavity will prevent the passage of light and the polymerization of the product. Completely remove any lining material when preparing the cavity.
5. When light curing the product, note the depth of cure in this Instructions for Use.
6. The emitting tip of the dental curing unit should be held as near and vertical to the resin surface as possible. If a large resin surface is to be light cured, it is advisable to divide the area into several sections and light-cure each section separately.
7. Low light intensity causes poor adhesion. Check the lamp for service life and the dental curing unit guide tip for contamination. It is advisable to check the dental curing light intensity using an appropriate light-evaluating device at appropriate intervals.
8. The paste of the product should be used as soon as practicable after being dispensed from the PLT tip. If the paste is to be left for a while before use, it should be covered with a Light blocking plate.
9. The product should be returned to room temperature before dispensing if it has been taken out of a refrigerator. Failure to do this might cause breakage of the PLT tip.
10. Squeeze the dispenser with a slow, steady pressure. Excessive force is not necessary.
11. The use of this product is restricted to dental professionals.

3. Storage precautions

1. The product must be used by the expiration date indicated on the package.
2. Keep away from extreme heat or direct sunlight.
3. The product must be stored at 2 - 25 °C / 36 - 77 °F when not in use.
4. The product must be stored in proper places where only dental practitioners can access it.

VI. SHADES SYSTEM AND COMPONENTS

1. Shades

CLEARFIL AP-X PLT is available in 8 shades separately ; 6 shades correspond on a best-match basis with VITA Shade Guide plus 2 shades for special needs.

Select an appropriate shade by matching the shade guide to the restoration site.

The XL and HO shades offer a strong, lighter white color. The HO shade is translucent white and is suitable for sites requiring extra white color.

Standard shades:	A2, A3, A3.5, B2, B3
Cervical shades:	A4
Translucent shades:	XL, HO

2. Components

Please see the outside of the package for contents and quantity.

3. Ingredients

Silanated barium glass filler (70-90%)/
Bisphenol A diglycidylmethacrylate (Bis-GMA) (5-15%)/
Silanated colloidal silica (1-7%)/
Triethyleneglycol dimethacrylate (TEGDMA) (< 5%)/
Silanated silica filler (< 3%)/Benzoyl peroxide (< 1%)/
dl-Camphorquinone (< 1%)

Units in parentheses are mass %.

The total amount of inorganic filler is approx. 68 vol%

The particle size of inorganic fillers ranges from 0,02 to 17 µm.

VII. CLINICAL PROCEDURES

1. Cleaning tooth structure

Be sure the cavity is adequately cleaned. An adequately cleaned cavity assures maximum adhesive performance.

2. Moisture control

Avoid contamination of the treatment area from saliva or blood to produce optimal results. A rubber dam is recommended to keep the tooth clean and dry.

3. Cavity preparations

Remove any infected dentin and prepare the cavity in the usual manner.

4. Shade selection

Select an appropriate shade using the VITA Lumin-VACUUM shade guide.

5. Pulp protection

Any actual or near pulp exposure should be covered with a hard setting calcium hydroxide material. However, usually a cement lining or basing is not necessary. Do not use eugenol materials for pulp protection.

6. Acid etching uncut enamel

If the resin restorative material will extend to uncut enamel, apply etching agent (e.g. K-ETCHANT GEL) to the enamel, let it stay for 10 seconds, wash with water, and then dry.

[CAUTION]

The use of a bonding agent alone does not adequately condition uncut enamel. Overfilling of resin onto unetched, uncut enamel could cause marginal discoloration.

7. Application of bonding agent

Tooth surface treatment and bonding should be performed according to Instructions for Use of the bonding system used (e.g. CLEARFIL SE BOND or CLEARFIL TRI-S BOND).

8. Placement and light curing of CLEARFIL AP-X PLT

1) Dispensing

Place a PLT tip of the product into the dispenser barrel according to the Instructions for Use for dispenser. The tip may be rotated to provide the proper angle for delivery into the cavity.

[CAUTION]

Squeeze the dispenser with a slow, steady pressure. Excessive force is not necessary. Discard the tip after use and sterilize the dispenser according to its' Instructions for Use.

2) Placement

Incremental placement and light-curing each increment is strongly recommended especially in deep cavities and in Class II cavities.

3) Curing

Light-cure the resin with a dental visible light curing unit* following the table. Hold the light tip as close to the resin as possible.

Table: Relation between curing time and depth of cure for each dental curing unit

Type	Curing time (sec)	Depth of cure (mm)	
		XL, A2, A3, A3.5, B2, B3	A4, HO
Conventional halogen*	20	2.0	1.5
	40	2.0	2.0
Fast halogen*	5	2.0	1.5
	10	2.0	2.0
Plasma arc*	5	2.0	1.5
	10	2.5	2.0
LED*	20	2.0	1.5
	40	2.0	2.0

*Dental curing unit

Type	Light source	Wavelength range and light intensity
Conventional halogen	Halogen lamp	Light intensity ²⁾ of 150 – 550 mW/cm ² in wavelength range from 400 – 515 nm
Fast halogen	Halogen lamp	Light intensity ²⁾ of more than 550 mW/cm ² in wavelength range from 400 – 515 nm
Plasma arc	Xenon lamp	Light intensity ³⁾ of more than 2000 mW/cm ² in wavelength range from 400 – 515 nm, and light intensity of more than 450 mW/cm ² in wavelength range from 400 – 430 nm
LED	Blue LED ¹⁾	Light intensity ²⁾ of more than 300 mW/cm ² in wavelength range from 400 – 515 nm

1) Peak of emission spectrum: 450 – 480 nm

2) Evaluated according to ISO 10650-1.

3) Wavelength distribution and light intensity values measured with a spectro-radiometer calibrated using an IEC or the NIST (National Institute of Standards and Technology) standard lamp

9. Finishing

Contour the restoration and adjust the occlusion and using a fine diamond point.
Polish with silicon rubber points or polishing discs.

[WARRANTY]

Kuraray Noritake Dental Inc. will replace any product that is proved to be defective. Kuraray Noritake Dental Inc. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

[NOTE]


If a serious incident attributable to this product occurs, report it to the manufacturer's authorized representative shown below and the regulatory authorities of the country in which the user/patient resides.

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"CLEARFIL", "CLEARFIL AP-X", "CLEARFIL TRI-S BOND" and "SE BOND" are registered trademarks or trademarks of KURARAY CO., LTD.
"VITA" is a trademark of VITA Zahnfabrik, H.Rauter GmbH & Co. KG.
„Kuraray Noritake Dental Inc.“ pakeičia visus produktus su įrodomais trūkumais.
„VITA“ yra registruotas VITA Zahnfabrik, H. Rauter GmbH & Co. KG. prekių ženklai.