DUAL-CURED DENTAL BONDING AGENT

CLEARFIL™ PHOTO BOND

2°C/36°F 8°C/46°F

**ENGLISH INSTRUCTIONS FOR USE****I. INTRODUCTION**

CLEARFIL PHOTO BOND is a dual curing dentin and enamel bonding agent developed for use with composite light-curing restorative materials. The general clinical benefit of CLEARFIL PHOTO BOND is to restore tooth function in the following INDICATIONS FOR USE.

II. INDICATIONS FOR USE

Enamel and dentin bonding agent for direct composite restorations

III. CONTRAINDICATIONS

Patients with a history of hypersensitivity to methacrylate monomers

IV. POSSIBLE SIDE EFFECTS

The oral mucosal membrane may turn whitish when contacted by BOND due to the coagulation of protein. This is usually a temporary phenomenon that will disappear in a few days. Instruct patients to avoid irritating the affected area while brushing.

V. INCOMPATIBILITIES

- [1] Do not use eugenol-containing materials for pulp protection or temporary sealing since the eugenol can retard the curing process.
- [2] Do not use hemostatics containing ferric compounds since these materials may impair adhesion and may cause discoloration of the tooth margin or surrounding gingiva due to remaining ferric ions.
- [3] Do not use a hydrogen peroxide solution for cleaning cavities since it may weaken the bond strength to the tooth structure.

VI. PRECAUTIONS**1. Safety precautions**

- This product contains substances that may cause allergic reactions. Avoid use of the product in patients with known allergies to methacrylate monomers or any other components.
- If the patient demonstrates a hypersensitivity reaction, such as rash, eczema, features of inflammation, ulcer, swelling, itching or numbness, discontinue use of the product and seek medical attention.
- Use caution to prevent the product from coming in contact with the skin or getting into the eye. Before using the product, cover the patient's eyes with a towel to protect them in the event of splashing material.
- If the product comes in contact with human body tissues, take the following actions:
-If the product gets in the eye> Immediately wash the eye with copious amounts of water and consult a physician.
-If the product comes in contact with the skin or the oral mucosa> Immediately wipe the area with a cotton pellet or a gauze moistened with alcohol, and rinse with copious amounts of water.
- Use caution to prevent the patient from accidentally swallowing the product.
- Avoid using the same product dispensed into a well of the dispensing dish and the same disposable brush tip and the same sponge pledget for different patients to prevent cross-contamination infection. Discard the disposable brush tip and the sponge pledget after use.
- Wear gloves or take other appropriate protective measures to prevent the occurrence of hypersensitivity that may result from contact with methacrylate monomers or any other components.
- Dispose of this product as a medical waste to prevent infection.

2. Handling and manipulation precautions**[Common precautions]**

- The product must not be used for any purposes other than specified in [II.INDICATIONS FOR USE]
- The use of this product is restricted to dental professionals.
- Use a pulp capping agent in a cavity close to the pulp or in the event of accidental pulp exposure.
- To prevent poor performance or poor handling characteristics, observe the specified light-curing times and other handling requirements.
- Clean the cavity sufficiently to prevent poor bonding. If the adherend surface is contaminated with saliva or blood, wash it thoroughly and dry before bonding.
- When using hemostatics containing aluminum chloride, minimize its quantity, and use caution to prevent contact with the adherend surface. Failure to do so might weaken the bond strength to the tooth structure.
- Do not mix the product with any other dental materials.
- Use a rubber dam to prevent contamination and to control moisture.

[BONDING AGENT]

- Use the light blocking plate to avoid exposing the material to an operating light or natural light, and use within 2 minutes after dispensing. As the volatile ethanol which is contained in the BOND evaporates, the viscosity increases, thereby making it difficult to apply.
- Universal Liquid contains a flammable substance. Do not use it near an open flame.

[Dental light-curing unit]

- Do not look directly at the light source. Protective glasses are recommended.
- Low light intensity causes poor adhesion. Check the lamp for service life and the dental curing light guide tip for contamination. It is advisable to check the dental curing light intensity using an appropriate light evaluating device at appropriate intervals.
- The emitting tip of the dental curing unit should be held as near and vertical to the resin surface as possible. If a large resin surface is to be light-cured, it is advisable to divide the area into several sections and light-cure each section separately.

3. Storage precautions

- The product must be used by the expiration date indicated on the package.
- This product and BONDING AGENT must be refrigerated (2-8°C / 36-46°F) when not in use.
- The container of BONDING AGENT should be tightly capped immediately after use to reduce the evaporation of the volatile solvent (ethanol contained in Universal Liquid).
- The product must be kept away from extreme heat, direct sunlight or a flame.
- The product must be stored in a proper place where only dental practitioners can access.

VII. COMPONENTS

- Components
See outer package for the items included in the package.
1) BONDING AGENT
(1) Catalyst Liquid
(2) Universal Liquid
- Principal ingredients
1) BONDING AGENT
(1) Catalyst Liquid
• Bisphenol A diglycidylmethacrylate (Bis-GMA) (20-40%)
• 10-Methacryloyloxydecyl dihydrogen phosphate (MDP) (20-40%)
• 2-Hydroxyethyl methacrylate (HEMA) (5-25%)
• Hydrophobic aliphatic dimethacrylate (15-35%)
• Benzoyl peroxide (< 2.5%)
• di-Camphorquinone (< 1%)
(2) Universal Liquid
• Ethanol (> 90%)
• Catalysts (< 5%)
• Accelerators (< 3%)

Units in parentheses are mass %.

VIII. CLINICAL PROCEDURES**1. Isolation and Moisture control**

Avoid contamination of the treatment area from saliva or blood to produce optimal results. A rubber dam is recommended to keep the tooth clean and dry.

2. Cavity preparations

Remove any infected dentin and prepare the cavity in the usual manner.

3. Pulp protection

Any actual or near pulp exposure could be covered with a hard setting calcium hydroxide material or other indicated materials. Do not use eugenol-containing materials for pulp protection.

4. Acid etching

Dispense one drop of an etching agent(e.g. K-ETCHANT GEL) into a well of mixing dish and apply it to the entire cavity surface with a small brush. Leave in place for 40 seconds, wash with water, and then dry. [NOTE]

Care should be taken to prevent saliva from coming in contact with the etched tooth surface. Any contaminated tooth surface must be re-etched.

5. Application of BOND

5-1. Dispense equal amounts of Catalyst Liquid and Universal Liquid into a well of a mixing dish and mix immediately with a sponge pledget. [NOTE]

Use the mixture of Catalyst Liquid and Universal Liquid as soon as possible after mixing. The mixture must be used within 2 minutes after mixing.

5-2. Apply mixed BOND to the entire cavity surface with the sponge pledget. [NOTE]

Use caution not to allow saliva or exudates to contact the treated surfaces.

5-3. Dry the entire cavity surface sufficiently by blowing mild air until BOND does not move. Use a vacuum aspirator to prevent BOND from scattering.

5-4. When a composite resin restorative has a thickness of greater than 2 mm, BOND must be light-cured before the composite resin is placed into the cavity. Light-cure BOND for 10 seconds with a dental curing unit (see table "Dental curing unit") [NOTE]

If the thickness of composite resin is less than 2 mm, BOND does not need to be light-cured prior to placement of the composite resin.

Table: Dental curing unit

Type	Light source	Wavelength range and light intensity
Halogen	Halogen lamp	Light intensity ³⁾ of more than 300 mW/cm ² in wavelength range from 400 - 515 nm
LED	Blue LED ¹⁾	Light intensity ³⁾ of more than 300 mW/cm ² in wavelength range from 400 - 515 nm

1) Peak of emission spectrum: 450 - 480 nm
2) Evaluated according to ISO 10650-1.
3) Evaluated according to ISO 10650-2.

6. Place composite resin restorative

Apply composite resins (e.g. CLEARFIL MAJESTY Esthetic, CLEARFIL MAJESTY Posterior, CLEARFIL MAJESTY Flow) into the cavity, light-cure, finish and polish according to the manufacturer's instructions. [NOTE]

If BOND has been light-cured, composite resins must be placed into the cavity within 1 minute after applying BOND.

[WARRANTY]

Kuraray Noritake Dental Inc. will replace any product that is proven to be defective. Kuraray Noritake Dental Inc. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

[NOTE]

- If a serious incident attributable to this product occurs, report it to the manufacturer through the EU Importer, and to the regulatory authorities of the country in which the user/patient resides.
- The Summary of Safety and Clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed>. Basic UDI-DI: 458211076Bonding0233W

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"CLEARFIL" and "CLEARFIL MAJESTY" are registered trademarks or trademarks of KURARAY CO., LTD.