CLEARFIL™ PHOTO CORE PLT

EN	ENGLISH	English	LIGHT-CURE CORE BUILD-UP COMPOSITE MATERIAL
FR	FRANÇAIS	French	MATÉRIAU COMPOSITE DE RECONSTITUTION CORONAIRE PHOTOPOLYMERISABLE
ES	ESPAÑOL	Spanish	MATERIAL COMPOSITE DE RESTAURACIÓN DE MUÑONES FOTOPOLIMERIZABLE
IT	ITALIANO	Italian	MATERIALE COMPOSITO FOTOPOLIMERIZZABILE PER RICOSTRUZIONE
NL	NEDERLANDS	Dutch	LICHTUITHARDEND COMPOSITET STOMPBOUW MATERIAAL
DE	DEUTSCH	German	Lichthärtendes Stumpfaufbau-Komposit
SV	SVENSKA	Swedish	LJUSHÄRDANDE PELARUPPBÄGGNAD KOMPOSITMATERIAL
NO	NORSK	Norwegian	LYSHERDET KOMPOSITTMATERIALE TIL KJERNEOPPBÆGGING
FI	SUOMI	Finnish	VALOKOVETTEINEN YHDISTELMÄMUOVI PILARIMATERIAALI
DA	DANSK	Danish	LYSHÆRDENDE KERNOPBÆGNINGS-KOMPOSITMATERIALE
PT	PORTUGUÊS	Portuguese	MATERIAL COMPÓSITO FOTOPOLIMERIZÁVEL PARA CONSTRUÇÃO DE NÚCLEOS
EL	ΕΛΛΗΝΙΚΑ	Greek	ΦΩΤΟΠΟΛΥΜΕΡΙΖΟΜΕΝΟ ΥΛΙΚΟ ΑΠΟΚΑΤΑΣΤΑΣΗΣ ΚΟΛΟΒΩΜΑΤΩΝ ΣΥΝΘΕΤΗΣ ΡΗΤΙΝΗΣ
TR	TÜRKÇE	Turkish	IŞIKLA SERTLEŞEN KOR YAPIMI KOMPOZİT MATERYALİ

Kuraray Noritake Dental Inc.

1621 Sakazu, Kurashiki, Okayama 710-0801, Japan

EMERGO EUROPE

Westervoortsedijk 60,
6827 AT Arnhem, The Netherlands


Kuraray Europe GmbH (EU Importer)

Philipp-Reis-Str. 4,
65795 Hattersheim am Main, Germany
Phone: +49 (0)69 305 35835 Fax: +49 (0)69 305 98 35835
URL: <https://www.kuraraynoritake.eu>

1561R518R

CLEARFIL™ PHOTO CORE PLT

I. INTRODUCTION

CLEARFIL PHOTO CORE PLT is a light-cure, radiopaque core build-up composite material with a special dispensing system for quick placement. The general clinical benefit of this product is to restore tooth function in the following INDICATIONS FOR USE.

II. INDICATIONS FOR USE

Core build-up of vital or non-vital tooth

III. CONTRAINDICATION

Patients with a history of hypersensitivity to methacrylate monomers

IV. INCOMPATIBILITY

Do not use eugenol containing materials for pulp protection or temporary sealing, since the eugenol could retard the bonding system curing process.

V. PRECAUTIONS

1. Safety precautions

- This product contains substances that may cause allergic reactions. Avoid use of the product in patients with known allergies to methacrylate monomers or any other components.
- If any hypersensitivity occurs, such as dermatitis, discontinue use of the product and consult a physician.
- Wear gloves or take other appropriate protective measures to prevent the occurrence of hypersensitivity that may result from contact with methacrylate monomers.
- Use caution to prevent the product from coming in contact with the skin or getting into the eye. Before using the product, cover the patient's eyes with a towel to protect the patient's eyes from splashing material.
- If the product comes in contact with human body tissues, take the following actions:
 - <If the product gets in the eye> Immediately wash the eye with copious amounts of water and consult a physician.
 - <If the product comes in contact with the skin> Immediately wipe it off with a cotton pledget moistened with alcohol or gauze and wash with copious amounts of water.
- Use caution to prevent the patient from accidentally swallowing the product.
- Avoid looking directly at the dental curing unit when curing the product.
- For infection control reasons, PLT tips are for single use only.
- Dispose of this product as a medical waste to prevent infection.

2. Handling and manipulation precautions

- Do not use the product in conjunction with the other composite resins. Mixing materials may cause a change in physical properties, possibly a decrease, from the properties expected.
- When light curing the product, note the depth of cure in this instructions for use.
- The product polymerize under an operating light or natural light (sunlight from windows). Avoid operating light or natural light, and use within 5 minutes from dispensing.
- The emitting tip of the dental curing unit should be held as near and vertical to the resin surface as possible. If a large resin surface is to be light cured, it is advisable to divide the area into several sections and light-cure each section separately.
- Low intensity of light causes poor adhesion. Check the lamp for service life and the dental curing unit guide tip for contamination. It is advisable to check the dental curing light intensity using an appropriate light evaluating device at appropriate intervals.
- Squeeze the dispenser with a slow, steady pressure. Excessive force is not necessary.
- The use of this product is restricted to dental professionals.

3. Storage precautions

- The product must be used by the expiration date indicated on the package.
- Keep away from extreme heat or direct sunlight.
- The product must be stored at 10- 25°C /50-77°F when not in use.
- Do not refrigerate or store at low temperature. If the product is refrigerated or stored at a low temperature, it will become excessively soft.
- The product must be stored in proper places where only dental practitioners can access it.

VI. COMPONENTS

1) Components

See outer package for the items included in the package.

2) Ingredients

Principal ingredients:

- Silanated silica filler (30-50%)
- Silanated barium glass filler (20-60%)
- Bisphenol A diglycidylmethacrylate (Bis-GMA) (< 15%)
- Triethyleneglycol dimethacrylate (TEGDMA) (< 6%)
- Benzoyl peroxide (< 1%)
- dl-Camphorquinone (< 0.2%)

Units in parentheses are mass %.

The total amount of inorganic filler is approx. 68vol%.

The particle size of inorganic fillers ranges from 0.49 to 75 µm.

VII. CLINICAL PROCEDURES

A. CORE BUILD-UP OF VITAL TOOTH

A-1. Cleaning tooth structure

Be sure the cavity is adequately cleaned. An adequately cleaned cavity assures maximum adhesive performance.

A-2. Moisture control

Avoid contamination of the treatment area from saliva or blood to produce optimal results. A rubber dam is recommended to keep the tooth clean and dry.

A-3. Cavity preparations

Remove any infected dentin and prepare the cavity in the usual manner.

A-4. Pulp protection

Any actual or near pulp exposure could be covered with a hard setting calcium hydroxide material. There is no need for cement lining or basing. Do not use eugenol materials for pulp protection.

A-5. Applying a matrix strip and wedges

Apply a matrix strip and wedges as needed.

A-6. Tooth surface treatment and bonding

Tooth surface treatment and bonding should be performed according to the Instructions for Use of the bonding system used (e.g. CLEARFIL TRI-S BOND or CLEARFIL SE BOND).

A-7. Placement and light curing of CLEARFIL PHOTO CORE PLT

1. Dispensing

Place a PLT tip of the product into the dispenser barrel according to the Instructions for Use of the dispenser. The tip may be rotated to provide the proper angle for delivery into the cavity.

[CAUTION]

Squeeze the dispenser with a slow, steady pressure. Excessive force is not necessary.

Discard the tip after use and sterilize the dispenser according to the dispenser instructions for use. The product polymerize under an operating light or natural light (sunlight from windows). Avoid operating light or natural light, and use within 5 minutes from dispensing.

2. Placement

Incremental placement and light-curing each increment is strongly recommended especially in deep cavities.

3. Curing

Cure the resin with a dental curing unit. Hold the emitting tip as close to the resin as possible. The light curing time and depth of cure for each dental curing unit is used is shown below.

Light curing time and depth of cure

	Dental curing unit	Light curing time			
		5 sec.	10 sec.	20 sec.	40 sec.
depth of cure	Conventional halogen*	—	—	3.5mm	4.0mm
	Fast halogen*	4.0mm	4.5mm	—	—
	Plasma arc*	5.0mm	6.0mm	—	—
	LED*	—	—	4.0mm	4.5mm

*Dental curing unit

Type	Light source	Wavelength range and light intensity
Conventional halogen	Halogen lamp	Light intensity ²⁾ of 300 - 550 mW/cm ² in wavelength range from 400 - 515 nm
Fast halogen	Halogen lamp	Light intensity ²⁾ of more than 550 mW/cm ² in wavelength range from 400 - 515 nm
Plasma arc	Xenon lamp	Light intensity ³⁾ of more than 2000 mW/cm ² in wavelength range from 400 - 515 nm, and light intensity of more than 450 mW/cm ² in wavelength range from 400 - 430 nm.
LED	Blue LED ¹⁾	Light intensity ²⁾ of more than 300 mW/cm ² in wavelength range from 400 - 515 nm

1) Peak of emission spectrum: 450 - 480 nm

2) Evaluated according to ISO 10650-1.

3) Wavelength distribution and light intensity values measured with a spectro-radiometer calibrated using an IEC or the NIST (National Institute of Standards and Technology) standard lamp

A-8. Preparation of the abutment tooth

After curing the composite resin, prepare the abutment tooth.

B. CORE BUILD-UP OF NON-VITAL TOOTH

B-1. Cleaning tooth structure

Be sure the cavity is adequately cleaned. An adequately cleaned cavity assures maximum adhesive performance.

B-2. Moisture control

Avoid contamination of the treatment area from saliva or blood to produce optimal results. A rubber dam is recommended to keep the tooth clean and dry.

B-3. Cleaning cavity and preparing root canal

Remove the temporary sealing material from the cavity and filling material from the root canal. Using a Pizo reamer, prepare and clean the root canal opening.

B-4. Try-in of the post

A post of the proper diameter is fitted to the prepared root canal and the length of the post adjusted.

Sand blast the post as needed.

B-5. Either B-5a or B-5b

B-5a. When PANAVIA F 2.0 is used.

Cement the post into the root canal according to the Instructions for Use of PANAVIA F 2.0.

[NOTE]

Do not use CLEARFIL SE BOND when PANAVIA F 2.0 is used. Since PANAVIA F 2.0 plays the role as adhesive bonding agent, apply excess PANAVIA F 2.0 paste to the remaining tooth structure on which ED PRIMER II, a component of PANAVIA F 2.0 is applied and the surface of the post.

B-5b. When another resin luting cement is used.

Cement the post into the root canal according to the Instructions for Use of the resin luting cement.

B-6. Placement and light curing of CLEARFIL PHOTO CORE PLT

1. Dispensing

Place a PLT tip of the product into the dispenser barrel according to the Instructions for Use of the dispenser. The tip may be rotated to provide the proper angle for delivery into the cavity.

[CAUTION]

Squeeze the dispenser with a slow, steady pressure. Excessive force is not necessary.

Discard the tip after use and sterilize the dispenser according to the instructions for use of the dispenser.

The product polymerize under an operating light or natural light (sunlight from windows). Avoid operating light or natural light, and use within 5 minutes from dispensing.

2. Build-up

Using appropriate instruments, build-up the resin to the appropriate anatomy being careful to avoid incorporation of air bubbles.

Incremental placement and light-curing each increment is strongly recommended.

3. Curing

Refer to "Light curing time and depth of cure" in A-7.

B-7. Preparation an abutment tooth

After curing the composite resin, prepare the abutment tooth.

[WARRANTY]

Kuraray Noritake Dental Inc. will replace any product that is proved to be defective. Kuraray Noritake Dental Inc. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.


[NOTE]

- If a serious incident attributable to this product occurs, report it to the manufacturer through the EU Importer, and to the regulatory authorities of the country in which the user/patient resides.
- The Summary of Safety and Clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed>.
Basic UDI-DI: 458211076Core1013E

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"CLEARFIL", "CLEARFIL PHOTO CORE", "SE BOND", "CLEARFIL TRI-S BOND" and "PANAVIA" are registered trademarks or trademarks of KURARAY CO., LTD.