CLEARFIL™ SE BOND 2 [Bottle]

EN	ENGLISH	English	Dental Universal Self-Etch Adhesive
FR	FRANÇAIS	French	Adhésif dentaire universel auto-mordant
ES	ESPAÑOL	Spanish	Adhesivo autograbante universal de uso dental
IT	ITALIANO	Italian	Adesivo dentale universale automordenzante (self-etch)
NL	NEDERLANDS	Dutch	Universeel dentaal zelfetsend adhesief
DE	DEUTSCH	German	Dentales selbstätzendes Universal Adhäsiv
SV	SVENSKA	Swedish	Dental universell självetsande adhesiv
NO	NORSK	Norwegian	Dentalt universal selvetsend (Self-Etch)-adhesiv
FI	SUOMI	Finnish	Hampaan itse-etsaava yleissidosaine
DA	DANSK	Danish	Dentalt universalt selvætsende adhæsiv
PT	PORTUGUÊS	Portuguese	Adesivo dentário auto-condicionante universal
EL	ΕΛΛΗΝΙΚΑ	Greek	Οδοντικό Αuto-αδροποιούμενο Συγκολλητικό Γενικής Χρήσεως
TR	TÜRKÇE	Turkish	Dental Üniversal Kendinden Asitli Adeziv
PL	POLSKI	Polish	Uniwersalny dentystyczny samowytrawiający materiał adhezyjny
RO	ROMÂNĂ	Romanian	Adeziv dentar universal autogravant
HR	HRVATSKI	Croatian	Dentalni univerzalni samojetkajući adheziv
HU	MAGYAR	Hungarian	Fogászati univerzális önmarató ragasztó
SL	SLOVENŠČINA	Slovenian	Dentalno samojedkalno univerzalno lepilo
CS	ČESKY	Czech	Dentální univerzální samoleptací adhezivum
SK	SLOVENSKY	Slovak	Dentálne samoleptacie univerzálne adhezívum
BG	БЪЛГАРСКИ	Bulgarian	Стоматологичен универсален самоецващ адхезив
UK	УКРАЇНСЬКА	Ukrainian	БОНДИНГОВА СИСТЕМА
ET	EESTI	Estonian	Dentaalne universaalne isesöövivut adhesiiv
LV	LATVIEŠU	Latvian	Zobārstniecības universālais paškodinošais adhezīvs
LT	LIETUVIŠKAI	Lithuanian	Odontologinis universalus savaime ésdinantis adhezyvas

Kuraray Noritake Dental Inc.

1621 Sakazu, Kurashiki, Okayama 710-0801, Japan

EMERGO EUROPE

Westervoortsedijk 60,
6827 AT Arnhem, The Netherlands

Kuraray Europe GmbH (EU Importer)

Philipp-Reis-Str. 4,
65795 Hattersheim am Main, Germany


Phone:+49 (0)69 305 35835 Fax:+49 (0)69 305 98 35835

URL:<https://www.kuraraynoritake.eu>

1562R052R

CLEARFIL™ SE BOND 2 [Bottle]

I. INTRODUCTION

"CLEARFIL SE BOND 2" [Bottle] is a two component, self-etch, light-cured bonding agent, which is intended for universal use for both direct and indirect restorations. The PRIMER allows simultaneous treatment of both dentin and enamel using one liquid. The product is activated by a dual-cure mechanism when the BOND is mixed with CLEARFIL DC Activator. This allows it to be used with dual-cure or self-cure composite filling materials, cements, or core build-up materials. The general clinical benefit of this product is to restore tooth function in the following INDICATIONS FOR USE.

II. INDICATIONS FOR USE

"CLEARFIL SE BOND 2" [Bottle] is indicated for the following uses:

- [1] Direct restorations using light-cured composite resin
- [2] Cavity sealing as a pretreatment for indirect restorations
- [3] Treatment of exposed root surfaces
- [4] Treatment of hypersensitive teeth
- [5] Intraoral repairs of fractured restorations
- [6] Post cementation using a dual- or self-cured composite resin
- [7] Core build-ups using a light-, dual- or self-cured core material
- [8] Cementing inlays, onlays, crowns, bridges and veneers using a composite resin cement

III. CONTRAINDICATIONS

Patients with a history of hypersensitivity to methacrylate monomers

IV. POSSIBLE SIDE EFFECTS

The oral mucosal membrane may turn whitish when contacted by the product due to the coagulation of protein. This is usually a temporary phenomenon that will disappear in a few days. Instruct patients to avoid irritating the affected area while brushing.

V. INCOMPATIBILITIES

- [1] Do not use eugenol-containing materials for pulp protection or temporary sealing, since the eugenol can retard the curing process.
- [2] Do not use hemostatics containing ferric compounds since these materials may impair adhesion and may cause discoloration of the tooth margin or surrounding gingiva due to remaining ferric ions.
- [3] When using hemostatics containing aluminum chloride, minimize the quantity; use caution to prevent contact with the adherend surface. Failure to do so might weaken the bond strength to the tooth structure.

VI. PRECAUTIONS

1. Safety precautions

1. This product contains substances that may cause allergic reactions. Avoid use of the product in patients with known allergies to methacrylate monomers or any other components.
2. If the patient demonstrates a hypersensitivity reaction, such as rash, eczema, features of inflammation, ulcer, swelling, itching or numbness, discontinue use of the product and seek medical attention.
3. Use caution to prevent the product from coming in contact with the skin or getting into the eye. Before using the product, cover the patient's eyes with a towel to protect them in the event of splashing material.
4. If the product comes in contact with human body tissues, take the following actions:
 - <If the product gets in the eye>
Immediately wash the eye with copious amounts of water and consult a physician.
 - <If the product comes in contact with the skin or the oral mucosa>
Immediately wipe the area with a cotton pellet or a gauze pad moistened with alcohol, and rinse with copious amounts of water.
5. Use caution to prevent the patient from accidentally swallowing the product.
6. Avoid looking directly at the dental curing light when curing the product.
7. Avoid using the same product dispensed into a well of the dispensing dish and the same applicator brush for different patients to prevent cross contamination. The applicator brush is single-use only. Discard the applicator brush after use.
8. Wear gloves or take other appropriate protective measures to prevent the occurrence of hypersensitivity that may result from contact with methacrylate monomers or any other components.
9. Dispose of this product as a medical waste to prevent infection.

2. Handling and manipulation precautions

1. The product must not be used for any purposes other than specified in [II.INDICATIONS FOR USE].
2. The use of this product is restricted to dental professionals.
3. Use a pulp capping agent in a cavity close to the pulp or in the event of accidental pulp exposure.
4. To prevent poor performance or poor handling characteristics, observe the specified light-curing times and other handling requirements.
5. Clean the cavity sufficiently to prevent poor bonding. If the adherend surface is contaminated with saliva or blood, wash it thoroughly and dry before bonding.
6. Use the light blocking plate to avoid exposing the material to an operating light or natural light, and use within 3 minutes after dispensing. PRIMER and BOND will set to a gel if left under an operating light or natural light (sunlight from windows).
7. PRIMER is applied to the adherend surface and should be left in place for 20 seconds. Move the operating light away from the mouth or turn off the light to prevent the applied PRIMER and BOND from being exposed to the operating light (the BOND may gelate when left under the light).
8. Dry sufficiently with mild air after treatment with PRIMER for 20 seconds.
9. Do not rinse after applying PRIMER.
10. Do not use BOND in conjunction with other bonding agents except CLEARFIL DC Activator.
11. Do not use PRIMER in conjunction with other bonding agents and silane coupling agents except CLEARFIL PORCELAIN BOND ACTIVATOR.
12. When using with CLEARFIL DC Activator, use the light blocking plate to avoid exposing the material to an operating light or natural light, and use within 90 seconds after dispensing.

13. CLEARFIL DC Activator contains ethanol. Dry the entire adherend surface sufficiently with mild air for more than 5 seconds until the mixture does not move. Adjust the air pressure according to the shape and size of the cavity and the prosthetic appliance. Use a vacuum aspirator to prevent the bond liquid from scattering.
14. If the treated surface is contaminated, wash it with water, dry, or clean with alcohol, and treat again with the product.

[Dental light-curing unit]

1. Low light intensity causes poor adhesion. Check the lamp for service life and the dental curing light guide tip for contamination. It is advisable to check the dental curing light intensity using an appropriate light evaluating device at appropriate intervals.
2. The emitting tip of the dental curing unit should be held as near and vertical to the resin surface as possible. If a large resin surface has to be light-cured, it is advisable to divide the area into several sections and light-cure each section separately.

[Outer case]

1. Do not hold the press points when removing the cap.
2. When dispensing PRIMER or BOND, hold the outer case and point it vertically downward.
3. Do not store the outer case with the bottom up.

3. Storage precautions

1. The product must be used by the expiration date indicated on the package.
2. The product must be refrigerated (2-8°C/ 36-46°F) when not in use, and should be brought to room temperature for more than 15 minutes before using. After taking out of the refrigerator, the product must be left standing until it comes to room temperature.
3. The container should be tightly capped immediately after use.
4. The product must be kept away from extreme heat or direct sunlight.
5. The product must be stored in a proper place where only dental practitioners can access.

VII. COMPONENTS

See outer package for the items included in the package.

1) PRIMER

Principal ingredients:

- 10-Methacryloyloxydecyl dihydrogen phosphate (MDP) (5-25%)
- 2-Hydroxyethyl methacrylate (HEMA) (20-40%)
- Hydrophilic aliphatic dimethacrylate (5-15%)
- Accelerators (< 6%)
- dl-Camphorquinone (< 1%)
- Water (40-60%)

2) BOND

Principal ingredients:

- 10-Methacryloyloxydecyl dihydrogen phosphate (MDP) (3-7%)
- Bisphenol A diglycidylmethacrylate (Bis-GMA) (25-45%)
- 2-Hydroxyethyl methacrylate (HEMA) (20-40%)
- Hydrophobic aliphatic dimethacrylate (10-30%)
- dl-Camphorquinone (< 1%)
- Phenyl bis(2,4,6-trimethylbenzoyl)-phosphine oxide (< 3%)
- Accelerators (< 3%)
- Colloidal silica (5-15%)

Units in parentheses are mass %.

3) Accessories

- Applicator brush (fine <silver>)
- Mixing dish*
- Light blocking plate*
- Outer case*
- *Consumables

VIII. CLINICAL PROCEDURES

A. Standard procedure I

- [1] Direct restorations using light-cured composite resin
- [2] Cavity sealing as a pretreatment for indirect restorations
- [3] Treatment of exposed root surfaces
- [4] Treatment of hypersensitive teeth

A-1. Isolation and Moisture control

Avoid contamination of the treatment area from saliva or blood to produce optimum results. A rubber dam is recommended to keep the tooth clean and dry.

A-2. Cavity preparations

Remove any infected dentin and prepare the cavity in the usual manner.

A-3. Pulp protection

Any actual or near pulp exposure can be covered with a hard setting calcium hydroxide material. There is no need for a cement liner or base. Do not use eugenol materials for pulp protection.

A-4. Acid-etching of enamel

If there is an adherend surface including uncut enamel, apply a phosphoric acid-etching gel (e.g. K-ETCHANT Syringe) to the uncut enamel, leave it in place for 10 seconds, then rinse and dry. If acid-etching to the cut enamel is preferred, it is optional to apply a phosphoric acid to the cut enamel in the same manner as described above.

A-5. Application of PRIMER

1. Dispense the necessary amount of PRIMER into a well of the dispensing dish immediately before application.

[CAUTION]

Use the light blocking plate to avoid exposing the material to an operating light or natural light, and use within 3 minutes after dispensing.

2. Apply PRIMER to the entire cavity wall with an applicator brush. Leave it in place for 20 seconds. Use caution not to allow saliva or exudate to contact the treated surfaces.

3. Dry the entire cavity wall sufficiently with mild air for more than 5 seconds until the PRIMER does not move. Use a vacuum aspirator to prevent the PRIMER from scattering.

[CAUTION]

Do not rinse after applying PRIMER.

A-6. Application of BOND

1. Dispense the necessary amount of BOND into a well of the dispensing dish immediately before application.

[CAUTION]

Use the light blocking plate to avoid exposing the material to an operating light or natural light, and use within 3 minutes after dispensing.

2. Apply BOND to the entire cavity wall with an applicator brush. Use caution not to allow saliva or exudate to contact the treated surfaces.
3. Make a uniform bond film using a gentle air flow.
4. Light-cure BOND with a dental curing unit (see table "Dental curing unit and curing time").

Table: Dental curing unit and curing time

Type	Light source	Light Intensity	Light-curing time
Halogen	Halogen lamp	More than 400 mW/cm ²	10 seconds
LED	Blue LED*	800 – 1400 mW/cm ²	10 seconds
		More than 1500 mW/cm ²	5 seconds

The effective wavelength range of each dental curing unit must be 400 - 515 nm.

*Peak of emission spectrum: 450 - 480 nm

A-7. Placement of composite resin restorative or Treatment of hypersensitive teeth

A-7a. Direct restorations using light-cured composite resin

Apply composite resins (e.g. CLEARFIL MAJESTY ES-2, CLEARFIL MAJESTY ES Flow) into the cavity, light-cure, finish and polish according to the manufacturer's instructions.

A-7b. Cavity sealing and treatment of exposed root surfaces

Apply a thin coat of composite resins (e.g. CLEARFIL MAJESTY ES Flow) onto the tooth, and light-cure according to the manufacturer's instructions. Remove unpolymerized resin with a cotton pellet or a gauze moistened with alcohol.

A-7c. Treatment of hypersensitive teeth

Remove the unpolymerized layer of BOND with a cotton pellet or a gauze moistened with alcohol.

B. Standard procedure II

[5] Intraoral repairs of fractured restorations

B-1. Preparation of adherend surfaces

Roughen the adherend surfaces using a diamond point, and place a bevel at the marginal area.

B-2. Phosphoric acid treatment of adherend surfaces

Apply a phosphoric acid (e.g. K-ETCHANT Syringe) to the adherend surface (including enamel if present). Leave in place for 5 seconds, then rinse and dry.

B-3. Precious metal surface

When precious metal is used, apply a metal primer (e.g. ALLOY PRIMER) according to the manufacturer's instructions.

B-4. Silane treatment

1. Dispense one drop each of PRIMER and CLEARFIL PORCELAIN BOND ACTIVATOR into a well of the dispensing dish, and mix.

[NOTE]

PRIMER and CLEARFIL PORCELAIN BOND ACTIVATOR should be mixed immediately before application.

[CAUTION]

Use the light blocking plate to avoid exposing the material to an operating light or natural light, and use within 3 minutes after dispensing.

2. Apply mixture to the adherend surface with an applicator brush.

3. Leave it in place for 5 seconds. Use caution not to allow saliva or exudate to contact the treated surfaces. If the adherend surface extends to the tooth, leave mixture in place for 20 seconds.

4. Dry the entire cavity wall sufficiently with mild air for more than 5 seconds until the mixture does not move. Use a vacuum aspirator to prevent the liquid from scattering.

[CAUTION]

Do not wash treated surface.

B-5. Bonding

Apply BOND to the entire adherend surface and light-cure. Refer to section A-6.

B-6. Place composite resin restorative

Apply composite resin (e.g. CLEARFIL MAJESTY ES-2) into the cavity, light-cure, finish and polish according to the manufacturer's instructions.

[NOTE]

Use an opaque resin (e.g. CLEARFIL ST OPAQUER) to mask metal color.

C. Standard procedure III

[6] Post cementation using a dual- or self-cured composite resin

C-1. Isolation and Moisture control

Avoid contamination of the treatment area from saliva or blood to produce optimal results. A rubber dam is recommended to keep the tooth clean and dry.

C-2. Preparing the root canal

Prepare and clean the root canal opening in the usual manner.

C-3. Post preparation

For glass fiber posts

1. Phosphoric acid treatment

Apply a phosphoric acid (e.g. K-ETCHANT Syringe) to the post surface. Leave it in place for 5 seconds, then rinse and dry.

2. Post surface treatment

Apply mixture of PRIMER and CLEARFIL PORCELAIN BOND ACTIVATOR. Refer to section B-4.

[CAUTION]

- Do not blast glass fiber posts with alumina powder, or the posts could be damaged.
- Avoid all contamination of the surfaces to be treated during pre-treatment and until the final core build-up.

For metal posts

1. Surface preparation

Roughen the post surface as necessary, according to the manufacturer's instructions.

2. Surface treatment

Apply PRIMER to the post surface and leave it in place for 5 seconds. Then dry the entire post surface thoroughly with mild air for more than 5 seconds, until the liquid does not move.

C-4. Application of PRIMER

1. Dispense the necessary amount of PRIMER into a well of the dispensing dish immediately before application.

[CAUTION]

Use the light blocking plate to avoid exposing the material to an operating light or

natural light, and use within 3 minutes after dispensing.

2. Apply PRIMER to the entire cavity wall with an applicator brush. Leave it in place for 20 seconds. Use caution not to allow saliva or exudate to contact the treated surfaces.
3. Dry the entire cavity wall sufficiently with mild air for more than 5 seconds until the PRIMER does not move. Use a vacuum aspirator to prevent the PRIMER from scattering. [CAUTION]

• Do not rinse after applying PRIMER.

• Remove excess PRIMER with a paper point inserted into the root canal during air-drying. After wiping away the excess PRIMER, dry the adherend surface again, if necessary.

C-5. Application of BOND

1. Dispense one drop each of BOND and CLEARFIL DC Activator into a well of the dispensing dish, and mix well. [CAUTION]

Use the light blocking plate to avoid exposing the material to an operating light or natural light, and use within 90 seconds after dispensing.

2. Apply the mixture to the entire cavity wall with the applicator brush. Use caution not to allow saliva or exudate to contact the treated surfaces.
3. Dry the entire adherend surface thoroughly with mild air for more than 5 seconds, until the bond does not move. Use a vacuum aspirator to prevent the bond liquid from scattering. Remove the excess mixture with a paper point. After wiping away the excess mixture, dry the adherend surface again, if necessary.

[NOTE]

To dry thoroughly, adjust the air pressure according to the shape and size of the adherend surface.

4. Light-cure BOND with a dental curing unit (see table "Dental curing unit and curing time" in A-6.)

[NOTE]

Light-curing the mixture is an option for optimum performance.

C-6. Post cementation

Place the post and cement it using a composite resin according to the manufacturer's instructions.

D. Standard procedure IV

[7] Core build-ups using light-, dual- or self-cured core material

D-1. Moisture control and preparing the root canal

Refer to C-1 and C-2.

D-2. Post cementation

Place the post using a luting material according to the manufacturer's instructions.

D-3. Application of PRIMER and BOND

When using with a light-cured core material or CLEARFIL DC CORE PLUS

Refer to the sections from A-5 to A-6.

When using with an other dual- or self-cured core material

Use with CLEARFIL DC Activator. Refer to sections from C-4 to C-5.

D-4. Core build-up

Perform a core build-up using a material (e.g. CLEARFIL PHOTO CORE, CLEARFIL DC CORE PLUS) according to the manufacturer's instructions.

E. Standard procedure V

[8] Cementing inlays, onlays, crowns, bridges and veneers using a composite resin cement

E-1. Conditioning the cavity and core (tooth, metal, composite) surfaces

- (1) Remove the temporary sealing material and temporary cement in the usual manner, and clean the cavity or core using moisture control.
- (2) Trial fit the prosthetic restoration to check its fit in the cavity or core. When using a try-in paste to check the color, follow the manufacturer's instructions.

E-2. Surface preparation of prosthetic restorations

Follow the Instructions for Use of the restorative material. In the absence of specific instructions, we recommend the following procedure:

For silica-based ceramics (e.g. conventional porcelain, lithium disilicate)

1. Etch glass ceramic surfaces with hydrofluoric acid in accordance with the manufacturer's instructions, or apply a phosphoric acid (e.g. K-ETCHANT Syringe), leave it in place for 5 seconds, then rinse and dry.

2. Silane treatment

Apply a mixture of PRIMER and CLEARFIL PORCELAIN BOND ACTIVATOR. Refer to section B-4.

For metal oxides (e.g. zirconia) or metals

1. Roughen the adherend surface by blasting with 30 to 50 µm alumina powder at an air pressure of 0.1-0.4MPa (1.4-5.8 PSI/ 1-4 kgf/cm²). The air pressure should be properly adjusted to suit the material and/ or shape of the prosthetic restoration, using caution to prevent chipping. After blasting, clean the prosthetic restoration by using ultrasound for 2 minutes followed by drying it with an air stream.

2. Apply PRIMER to the adherend surface, leave it for 5 seconds.

3. Dry the entire surface thoroughly with mild air for more than 5 seconds, until the liquid does not move.

[Note]

For precious metals, use a metal primer (e.g. ALLOY PRIMER) according to the manufacturer's instructions.

For composite resins

1. Roughen the adherend surface by blasting. Refer to the section E-2 "For metal oxides".
2. Apply a phosphoric acid (e.g. K-ETCHANT Syringe), leave it in place for 5 seconds, then rinse and dry.

3. Apply a mixture of PRIMER and CLEARFIL PORCELAIN BOND ACTIVATOR. Refer to section B-4.

[CAUTION]

Avoid all contamination of the surfaces to be treated during pre-treatment and until final cementation.

E-3. Acid-etching of enamel

If there is an adherend surface including uncut enamel, apply a phosphoric acid-etching gel (e.g. K-ETCHANT Syringe) to the uncut enamel, leave it in place for 10 seconds, then rinse and dry. If acid-etching to the cut enamel is preferred, it is optional to apply a phosphoric acid to the cut enamel in the same manner as described above.

[CAUTION]

When cementing adhesion bridges or veneers, apply a phosphoric acid (e.g. K-ETCHANTE Syringe) to the enamel surface.

E-4. Pretreatment of an abutment tooth

Apply PRIMER to the adherend surface. Refer to section C-4.

E-5. Bonding to an abutment tooth

Mix one drop of BOND and CLEARFIL DC Activator, then apply to the adherend surface. Refer to section C-5.

E-6. Cementing

Cement the prosthetic restoration using a resin cement according to the manufacturer's instructions.

[WARRANTY]

Kuraray Noritake Dental Inc. will replace any product that is proven to be defective. Kuraray Noritake Dental Inc. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.


[NOTE]

- If a serious incident attributable to this product occurs, report it to the manufacturer through the EU Importer, and to the regulatory authorities of the country in which the user/patient resides.
- The Summary of Safety and Clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed>. Basic UDI-DI: 458211076Bonding1954U

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"CLEARFIL", "SE BOND", "CLEARFIL MAJESTY", "CLEARFIL ST", "CLEARFIL DC CORE PLUS" and "CLEARFIL PHOTO CORE" are registered trademarks or trademarks of KURARAY CO., LTD.