CLEARFIL™ SE BOND 2 [Unit Dose]

EN	ENGLISH	English	Dental Universal Self-Etch Adhesive
FR	FRANÇAIS	French	Adhésif dentaire universel auto-mordantant
ES	ESPAÑOL	Spanish	Adhesivo autograbante universal de uso dental
IT	ITALIANO	Italian	Adesivo dentale universale automordenzante (self-etch)
NL	NEDERLANDS	Dutch	Universeel dentaal zelfetsend adhesief
DE	DEUTSCH	German	Dentales selbstätzendes Universal Adhäsiv
SV	SVENSKA	Swedish	Dental universell självetsande adhesiv
NO	NORSK	Norwegian	Dentalt universal selvetsend (Self-Etch)-adhesiv
FI	SUOMI	Finnish	Hampaan itse-etsaava yleissidosaine
DA	DANSK	Danish	Dentalt universalt selvætsende adhæsiv
PT	PORTUGUÊS	Portuguese	Adesivo dentário auto-condicionante universal
EL	ΕΛΛΗΝΙΚΑ	Greek	Οδοντικό Αuto-αδροποιούμενο Συγκολλητικό Γενικής Χρήσεως
TR	TÜRKÇE	Turkish	Dental Üniversal Kendinden Asitli Adeziv
PL	POLSKI	Polish	Uniwersalny dentystyczny samowytrawiający materiał adhezyjny
RO	ROMÂNĂ	Romanian	Adeziv dentar universal autogravant
HR	HRVATSKI	Croatian	Dentalni univerzalni samojetkajući adheziv
HU	MAGYAR	Hungarian	Fogászati univerzális önmarató ragasztó
SL	SLOVENŠČINA	Slovenian	Dentalno samojedkalno univerzalno lepilo
CS	ČESKY	Czech	Dentální univerzální samoleptací adhezivum
SK	SLOVENSKY	Slovak	Dentálne samoleptacie univerzálne adhezívum
BG	БЪЛГАРСКИ	Bulgarian	Стоматологичен универсален самоецващ адхезив
UK	УКРАЇНСЬКА	Ukrainian	Стоматологічний універсальний самопротравлювальний адгезив
ET	EESTI	Estonian	Dentaalne universaalne isesöövituv adhesiiv
LV	LATVIEŠU	Latvian	Zobārstniecības universālais paškodinošais adhezivs
LT	LIETUVIŠKAI	Lithuanian	Odontologinis universalus savaime išdinantis adhezyvas

Kuraray Noritake Dental Inc.

1621 Sakazu, Kurashiki, Okayama 710-0801, Japan

EMERGO EUROPE

Westervoortsedijk 60,
6827 AT Arnhem, The Netherlands


Kuraray Europe GmbH (EU Importer)

Philipp-Reis-Str. 4,
65795 Hattersheim am Main, Germany
Phone:+49 (0)69 305 35835 Fax:+49 (0)69 305 98 35835
URL:<https://www.kuraraynoritake.eu>

CLEARFIL™ SE BOND 2 [Unit Dose]

I. INTRODUCTION

CLEARFIL SE BOND 2 [Unit Dose] is a two component, self-etch, light-cured bonding agent. The PRIMER allows simultaneous treatment of dentin, enamel, non-precious metals and metal oxides using one liquid. The product is provided for single-use. The general clinical benefit of this product is to restore tooth function for the following INDICATIONS FOR USE.

II. INDICATIONS FOR USE

CLEARFIL SE BOND 2 [Unit Dose] is indicated for the following uses:

- [1] Direct restorations using light-cured composite resin
- [2] Cavity sealing as a pretreatment for indirect restorations
- [3] Treatment of exposed root surfaces
- [4] Treatment of hypersensitive teeth
- [5] Intraoral repairs of fractured restorations
- [6] Core build-ups using a light- or dual-cured core material
- [7] Cementing inlays, onlays, crowns, bridges and veneers using a composite resin cement

III. CONTRAINDICATIONS

Patients with a history of hypersensitivity to methacrylate monomers

IV. POSSIBLE SIDE EFFECTS

The oral mucosal membrane may turn whitish when contacted by the product due to the coagulation of protein. This is usually a temporary phenomenon that will disappear in a few days. Instruct patients to avoid irritating the affected area while brushing.

V. INCOMPATIBILITIES

- [1] Do not use eugenol-containing materials for pulp protection or temporary sealing, since the eugenol can retard the curing process.
- [2] Do not use hemostatics containing ferric compounds since these materials may impair adhesion and may cause discoloration of the tooth margin or surrounding gingiva due to remaining ferric ions.
- [3] When using hemostatics containing aluminum chloride, minimize the quantity; use caution to prevent contact with the adherend surface. Failure to do so might weaken the bond strength to the tooth structure.

VI. PRECAUTIONS

1. Safety precautions

1. This product contains substances that may cause allergic reactions. Avoid use of the product in patients with known allergies to methacrylate monomers or any other components.
2. If the patient demonstrates a hypersensitivity reaction, such as rash, eczema, features of inflammation, ulcer, swelling, itching or numbness, discontinue use of the product and seek medical attention.
3. Use caution to prevent the product from coming in contact with the skin or getting into the eye. Before using the product, cover the patient's eyes with a towel to protect them in the event of splashing material.
4. If the product comes in contact with human body tissues, take the following actions:
 - <If the product gets in the eye>
Immediately wash the eye with copious amounts of water and consult a physician.
 - <If the product comes in contact with the skin or the oral mucosa>
Immediately wipe the area with a cotton pellet or a gauze pad moistened with alcohol, and rinse with copious amounts of water.
5. Use caution to prevent the patient from accidentally swallowing the product.
6. Avoid looking directly at the dental curing light when curing the product.
7. Single use only. Avoid using the same applicator brush on different patients to prevent cross contamination. Discard the container and applicator brush after each use.
8. Wear gloves or take other appropriate protective measures to prevent the occurrence of hypersensitivity that may result from contact with methacrylate monomers or any other components.
9. Dispose of this product as a medical waste to prevent infection.

2. Handling and manipulation precautions

1. The product must not be used for any purposes other than specified in [II.INDICATIONS FOR USE].
2. The use of this product is restricted to dental professionals.
3. Use a pulp capping agent in a cavity close to the pulp or in the event of accidental pulp exposure.
4. To prevent poor performance or poor handling characteristics, observe the specified light-curing times and other handling requirements.
5. Clean the cavity sufficiently to prevent poor bonding. If the adherend surface is contaminated with saliva or blood, wash it thoroughly and dry before bonding.
6. PRIMER is applied to the adherend surface and should be left in place for 20 seconds. Move the operating light away from the mouth or turn off the light to prevent the applied PRIMER and BOND from being exposed to the operating light (the BOND may gelate when left under the light).
7. Dry sufficiently with mild air after treatment with PRIMER for 20 seconds.
8. Do not rinse after applying PRIMER.
9. If the treated surface is contaminated, wash it with water, dry, or clean with alcohol, and treat again with the product.
10. Apply PRIMER/BOND with the applicator brush immediately after snapping off the container cap.
11. Do not use the same applicator brush for PRIMER and BOND.

12. When snapping off the container cap, do not tilt the container to avoid spilling PRIMER/BOND.
13. Do not mix the product with other materials.

[Dental light-curing unit]

1. Low light intensity causes poor adhesion. Check the lamp for service life and the dental curing light guide tip for contamination. It is advisable to check the dental curing light intensity using an appropriate light evaluating device at appropriate intervals.
2. The emitting tip of the dental curing unit should be held as near and vertical to the resin surface as possible. If a large resin surface has to be light-cured, it is advisable to divide the area into several sections and light-cure each section separately.

3. Storage precautions

1. The product must be used by the expiration date indicated on the package.
2. The product must be refrigerated (2-8°C/ 36-46°F) when not in use, and should be brought to room temperature for more than 15 minutes before using. After taking out of the refrigerator, the product must be left standing until it comes to room temperature.
3. The product must be kept away from extreme heat or direct sunlight.
4. The product must be stored in a proper place where only dental practitioners can access.

VII. COMPONENTS

Please see the outside of the package for contents and quantity.

1) PRIMER (Yellow container)

Principal ingredients:

- 10-Methacryloyloxydecyl dihydrogen phosphate (MDP) (5-25%)
- 2-Hydroxyethyl methacrylate (HEMA) (20-40%)
- Hydrophilic aliphatic dimethacrylate (5-15%)
- Accelerators (< 6%)
- dl-Camphorquinone (< 1%)
- Water (40-60%)

2) BOND (Black container)

Principal ingredients:

- 10-Methacryloyloxydecyl dihydrogen phosphate (MDP) (3-7%)
- Bisphenol A diglycidylmethacrylate (Bis-GMA) (25-45%)
- 2-Hydroxyethyl methacrylate (HEMA) (20-40%)
- Hydrophobic aliphatic dimethacrylate (10-30%)
- dl-Camphorquinone (< 1%)
- Phenyl bis(2,4,6-trimethylbenzoyl)-phosphine oxide (< 3%)
- Accelerators (< 3%)
- Colloidal silica (5-15%)

Units in parentheses are mass %.

3) Accessory

- Applicator brush (fine <silver>)

VIII. CLINICAL PROCEDURES

A. Standard procedure I

- [1] Direct restorations using light-cured composite resin
- [2] Cavity sealing as a pretreatment for indirect restorations
- [3] Treatment of exposed root surfaces
- [4] Treatment of hypersensitive teeth

A-1. Isolation and Moisture control

Avoid contamination of the treatment area from saliva or blood to produce optimum results. A rubber dam is recommended to keep the tooth clean and dry.

A-2. Cavity preparations

Remove any infected dentin and prepare the cavity in the usual manner.

A-3. Pulp protection

Any actual or near pulp exposure can be covered with a hard setting calcium hydroxide material. There is no need for a cement liner or base. Do not use eugenol materials for pulp protection.

A-4. Acid-etching of enamel

If there is an adherend surface including uncut enamel, apply a phosphoric acid-etching gel (e.g. K-ETCHANT Syringe) to the uncut enamel, leave it in place for 10 seconds, then rinse and dry. If acid-etching to the cut enamel is preferred, it is optional to apply a phosphoric acid to the cut enamel in the same manner as described above.

A-5. Application of PRIMER

1. Snap-off the container cap and apply PRIMER to the entire cavity wall with the applicator brush; leave it in place for 20 seconds. Use caution not to allow saliva or exudate to contact the treated surfaces.
2. Dry the entire cavity wall sufficiently with mild air for more than 5 seconds until the PRIMER does not move. Use a vacuum aspirator to prevent the PRIMER from scattering.
 - [CAUTION]
 - Do not rinse after applying PRIMER.

A-6. Application of BOND

1. Snap-off the container cap and apply BOND to the entire cavity wall with an applicator brush. Use caution not to allow saliva or exudate to contact the treated surfaces.
2. Make a uniform bond film using a gentle air flow.
3. Light-cure BOND with a dental curing unit (see table "Dental curing unit and curing time").

Table: Dental curing unit and curing time

Type	Light source	Light Intensity	Light-curing time
Halogen	Halogen lamp	More than 400 mW/cm ²	10 seconds
LED	Blue LED*	800 – 1400 mW/cm ²	10 seconds
		More than 1500 mW/cm ²	5 seconds

The effective wavelength range of each dental curing unit must be 400 - 515 nm.

*Peak of emission spectrum: 450 - 480 nm

A-7. Placement of composite resin restorative or Treatment of hypersensitive teeth

A-7a. Direct restorations using light-cured composite resin

Apply composite resins (e.g. CLEARFIL MAJESTY ES-2, CLEARFIL MAJESTY ES Flow) into the cavity, light-cure, finish and polish according to the manufacturer's instructions.

A-7b. Cavity sealing and treatment of exposed root surfaces

Apply a thin coat of composite resins (e.g. CLEARFIL MAJESTY ES Flow) onto the tooth, and light-cure according to the manufacturer's instructions. Remove unpolymerized resin with a cotton pellet or a gauze moistened with alcohol.

A-7c. Treatment of hypersensitive teeth

Remove the unpolymerized layer of BOND with a cotton pellet or a gauze moistened with alcohol.

B. Standard procedure II

[5] Intraoral repairs of fractured restorations

B-1. Preparation of adherend surfaces

Roughen the adherend surfaces using a diamond point, and place a bevel at the marginal area.

B-2. Phosphoric acid treatment of adherend surfaces

Apply a phosphoric acid (e.g. K-ETCHANT Syringe) to the adherend surface (including enamel if present). Leave in place for 5 seconds, then rinse and dry.

B-3. Precious metal treatment

When the adherend surface includes precious metal, apply a metal primer (e.g. ALLOY PRIMER) according to the manufacturer's instructions.

B-4. Silane treatment

When the adherend surface includes silica-based ceramics (e.g. conventional porcelain, lithium disilicate) or composite resins, apply a silane coupling agent (e.g. CLEARFIL CERAMIC PRIMER), according to the manufacturer's instructions.

B-5. Application of PRIMER and BOND

1. Snap-off the container cap and apply PRIMER to the adherend surface and leave it in place for 5 seconds (20 seconds for tooth surface). Then dry the entire surface thoroughly with mild air for more than 5 seconds, until the liquid does not move. Use a vacuum aspirator to prevent the PRIMER from scattering. Use caution not to allow saliva or blood to contact the treated surfaces.

[CAUTION]

Do not rinse after applying PRIMER.

2. Apply BOND to the entire adherend surface and light-cure. Refer to section A-6.

B-6. Place composite resin restorative

Apply composite resin (e.g. CLEARFIL MAJESTY ES-2) into the cavity, light-cure, finish and polish according to the manufacturer's instructions.

[NOTE]

Use an opaque resin (e.g. CLEARFIL ST OPAQUER) to mask metal color.

C. Standard procedure III

[6] Core build-ups using a light- or dual-cured core material

C-1. Isolation and Moisture control

Avoid contamination of the treatment area from saliva or blood to produce optimal results. A rubber dam is recommended to keep the tooth clean and dry.

C-2. Preparing the root canal

Prepare and clean the root canal opening in the usual manner.

C-3. Post cementation

Place the post using a luting material according to the manufacturer's instructions.

C-4. Application of PRIMER and BOND

Apply PRIMER and BOND to the entire adherend surface. Refer to the sections from A-5 to A-6.

C-5. Core build-up

Perform a core build-up using a material (e.g. CLEARFIL PHOTO CORE, CLEARFIL DC CORE PLUS) according to the manufacturer's instructions.

[CAUTION]

When using a dual-cure core material, the first increment should be light-cured within the depth of cure of the material.

D. Standard procedure IV

[7] Cementing inlays, onlays, crowns, bridges and veneers using a composite resin cement

D-1. Conditioning the cavity and core (tooth, metal, composite) surfaces

1. Remove the temporary sealing material and temporary cement in the usual manner, and clean the cavity or core using moisture control.
2. Trial fit the prosthetic restoration to check its fit in the cavity or core. When using a try-in paste to check the color, follow the manufacturer's instructions.

D-2. Surface preparation of prosthetic restorations

Follow the Instructions for Use of the restorative material. In the absence of specific instructions, we recommend the following procedure:

For metal oxides (e.g. zirconia) or non-precious metals

1. Roughen the adherend surface by blasting with 30 to 50 µm alumina powder at an air pressure of 0.1-0.4MPa (14-58 PSI/ 1-4 kgf/cm²). The air pressure should be properly adjusted to suit the material and/ or shape of the prosthetic restoration, using caution to prevent chipping. After blasting, clean the prosthetic restoration by using ultrasound for 2 minutes followed by drying it with an air stream.
2. Snap-off the container cap and apply PRIMER to the adherend surface; leave it for 5 seconds.
3. Dry the entire surface thoroughly with mild air for more than 5 seconds, until the liquid does not move. Use a vacuum aspirator to prevent the PRIMER from scattering. Use caution not to allow saliva or blood to contact the treated surfaces.
[Note]
 - For precious metals, use a metal primer (e.g. ALLOY PRIMER) according to the manufacturer's instructions.
 - For silica-based ceramics/composite resins, Use a silane coupling agent (e.g. CLEARFIL CERAMIC PRIMER), according to the manufacturer's instructions.

D-3. Cementing

Cement the prosthetic restoration using a resin cement (e.g. CLEARFIL ESTHETIC CEMENT) according to the manufacturer's instructions.

[WARRANTY]

Kuraray Noritake Dental Inc. will replace any product that is proven to be defective. Kuraray Noritake Dental Inc. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

[NOTE]


If a serious incident attributable to this product occurs, report it to the manufacturer's authorized representative shown below and the regulatory authorities of the country in which the user/patient resides.

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"CLEARFIL", "SE BOND", "CLEARFIL MAJESTY", "CLEARFIL ST", "CLEARFIL DC CORE PLUS" and "CLEARFIL PHOTO CORE" are registered trademarks or trademarks of KURARAY CO., LTD.