CZR PRESS LF

Rx Only

EN	ENGLISH	English	Dental Ceramic
FR	FRANÇAIS	French	Céramique Dentaire
ES	ESPAÑOL	Spanish	Cerámica Dental
IT	ITALIANO	Italian	Ceramica Dentale
NL	NEDERLANDS	Dutch	DENTAAL PORSELEIN
DE	DEUTSCH	German	Dentalkeramik
SV	SVENSKA	Swedish	Dentalkeramik
NO	NORSK	Norwegian	Dentalkeramikk
FI	SUOMI	Finnish	Hammaskeramiikka
DA	DANSK	Danish	Dental Keramik
PT	PORTUGUÊS	Portuguese	Cerâmica Dentária
EL	ΕΛΛΗΝΙΚΑ	Greek	Οδοντικό Κεραμικό Υλικό
TR	TÜRKÇE	Turkish	Dental Seramik
PL	POLSKI	Polish	Dentystyczny materiał ceramiczny
RO	ROMÂNĂ	Romanian	Ceramică dentară
HR	HRVATSKI	Croatian	Dentalna keramika
HU	MAGYAR	Hungarian	Fogászati kerámia
SL	SLOVENŠČINA	Slovenian	Dentalna keramika
CS	ČESKY	Czech	Dentální keramika
SK	SLOVENSKY	Slovak	Dentálna keramika
BG	БЪЛГАРСКИ	Bulgarian	Стоматологична керамика
UK	УКРАЇНСЬКА	Ukrainian	Стоматологічна кераміка
ET	EESTI	Estonian	Dentaalkeraamika
LV	LATVIEŠU	Latvian	Zobārstniecības keramika
LT	LIETUVIŠKAI	Lithuanian	Odontologinė keramika

UK REP	EN	UK Responsible Person	PL	Osoba odpowiedzialna w Wielkiej Brytanii
	FR	Personne responsable au Royaume-Uni	RO	Persoana responsabilă din Marea Britanie
	ES	Responsable en el Reino Unido	HR	Osoba odgovorna za Ujedinjeno Kraljevstvo (UK)
	IT	Responsabile per il Regno Unito	HU	Felelős személy az Egyesült Királyságban
	NL	Verantwoordelijke in het VK	SL	Odgovorna oseba v Združenem kraljestvu
	DE	Für UK verantwortliche Person	CS	Odpovědná osoba ve Spojeném království
	SV	Ansvarig person i Storbritannien	SK	Zodpovedná osoba v Spojenom kráľovstve
	NO	Ansvarlig person i Storbritannia	BG	Отговорно лице за Обединеното кралство
	FI	Vastuushenkilö Yhdistyneessä kuningaskunnassa	UK	Відповідальна особа у Великій Британії
	DA	Ansvarlig person i Det Forenede Kongerige	ET	Suurbritannia vastutav isik
	PT	Responsável no Reino Unido	LV	Arvienotās Karalistes atbildīgā persona
	EL	Υπεύθυνο Πρόσωπο για το Ηνωμένο Βασίλειο	LT	Atsakingas asmuo Jungtinėje Karalystėje
	TR	Birleşik Krallık Sorumlusu		

Kuraray Noritake Dental Inc.

300 Higashiyama, Miyoshi-cho, Miyoshi, Aichi 470-0293, Japan
TEL +81-561-32-8953 FAX +81-561-32-8976

EMERGO EUROPE

Westervoortsedijk 60, 6827 AT Arnhem, The Netherlands

Kuraray Europe GmbH (EU Importer)

Philipp-Reis-Str. 4,
65795 Hattersheim am Main, Germany
Phone:+49 (0)69 305 35835 Fax:+49 (0)69 305 98 35835
URL:https://www.kuraraynoritake.eu

CZR PRESS LF

I. Introduction

CZR PRESS LF is a specific type of porcelain for dental restorations. The system consists of glass powder and liquid to be applied utilizing the layering technique and baked onto a zirconia (ZrO₂) ceramic frame. The general clinical benefit of this product is to restore tooth function for the following INDICATIONS FOR USE.

II. Indications for Use

For use in prosthetic dentistry to create an all-ceramic prosthesis.
-Ceramic for coverage of a zirconia substructure

III. Components

For a complete list of the components (including shades and accessories) for this product, see "CZR PRESS LF Baking Schedule".

IV. Composition

Potassium-aluminosilicate glass (>70 %), inorganic pigments (0-30 %)
LF LIQUID: Water (>99 %), oxidizing agents (<1 %)
IS LIQUID: 2-Phenoxyethanol (>99 %)
ES LIQUID: Propane-1,2-diol (>99 %)

Units in parentheses are mass %.

V. Type and Class (ISO6872: 2015)

Type:1 / Class:1

VI. Physical Properties

See Table 2.

VII. Directions for Use

See Table 1 for the Baking schedule.

See "Comprehensive Instructions for Zirconia Solutions" for detailed information about this product's use, materials available, and color variations.

1. Preparation of zirconia framework

- Adjust the thickness of the zirconia framework in the usual manner, assuring that all areas for use on the anterior teeth to the canines are at least 0.4 mm thick, and those for use on premolars and molars are at least 0.5 mm thick.
- When using in combination with CZR PRESS Ingot, the margins should be finished with a perfectly smooth edge, maintaining a minimum depth of 0.8 mm so that press ingot material can be applied.

2. Checking for cracks and cleaning

- Check the zirconia framework for cracks.
- After making sure the framework does not have any cracks, sandblast the surface with aluminum oxide (50µm) at a pressure of 0.2MPa or less, followed by about 10 minutes of ultrasonic cleaning.
- If any cracks are found, do not use that zirconia framework.

3. 1st application and baking of Shade Base Stain

- Mix the desired shade of CZR PRESS LF Shade Base Stain with IS LIQUID only.
- The viscosity of the mixture should have a consistency like "Maple Syrup". Apply a thin and even layer of the mixture, covering the zirconia framework to a thickness of 0.15 mm.
- Bake it according to Table 1.

4. 2nd application and baking of Shade Base Stain

- As in the case of the 1st application and baking, apply the desired shade of CZR PRESS LF Shade Base Stain as necessary before baking.

5. After baking the Shade Base Stain, wax-up the dentinal form or full-contour crown and then press the ingot material to the zirconia framework following the Instructions for Use supplied with CZR PRESS.

6. If CERABIEN™ ZR is used, apply the porcelain to form the dentin contour, followed by baking and defining the shape, according to the Instructions for Use supplied with CERABIEN™ ZR.

7. If the dentin structure requires characterization staining, apply LF Internal Stain and bake it according to Table 1.

8. Application of LF Enamel

- Mix LF Enamel with LF LIQUID.
- Apply the mixture to the incisal one-third of the dentinal form of the crown.

9. Application of LF Translucent and/or LF Luster

- Apply LF Translucent and/or LF Luster to cover the whole surface of the crown, extending it slightly beyond the target contour of the crown to allow for the possibility of shrinkage during baking.

10. Baking the tooth component layered porcelains

- Bake the porcelains applied to the restoration according to Table 1.
- After baking, the surface of a correctly baked restoration should have a slight luster or sheen.

11. Adjusting for morphological corrections, glazing, and completion

- Correct and contour the morphological shape of the restoration, followed by cleaning and glazing. If necessary, apply LF External Stain or LF External Stain Glaze and bake it according to Table 1.

12. Additional application of LF Add-On

- Any shortage or contact additions can be corrected by adding LF Add-On porcelain after glazing. Bake the LF Add-On porcelain according to Table 1. After baking, finish the restoration by polishing in the usual manner.

VIII. Remarks on Handling

Contraindications:

If the patient is hypersensitive to potassium-aluminosilicate glass or any other components, this product must not be used.

Caution and Warning:

- If the patient or the dental professional demonstrates a hypersensitivity reaction, such as rash, dermatitis, etc., discontinue use of the product and seek medical attention immediately.
- When mixing materials, or grinding a restoration being fabricated, use an approved dust mask and vacuum with an air filter to protect your lungs from inhaling either dust or liquid.

- When mixing materials or grinding a restoration being fabricated, use safety glasses to prevent dust or liquid from getting into your eyes. If dust or liquid gets into your eyes, immediately rinse with copious amounts of water and seek medical attention.
- Do not use for any purposes except for fabricating dental restorations. This product is for dental application only.
- The restoration is very hot immediately after baking. Do not touch it with your bare hands for a considerable amount of time after removing it from the furnace.
- IS LIQUID and ES LIQUID are flammable. Do not use them near fire or high-temperature objects.
- IS LIQUID is a stimulant liquid. When using it, have good ventilation and if necessary, wear a protective mask, safety glasses and gloves.
- Stains and plaque can accumulate on a restoration after it has been placed in the patient's oral cavity, depending on the patient's eating habits or food choices. Instruct patients about the importance of cleaning their restorations.

Precautions in use:

- This product should only be used on frameworks made of KATANA™ Zirconia or zirconia (3YTZP type) in which a small amount of yttria formulation is utilized.
- This product's coefficient of thermal expansion is different from that of the porcelain materials used for dental metal ceramic restorations and aluminum oxide porcelain materials (including Stain). Do not mix this product with any other dental porcelain materials or use it in unauthorized combinations.
- The correct baking temperature for this product varies with different furnaces. Perform a test baking trial in advance to identify the appropriate baking temperature for this product.
- The correct baking temperature for this product varies with the size of the restorations to be baked. Identify the appropriate baking temperature by referring to the baking schedule covered in the Technical Instructions.
- The Stain porcelain should be mixed using the specified liquid. (Use of liquids other than those specified may result in a darker restoration color than desired, after baking.)
- For staining and glazing the crown surface, use LF External Stain. (Do not use LF Internal Stain for this purpose.)
- This product is a low-temperature baking porcelain. To prevent discoloration of the porcelain, observe the following conditions:
 - Only use Magic Separator when separating the porcelain from the gypsum die.
 - If tissue paper fibers get incorporated into the porcelain slurry during water sorption in the condensation procedure, the fibers will not burn off completely. After drying, check to make sure that no residual tissue fibers remain before baking.
 - Always use freshly mixed porcelain. Do not use other liquids than those specified for mixing with this product.
 - Porcelain furnace temperature variations may be significant in the lower ranges. Determine the exact baking schedule by making a preliminary test before baking an actual case restoration.
- Keep the interior of the porcelain furnace clean by periodically firing it above 1,000 °C (1,832°F) while empty.
- Dispose of this product as a medical waste to prevent infection.

Storage:

- The product must be stored in a cool and dry place. Keep away from direct sunlight.
- The product should be stored at 1-30°C (33.8-86°F).
- The product must be stored in an appropriate place where only dental personnel have access.
- The product must be used by the expiration date indicated on the package.
- Do not store IS LIQUID and ES LIQUID near fire or high-temperature objects.

[WARRANTY]

Kuraray Noritake Dental Inc. will replace any product that is proven to be defective. Kuraray Noritake Dental Inc. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

[NOTE]

- If a serious incident attributable to this product occurs, report it to the manufacturer's authorized representative shown below and the regulatory authorities of the country in which the user/patient resides.
- Refill does not include instructions for use (IFU). We provide the IFU in electronic format on the dedicated eIFU website.
- "CERABIEN" and "KATANA" are registered trademarks or trademarks of NORITAKE CO., LIMITED.

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