EN	ENGLISH	English	Dental Ceramic
FR	FRANÇAIS	French	Céramique Dentaire
ES	ESPAÑOL	Spanish	Cerámica Dental
IT	ITALIANO	Italian	Ceramica Dentale
NL	NEDERLANDS	Dutch	DENTAAL PORSELEIN
DE	DEUTSCH	German	Dentalkeramik
SV	SVENSKA	Swedish	Dentalkeramik
NO	NORSK	Norwegian	Dentalkeramikk
FI	SUOMI	Finnish	Hammaskeramiikka
DA	DANSK	Danish	Dental Keramik
PT	PORTUGUÊS	Portuguese	Cerâmica Dentária
EL	ΕΛΛΗΝΙΚΑ	Greek	Οδοντικό Κεραμικό Υλικό
TR	TÜRKÇE	Turkish	Dental Seramik
PL	POLSKI	Polish	Dentystyczny materiał ceramiczny
RO	ROMÂNĂ	Romanian	Ceramică dentară
HR	HRVATSKI	Croatian	Dentalna keramika
HU	MAGYAR	Hungarian	Fogászati kerámia
SL	SLOVENŠČINA	Slovenian	Dentalna keramika
CS	ČESKY	Czech	Dentální keramika
SK	SLOVENSKY	Slovak	Dentálna keramika
BG	БЪЛГАРСКИ	Bulgarian	Стоматологична керамика
UK	УКРАЇНСЬКА	Ukrainian	Стоматологічна кераміка
ET	EESTI	Estonian	Dentaalkeraamika
LV	LATVIEŠU	Latvian	Zobārstniecības keramika
LT	LIETUVIŠKAI	Lithuanian	Odontologinė keramika

MD	EN	Medical Device	DA	Medicinsk udstyr	CS	Zdravotnický prostředek
	FR	Dispositif médical	PT	Dispositivo médico	SK	Zdravotnícka pomôcka
	ES	Producto sanitario	EL	Ιατροτεχνολογικό προϊόν	BG	Медицинско изделие
	IT	Dispositivo medico	TR	Tıbbi cihaz	UK	Виріб медичного призначення
	NL	Medisch hulpmiddel	PL	Wyrób medyczny	ET	Meditsiiniseade
	DE	Medizinprodukt	RO	Dispozitiv medical	LV	Medicīniska ierīce
	SV	Medicinteknisk produkt	HR	Medicinski proizvod	LT	Medicinos priemonė
	NO	Medisinsk utstyr	HU	Orvostechikai eszköz		
	FI	Lääkinnällinen laite	SL	Medicinski pripomoček		

Kuraray Noritake Dental Inc.

300 Higashiyama, Miyoshi-cho, Miyoshi, Aichi 470-0293, Japan
TEL +81-561-32-8953 FAX +81-561-32-8976

EMERGO EUROPE

Prinsessegracht 20
2514 AP The Hague, The Netherlands

Kuraray Europe GmbH

Philipp-Reis-Str. 4,
65795 Hattersheim am Main, Germany
Phone: +49 (0)69 305 35835 Fax: +49 (0)69 305 98 35835
URL: <https://www.kuraraynoritake.eu>

CZR PRESS

CE 1639

I. Introduction

CZR PRESS is a specific type of porcelain for dental restorations. The system consists of glass ingots to be applied utilizing the pressing technique onto a zirconia (ZrO₂) ceramic frame.

II. Indications for Use

- For use in prosthetic dentistry to create an all-ceramic prosthesis.
- Ceramic for coverage of a zirconia substructure

III. Components

For the components (including shades and accessories) of this product, see Table 1.

IV. Composition

Potassium-aluminosilicate glass, inorganic pigment, etc.

V. Type and Class (ISO6872: 2015)

Type:II / Class:1

VI. Physical Properties

See Table 2.

VII. Directions for Use

See "Comprehensive Instructions for Zirconia Solutions" for detailed information about this product's use, materials available, and color variations.

1. Preparation of zirconia framework

- Adjust the thickness of the zirconia framework in the usual manner, assuring that all areas for use on the anterior teeth to the canines are at least 0.4 mm thick, and those for use on premolars and molars are at least 0.5 mm thick.
 - The margins should be finished with a perfectly smooth edge, maintaining a minimum depth of 0.8 mm so that press ingot material can be applied.
- ### 2. Checking for cracks and cleaning
- Check the zirconia framework for cracks.
 - After making sure the framework does not have any cracks, sandblast the surface with aluminum oxide (50µm) at a pressure of 0.2MPa or less, followed by about 10 minutes of ultrasonic cleaning.
 - If any cracks are found, do not use that zirconia framework.

3. 1st application and baking of Shade Base Stain

- Mix the desired shade of CZR PRESS LF Shade Base Stain with IS LIQUID only.
- The viscosity of the mixture should have a consistency like "Maple Syrup". Apply a thin and even layer of the mixture, covering the zirconia framework to a thickness of 0.15 mm.
- Bake it according to the Instructions for Use supplied with CZR PRESS LF.

4. 2nd application and baking of Shade Base Stain

- As in the case of the 1st application and baking, apply the desired shade of CZR PRESS LF Shade Base Stain as necessary before baking.

5. Wax-up

- Place the zirconia framework in position on the abutment model and begin wax up of the restoration to the framework, starting at the margins.
- Avoid sharp angles or deep undercuts to prevent poor pressing. After waxing up, remove any excess wax from the inside of the framework. Make sure that there is no clearance or gaps between wax and framework.

[Stain method]

Wax-up the framework to full-contour form of the restoration needed.

[Layering method]

Wax-up the framework to form the complete shape of the dentin. (The mamelon structures will be defined after pressing ingot material)

6. Preparing for investment

- Cut a 3.3 mm diameter wax wire sprue to a length of 2 to 3 mm and attach it to the ring former.
- Attach the waxed restoration to the other end of the wax wire sprue utilizing the outer incline of the nonfunctional cusps of the wax pattern form, or to the marginal ridges if the wax pattern form is a molar.
- The angle for attaching the sprue wax line should be within the range of 30 to 60°
- The distance from the wax pattern to the inner wall of the ring must be not less than 8 mm and the distance to the top leveling cap of the ring must be not less than 10 mm. If multiple wax patterns are attached to the ring former, they must be spaced not less than 5 mm apart from each other, with their marginal lines perpendicular to the sprue line and all at the same height.
- It is advisable to apply petrolatum (petroleum jelly) in advance of investment to the ring, ring former, and ring gauge (leveling cap).

7. Investment and burn-out

Using one of the investments recommended in "Comprehensive Instructions for Zirconia Solutions", invest and burn-out the wax pattern in the specified manner.

8. Selecting of the press ingot

- Select H-ingot when using the stain method.
- Select L-ingot when using the layering method
- Select EW-ingot when matching the bleach EW shades. (Esthetic White)

9. Pressing

- Insert ingots and a plunger vertically into the investment ring. Place the investment ring in position inside the press furnace.
- Perform the pressing as directed in "Comprehensive Instructions for Zirconia Solutions".
- After pressing, remove the investment ring from the press furnace and let it slowly cool down to room temperature.

10. Divesting

- Locate and mark the position of the plunger in the investment ring. Using a cutting disc, cut a circumferential notch into the investment at the height at which the bottom of the plunger is located. Separate the investment ring using plaster nippers.
- Then, remove the bulk of the investment material from the restoration(s) by sandblasting with aluminum oxide (50µm) at a pressure of 0.4 to 0.6MPa.

- Once the pressed restoration is exposed to view, sandblast it using glass beads at a pressure of 0.2MPa along the long axis of the crown.

11. Cutting off the sprue

Cut off the sprue at a distance of more than 2 mm from the pressed restoration by using a diamond disc. Then, use a diamond point to remove the remainder of the sprue from the pressed restoration.

12. Adjust for morphological corrections

Put the pressed restoration on the working model and check the fit at the margins.

[Stain method]

Shape the pressed restoration to the contours you want to achieve checking contacts and occlusion.

[Layering method]

Form and define the mamelon structure using a diamond disc and point, taking care to ensure that the final thickness of the pressed restoration remains at least 0.8 mm. Sandblast the pressed restoration with aluminum oxide (50µm) at a pressure of 0.2MPa, followed by about 10 minutes of ultrasonic cleaning.

13. Finishing using the stain method

- Characterize by applying CZR PRESS LF External Stain to the surface of the crown restoration, followed by baking at the specified parameters.
- Perform the 1st glazing application. Mix CZR PRESS LF External Stain Glaze with IS LIQUID. The viscosity of the mixture should have a consistency like "Maple Syrup".
- Apply the mixture in a thin and even coat, covering the entire surface of the crown restoration to a thickness of 0.2 mm. Bake it according to the Instructions for Use supplied with CZR PRESS LF.
- After baking, adjust the contact area and shape of the restoration again, if necessary.
- Clean the crown restoration for about 10 minutes, using an ultrasonic cleaner.
- Mix CZR PRESS LF External Stain Glaze with IS LIQUID, using a slightly lower viscosity than that for the 1st glazing step. Apply the mixture to the entire surface of the crown restoration, followed by baking at the specified parameters. If there is insufficient gloss after baking, apply an additional layer of the mixture, followed by re-baking.

14. Finishing using the layering method

- Apply CZR PRESS LF onto the surface of the restoration and bake it according to the Instructions for Use supplied with CZR PRESS LF.
- If an internal stain is needed to add characterization, use CZR PRESS LF Internal Stain.
- It is difficult to develop a sufficient level of gloss by self-glazing if the restoration is made from an ingot. Apply CZR PRESS LF Translucent, covering the entire surface of the restoration or apply CZR PRESS LF External Stain Glaze.
- Correct the morphological shape in the usual manner.
- If surface stain characterization is needed, apply CZR PRESS LF External Stain as instructed in "Comprehensive Instructions for Zirconia Solutions".

VIII. Remarks on Handling

Contraindications:

If the patient is hypersensitive to potassium-aluminosilicate glass or any other components, this product must not be used.

Caution and Warning:

- If the patient or the dental professional demonstrates a hypersensitivity reaction, such as rash, dermatitis, etc., discontinue use of the product and seek medical attention immediately.
- When mixing materials, or grinding a restoration being fabricated, use an approved dust mask and vacuum with an air filter to protect your lungs from inhaling either dust or liquid.
- When mixing materials or grinding a restoration being fabricated, use safety glasses to prevent dust or liquid from getting into your eyes. If dust or liquid gets into your eyes, immediately rinse with copious amounts of water and seek medical attention.
- Do not use for any purposes except for fabricating dental restorations. This product is for dental application only.
- The restoration is very hot immediately after baking. Do not touch it with your bare hands for a considerable amount of time after removing it from the furnace.
- Stains and plaque can accumulate on a restoration, once it has been put in position in the patient's oral cavity depending on the patient's eating habits or food choices. Instruct patients about the importance of cleaning their restorations.
- Never use a cracked zirconia framework.

Precautions in use:

- This product should only be used on frameworks made of zirconia (3YTZP type) in which a small amount of yttria formulation is utilized.
- This product's coefficient of thermal expansion is different from that of the porcelain materials used for dental metal ceramic restorations and aluminum oxide porcelain materials (including Stain). Do not mix this product with any other dental porcelain materials or use it in unauthorized combinations.
- The correct press temperature for this product varies with different press furnaces. Perform a test baking trial in advance to identify the appropriate press temperature for this product.
- Use frameworks with the specified thicknesses. The margins of the frameworks should be finished with a perfectly smooth edge.
- Do not use dental porcelain materials other than those specified.
- Ingots must not be reused.
- The wax used for investment must be kept clean.
- Dispose of this product as a medical waste to prevent infection.

Storage:

- The product must be stored in a cool and dry place. Keep away from direct sunlight.
- The product should be stored at 1-30°C (33.8-86°F).
- The product must be stored in an appropriate place where only dental personnel have access.
- The product must be used by the expiration date indicated on the package.

[WARRANTY]

Kuraray Noritake Dental Inc. will replace any product that is proven to be defective. Kuraray Noritake Dental Inc. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

[NOTE]

1. If a serious accident attributable to this product occurs, report it to the manufacturer's authorized representative shown below and the regulatory authorities of the country in which the user/patient resides.
2. Refill does not include instructions for use (IFU). We provide the IFU in electronic format on the dedicated eIFU website.

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