|           |             |            |  |
|-----------|-------------|------------|--|
| <b>EN</b> | ENGLISH     | English    | Dental Zirconia                            |
| <b>FR</b> | FRANÇAIS    | French     | Zircone dentaire                           |
| <b>ES</b> | ESPAÑOL     | Spanish    | Óxido de zirconio dental                   |
| <b>IT</b> | ITALIANO    | Italian    | Ossido di zirconio dentale                 |
| <b>NL</b> | NEDERLANDS  | Dutch      | Dentaaal zirkoniumoxide                    |
| <b>DE</b> | DEUTSCH     | German     | Dentales Zirkonoxid                        |
| <b>SV</b> | SVENSKA     | Swedish    | Dental zirkoniumdioxid                     |
| <b>NO</b> | NORSK       | Norwegian  | Dental zirkoniumdioksid                    |
| <b>FI</b> | SUOMI       | Finnish    | Hammaslääketieteellinen zirkoniumdioksidi  |
| <b>DA</b> | DANSK       | Danish     | Dentalt zirkonia                           |
| <b>PT</b> | PORTUGUÊS   | Portuguese | Óxido de zircónio dentário                 |
| <b>EL</b> | ΕΛΛΗΝΙΚΑ    | Greek      | Οδοντιατρικό οξειδίο του ζirkονίου         |
| <b>TR</b> | TÜRKÇE      | Turkish    | Dental zirkonyum oksit                     |
| <b>PL</b> | POLSKI      | Polish     | Tlenek cyrkonu do aplikacji dentystycznych |
| <b>RO</b> | ROMÂNĂ      | Romanian   | Oxid de zirconiu dentar                    |
| <b>HR</b> | HRVATSKI    | Croatian   | Dentalni cirkonij                          |
| <b>HU</b> | MAGYAR      | Hungarian  | Fogászati cirkonoxid                       |
| <b>SL</b> | SLOVENŠČINA | Slovenian  | Dentalni cirkonij                          |
| <b>CS</b> | ČESKY       | Czech      | Dentální zirkonoxid                        |
| <b>SK</b> | SLOVENSKY   | Slovak     | Dentálny oxid zirkónu                      |
| <b>BG</b> | БЪЛГАРСКИ   | Bulgarian  | Стоматологичен цирконий                    |
| <b>UK</b> | УКРАЇНСЬКА  | Ukrainian  | Стоматологічний цирконій                   |
| <b>ET</b> | EESTI       | Estonian   | Dentaaalne tsirkoonium                     |
| <b>LV</b> | LATVIEŠU    | Latvian    | Zobārstniecības cirkonijs                  |
| <b>LT</b> | LIETUVIŠKAI | Lithuanian | Odontologinis cirkonio oksidas             |

|               |           |   |           |  |
|---------------|-----------|---|-----------|--|
| <b>UK REP</b> | <b>EN</b> | UK Responsible Person                         | <b>PL</b> | Osoba odpowiedzialna w Wielkiej Brytanii       |
|               | <b>FR</b> | Personne responsable au Royaume-Uni           | <b>RO</b> | Persoana responsabilă din Marea Britanie       |
|               | <b>ES</b> | Responsable en el Reino Unido                 | <b>HR</b> | Osoba odgovorna za Ujedinjeno Kraljevstvo (UK) |
|               | <b>IT</b> | Responsabile per il Regno Unito               | <b>HU</b> | Felelős személy az Egyesült Királyságban       |
|               | <b>NL</b> | Verantwoordelijke in het VK                   | <b>SL</b> | Odgovorna oseba v Združenem kraljestvu         |
|               | <b>DE</b> | Für UK verantwortliche Person                 | <b>CS</b> | Odpovědná osoba ve Spojeném království         |
|               | <b>SV</b> | Ansvarig person i Storbritannien              | <b>SK</b> | Zodpovedná osoba v Spojenom kráľovstve         |
|               | <b>NO</b> | Ansvarlig person i Storbritannia              | <b>BG</b> | Отговорно лице за Обединеното кралство         |
|               | <b>FI</b> | Vastuushenkilö Yhdistyneessä kuningaskunnassa | <b>UK</b> | Відповідальна особа у Великій Британії         |
|               | <b>DA</b> | Ansvarlig person i Det Forenede Kongerige     | <b>ET</b> | Suurbritannia vastutav isik                    |
|               | <b>PT</b> | Responsável no Reino Unido                    | <b>LV</b> | Apvienotās Karalistes atbildīgā persona        |
|               | <b>EL</b> | Υπεύθυνο Πρόσωπο για το Ηνωμένο Βασίλειο      | <b>LT</b> | Atsakingas asmuo Jungtinėje Karalystėje        |
|               | <b>TR</b> | Birleşik Krallık Sorumlusu                    |           |  |

**Kuraray Noritake Dental Inc.**

300 Higashiyama, Miyoshi-cho, Miyoshi, Aichi 470-0293, Japan  
TEL +81-561-32-8953 FAX +81-561-32-8976

**Kuraray Europe GmbH (EU Importer)**

Philipp-Reis-Str. 4,  
65795 Hattersheim am Main, Germany  
Phone:+49 (0)69 305 35835 Fax:+49 (0)69 305 98 35835  
URL:<https://www.kuraraynoritake.eu>


**EMERGO EUROPE**

Westervoortsedijk 60, 6827 AT Arnhem, The Netherlands

## HTML PLUS (High Translucent Multi Layered PLUS)

### I. Introduction

This IFU is for KATANA Zirconia HTML PLUS (High Translucent Multi Layered PLUS). KATANA Zirconia HTML PLUS is a pre-sintered zirconia disc 98.5 mm in diameter, which contains a plastic ring. This is designed for all milling systems using this generic-type disc. (Please refer to your milling system's technical instructions for correct machine operation.) KATANA Zirconia HTML PLUS has 3 available thicknesses: (14 mm, 18 mm and 22 mm) and has 14 shade variations: (A1, A2, A3, A3.5, A4, B1, B2, B3, C1, C2, C3, D2, D3, NW). KATANA Zirconia HTML PLUS consists of 4 graded shade layers. KATANA Zirconia HTML PLUS is recommended for use in fabricating FCZ (Full Contour Zirconia) restorations or the frameworks. The general clinical benefit of this product is to restore tooth function in the following INTENDED USE/INDICATIONS FOR USE.

### II. INTENDED USE/ INDICATIONS FOR USE

KATANA Zirconia is used for the fabrication of the all-ceramic restorations (frameworks, FCZ crowns, FCZ bridges, inlays, onlays and veneers).

[Intended user]

- Dentists (chairside fabrication of restorations; clinical workflow)
- Dental laboratory technicians (fabrication of restorations in the dental laboratory)

### III. Sintering program

#### Sintering program 1

| Temperature                  | Programming Rate        | Holding Time |
|------------------------------|-------------------------|--------------|
| Room Temp. — 1550°C (2822°F) | 10°C/min. (18°F/min.)   | —            |
| 1550°C (2822°F)              | —                       | 2 hrs        |
| 1550°C (2822°F) — Room Temp. | -10°C/min. (-18°F/min.) | —            |

#### Sintering program 2

| Temperature                       | Programming Rate        | Holding Time |
|-----------------------------------|-------------------------|--------------|
| Room Temp. — 1400°C (2552°F)      | 50°C/min. (90°F/min.)   | —            |
| 1400°C (2552°F) — 1500°C (2732°F) | 4°C/min. (7°F/min.)     | —            |
| 1500°C (2732°F) — 1560°C (2840°F) | 10°C/min. (18°F/min.)   | —            |
| 1560°C (2840°F)                   | —                       | 16 min       |
| 1560°C (2840°F) — Room Temp. (*)  | -50°C/min. (-90°F/min.) | —            |

\* The restorations may be removed at 800°C (1472°F) or less depending on the circumstances.

#### Sintering program 3

| Temperature                       | Programming Rate          | Holding Time |
|-----------------------------------|---------------------------|--------------|
| Room Temp. — 1450°C (2642°F)      | 120°C/min. (216°F/min.)   | —            |
| 1450°C (2642°F) — 1600°C (2912°F) | 10°C/min. (18°F/min.)     | —            |
| 1600°C (2912°F)                   | —                         | 20 min       |
| 1600°C (2912°F) — Room Temp. (*)  | -120°C/min. (-216°F/min.) | —            |

\* The restorations may be removed at 800°C (1472°F) or less depending on the circumstances.

### IV. Composition

ZrO<sub>2</sub> (80-95%), Y<sub>2</sub>O<sub>3</sub> (3-15%)

Units in parentheses are mass %.

### V. Type and Class (ISO6872:2015)

Type:II/ Class:5

### VI. Physical Properties

Coefficient of Thermal Expansion (25-500°C (77-932°F)): 10.1×10<sup>-6</sup>/K

### VII. Directions for Use

- Take the disc from the packaging and confirm that the disc does not have a crack or other damage.
- Place the disc into the milling machine; then begin the milling process following the milling systems technical instructions.
- After milling, remove the restorations from the disc with a diamond bur, etc.
- The cutting waste or dust, which is attached to the restorations, can be removed with a gentle air stream.
- Put the restorations into the refractory sagger tray and place them in the sintering furnace.
- Depending on the performance of the sintering furnace used, review the sintering schedule shown above (III. Sintering Program) before sintering the restorations.
- After sintering, adjust the restorations with a diamond bur as needed.
- Confirm that the restorations have no cracks.
- 1 FCZ:
  - Glaze baking: Create a high shine surface by polishing, especially on the contact areas, then apply the glaze on all surfaces in the usual manner.
  - Hand polishing: Create a high shine surface on the entire restoration by polishing without using the glaze. When the restoration is to be finished without the glaze, select a brighter shade than the final shade.
- 2 Frameworks: Build-up the dental porcelain (CERABIEN ZR or CZR PRESS LF) on the frameworks following the manufacturer's technical instructions. Check the coefficient of thermal expansion of the porcelain in the manufacturer's technical instructions to confirm compatibility.

### VIII. Remarks on Handling

#### Contraindications:

- If the patient is hypersensitive to zirconia or any other components, this product must not be used.

- Do not use this product to build fixed FCZ bridges or frameworks of more than 6 units and 2 consecutive pontics. (For Canada only)

#### Warning:

- If the patient or the dental professional demonstrates a hypersensitivity reaction, such as rash, dermatitis etc., discontinue use of the product and seek medical attention immediately.

#### Caution:

- This product should not be used when malocclusion, clenching or bruxism conditions are present.
- When milling the disc or cutting, grinding and polishing the restorations, use an approved dust mask and vacuum with air filter to protect your lungs from inhaling dust.
- When milling the disc or cutting, grinding and polishing the restorations, use safety glasses to prevent dust from getting into your eyes. If the dust gets into your eyes, immediately rinse with copious amounts of water and consult a physician.
- Do not use for any purposes except for dental restoration. This product is for dental application only.
- The use of this product is restricted to dental professionals.
- Do not touch the items heated in the furnace with your bare hands.
- Do not take the restorations out of the sintering furnace during high temperature, as the quenching causes the breaking. However, if an automatic opening type furnace is used under the sintering program 2 or 3, the restorations may be removed from the furnace at 800°C (1472°F) or less. When removed from the furnace, the restorations must be put on a tray made of ceramic fiber (for example: Noritake Porcelain Mat) to cool slowly.
- There may be a large difference between the setting temperature in the program and the temperature in the real furnace when using the furnace for sintering program 2 or 3. Please use the furnace after checking with the manufacturer that the furnace and temperature listed in sintering program 2 or 3 is compatible.
- Sintering program 2 or 3 is only recommended for frameworks (up to 3 units), FCZ crowns, FCZ bridges (up to 3 units), inlays, onlays and veneers.**
- The margins should be prepared with a deep chamfer and rounded shoulders, with cutting edges and corners rounded to eliminate sharp preparation corners. The angle of the axial surface should be within the range of 5 to 15 degrees.
- When preparing teeth, avoid the following: deep shoulders, J-margins, knife edges, serrated margins, non-tapered abutments, undercuts, guide grooves, the formation of retentive holes, and sharp corners.
- Keep the following thickness of this product for fabricating prosthetics:

| Location & indication     | Wall thickness               |
|---------------------------|------------------------------|
| Anterior crown or bridge  | 0.4 mm or more               |
| Veneer                    | 0.4 mm or more <sup>1)</sup> |
| Posterior crown or bridge | 0.5 mm or more               |
| Inlay or onlay            | 0.5 mm or more               |

1) 0.4 mm or more of this product is for full zirconia veneers. Keep thickness 0.8 mm or more, if it is used for combination with the porcelain.

- Use the following cross-sectional areas for connectors when fabricating bridges.

| Location & indication | Connector cross section  |                           |
|-----------------------|--------------------------|---------------------------|
| Anterior              | 2- or 3-unit bridges     | 7 mm <sup>2</sup> or more |
|                       | more than 3-unit bridges | 9 mm <sup>2</sup> or more |
| Posterior             | 2- or 3-unit bridges     | 9 mm <sup>2</sup> or more |
|                       | more than 3-unit bridges | 9 mm <sup>2</sup> or more |

- Keep the connected bridge pontics maximum 2 teeth in bridge constructions. When the pontic of 2 teeth continues, keep the connector cross section between the pontic at 12 mm<sup>2</sup>.
- Keep the cantilever bridge to 1 pontic tooth and the connector cross section at 12 mm<sup>2</sup>.
- Choose a shade color that is brighter than the intended color for a thick restoration, as it may look duller depending on the thickness of the restorations.
- When using a sintering furnace for the first time and changing a sintering condition, colors after sintering may vary. Sinter a small piece of Zirconia beforehand and confirm the color.

#### Caution in conjunction with usage:

- Do not use the disc if there is a crack noticed after removing it from the package.
  - Dry milling is recommended. If wet milling is used, the translucency of the disc may be reduced.
  - If you find a crack in the restorations, do not use.
  - The above sintering recommendation is only a guideline; some adjustments may be required depending on each individual furnace.
  - When milling the disc, use caution when approaching the milling bur with the plastic ring to prevent detachment of the disc. Mill the disc as if leaving a zirconia part contacting the internal side of the plastic ring at 2 mm.
  - Cut and remove the plastic ring prior to sintering the whole disc.
  - Dispose of this product as a medical waste to prevent infection.
- Storage:**
- Store in a cool and dry place. Keep away from direct sunlight.
  - The product should be stored at 10-30°C (50-86°F).
  - Do not remove the disc from its packaging during storage.
  - The disc is fragile, and requires care when handling.
  - The product must be stored in an appropriate place where only dental personnel have access.
  - The product must be used by the expiration date indicated on the package.

#### [WARRANTY]

Kuraray Noritake Dental Inc. will replace any product that is proven to be defective. Kuraray Noritake Dental Inc. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

#### [NOTE]

If a serious incident attributable to this product occurs, report it to the manufacturer through the EU importer, and to the regulatory authorities of the country in which the user/patient resides (if the product code ends with "EU" (e.g.#123-4567EU)).


"KATANA" and "CERABIEN" are registered trademarks or trademarks of NORITAKE CO., LIMITED.

---

005 TI-039-EN 08/2025

## HTML PLUS (High Translucent Multi Layered PLUS)

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