PANAVIA™ F 2.0

EN	ENGLISH	English	DENTAL DUAL-CURED ADHESIVE RESIN CEMENT
FR	FRANÇAIS	French	CIMENT RÉSINE ADHÉSIF DENTAIRE À PRISE DUALE
ES	ESPAÑOL	Spanish	CEMENTO DE RESINA ADHESIVO DE POLIMERIZACIÓN DUAL DE USO DENTAL
IT	ITALIANO	Italian	CEMENTO RESINA ADESIVO DENTALE A POLIMERIZZAZIONE DUALE
NL	NEDERLANDS	Dutch	DUAAL HARDEND COMPOSITCEMENT
DE	DEUTSCH	German	DENTALER dualhärtender Kompositzement
SV	SVENSKA	Swedish	DENTAL DUALHÄRDANDE VIDHÅFTANDE RESINCEMENT
NO	NORSK	Norwegian	DENTAL DOBBELTHERDET ADHESIV RESINSEMENT
FI	SUOMI	Finnish	HAMPAAN KAKSOISKOVETTEINEN YHDISTELMÄMUOVISEMENTTI
DA	DANSK	Danish	DENTAL DUALHÆRDENDE ADHÆSIV RESINCEMENT
PT	PORTUGUÊS	Portuguese	CIMENTO DE RESINA ADESIVA DE ACÇÃO DUAL PARA USO DENTÁRIO
EL	ΕΛΛΗΝΙΚΑ	Greek	ΟΔΟΝΤΙΚΗ ΡΗΤΙΝΩΔΗΣ ΣΥΓΚΟΛΛΗΤΙΚΗ ΚΟΝΙΑ ΔΙΠΛΟΥ ΠΟΛΥΜΕΡΙΣΜΟΥ
TR	TÜRKÇE	Turkish	DUAL SERTLEŞEN DENTAL ADEZİV REZİN SİMAN
PL	POLSKI	Polish	PODWÓJNIE UTWARDZANY ADHEZYJNY ŻYWICZNY CEMENT DENTYSTYCZNY
RO	ROMÂNĂ	Romanian	CIMENT DENTAR PE BAZĂ DE RĂȘINI ADEZIVE CU POLIMERIZARE DUBLĂ
HR	HRVATSKI	Croatian	DENTALNI DVOJNO STVRDNJAVAJUĆI ADHEZIVNI CEMENT NA BAZI SMOLE
HU	MAGYAR	Hungarian	FOGÁSZATI, DUÁL KÖTÉSŰ, MŰGYANTA ALAPÚ RAGASZTÓCEMENT
SL	SLOVENŠČINA	Slovenian	DENTALNI DVOJNO STRJUJOČI ADHEZIVNI KOMPOZITNI CEMENT
CS	ČESKY	Czech	DENTÁLNÍ DUÁLNĚ TUHNOUCÍ KOMPOZITNÍ UPEVŇOVACÍ CEMENT
SK	SLOVENSKY	Slovak	DENTÁLNY DUÁLNE VYTVRDZOVANÝ KOMPOZITNÝ ŽIVICOVÝ CEMENT
BG	БЪЛГАРСКИ	Bulgarian	СТОМАТОЛОГИЧЕН ДВОЙНОПОЛИМЕРИЗИРАЩ АДХЕЗИВЕН КОМПОЗИТЕН ЦИМЕНТ
UK	УКРАЇНСЬКА	Ukrainian	СТОМАТОЛОГІЧНИЙ АДГЕЗИВНИЙ КОМПОЗИТНИЙ ЦЕМЕНТ ПОДВІЙНОГО ТВЕРДІННЯ
ET	EESTI	Estonian	DENTAALNE KAKSIKÕVASTUV ADHESIIVNE KOMPOSIITSEMENT
LV	LATVIEŠU	Latvian	ZOBĀRSTNIECĪBAS DIVĒJĀDI CIETĒJOŠS ADHEZĪVS SVEĶU CEMENTS
LT	LIETUVIŠKAI	Lithuanian	ODONTOLOGINIS DVIGUBO KIETĖJIMO DERVINIS CEMENTAS FIKSAVIMUI

Kuraray Noritake Dental Inc.

1621 Sakazu, Kurashiki, Okayama 710-0801, Japan

EMERGO EUROPE

Westervoortsedijk 60,
6827 AT Arnhem, The Netherlands

Kuraray Europe GmbH (EU Importer)

Philipp-Reis-Str. 4,
65795 Hattersheim am Main, Germany
Phone:+49 (0)69 305 35835 Fax:+49 (0)69 305 98 35835
URL:<https://www.kuraraynoritake.eu>

1561R822R

PANAVIA™ F 2.0

I. INTRODUCTION

PANAVIA F 2.0 is a dual-cure (light- and/or self-cure), radiopaque resin-based cement system for metal, composite, and silanated porcelain restorations. The general clinical benefit of PANAVIA F 2.0 is to restore tooth function for the following INDICATIONS FOR USE.

PANAVIA F 2.0 consists of ED PRIMER II, PANAVIA F 2.0 Paste, ALLOY PRIMER and OXYGUARD II.

ED PRIMER II contains HEMA and 5-NMSA as well as MDP and consists of liquid A and liquid B.

PANAVIA F 2.0 Paste releases fluoride.

ALLOY PRIMER improves the bond strength for precious alloy and PANAVIA F 2.0 Paste. OXYGUARD II is an oxygen-blocking agent to allow the PANAVIA F 2.0 Paste to polymerize when not light-cured. Tin plating is not required.

II. INDICATIONS FOR USE

PANAVIA F 2.0 is indicated for the following applications:

- [1] Cementation of metal crowns and bridges, inlays and onlays.
- [2] Cementation of porcelain crowns, inlays, onlays and veneers.
- [3] Cementation of composite resin crowns, inlays, and onlays.
- [4] Cementation of adhesion bridges.
- [5] Cementation of endodontic cores and prefabricated posts.
- [6] Amalgam bonding.

[NOTE]

Use cement shades appropriate to individual cases.

The shades of adhesive cements and applicable cases:

Restoration	Cement shade	TC,	White	Opaque
		Light		
Metal inlays and onlays; metal crowns and bridges		⊙	⊙	⊙
Porcelain or composite inlays, onlays and veneers;		⊙	△	△
Preformed posts and cast metal cores		⊙	⊙	⊙
Adhesion bridges and splints	anterior	△	⊙	⊙
	posterior	⊙	⊙	⊙
Bonded amalgam restorations		⊙	⊙	⊙

⊙:RECOMMENDED ○:ADEQUATE △:NOT RECOMMENDED

III. CONTRAINDICATIONS ⚠

- [1] Patients with a history of hypersensitivity to methacrylate monomers. Patients known to be allergic to any of the ingredients contained in this product.
- [2] Patients with a history of hypersensitivity to acetone

IV. SIDE EFFECTS ⚠

The oral mucosal membrane may turn a whitish color when contacted by ED PRIMER II or ALLOY PRIMER due to the coagulation of protein. This is usually a temporary phenomenon that will disappear in a few days. In individual cases ulceration has been reported.

V. INCOMPATIBILITY ⚠

- [1] Do not use eugenol containing materials for pulp protection and temporary sealing since they retard the curing process.
- [2] Do not use hemostatics containing ferric compounds, since they may impair adhesion and the remaining ferric ion may cause discoloration in the margin or surrounding gingiva.
- [3] Do not use PANAVIA F 2.0 Paste with PANAVIA F Paste, or the mixed paste will not cure completely by light curing.

VI. PRECAUTIONS ⚠

1. Safety precautions

1. This product contains substances that may cause allergic reactions. Avoid use of the product in patients with known allergies to methacrylate monomers or any other components.
2. If the patient demonstrates a hypersensitivity reaction, such as rash, eczema, features of inflammation, ulcer, swelling, itching or numbness, discontinue use of the product and seek medical attention.
3. Use caution to prevent the product from coming in contact with the skin or getting into the eye. Before using the product, cover the patient's eyes with a towel to protect them in the event of splashing material.
4. If the product comes in contact with human body tissues, take the following actions:
 - <If the product gets in the eye>
 - Immediately wash the eye with copious amounts of water and consult a physician.
 - <If the product comes in contact with the skin or the oral mucosa>
 - Immediately wipe the area with a cotton pellet or a gauze moistened with alcohol, and rinse with copious amounts of water.
5. Use caution to prevent the patient from accidentally swallowing the product.
6. If a patient or operator feels sick from inhaling the acetone contained in the product, allow them to rest and breathe fresh air.
7. Any actual or near pulp exposure area should be covered with a hard setting calcium hydroxide material. Do not use eugenol materials for pulp protection.
8. When using with preformed stainless posts, the post should not contact metal restorations. Cover the post with composite resin.
9. Avoid looking directly at the curing light when curing the product.
10. Dispose of this product as a medical waste to prevent infection.

2. Handling and manipulation precautions

1. PANAVIA F 2.0 polymerizes by a dental curing light (irradiate wave length : 400-515nm). Use the light blocking plate to avoid exposing the material to an operating light or natural light (sunlight from windows).
2. Make sure the disposable nozzle or disposable brush tip is securely attached to prevent the patient from swallowing them.
3. After placing the restoration on the tooth, the cement could cure due to the operatory light. Use caution not to let the operatory light come too close to the patient.

3. Storage precautions

1. Do not use after the expiration date. Note expiration date on the outside of package.
2. ALLOY PRIMER is flammable. Keep away from flame.
3. The product must be refrigerated (2-8°C/36-46°F) when not in use; bring it to room temperature before using.
4. ALLOY PRIMER must be stored at 2-25°C/36-77°F when not in use.
5. Keep away from extreme heat or direct sunlight.
6. The bottle or syringe cap should be replaced as soon as the resin has been dispensed from the bottle or syringe. This prevents evaporation of volatile ingredients.
7. The product must be stored in a proper place where only dental practitioners can access.

VII. COMPONENTS

1. Shades

PANAVIA F 2.0 Paste is available in 4 shades; TC, Light, White or Opaque

2. Components

Please see the outside of the package for quantity.

- 1) PANAVIA F 2.0 Paste (A Paste / B Paste)
- 2) ED PRIMER II (Liquid A/Liquid B)
- 3) ALLOY PRIMER
- 4) OXYGUARD II
- 5) Accessories

- Spatula
- Mixing pad
- Disposable brush tips
- Brush tip handle
- Mixing dish
- Light blocking plate
- Disposable nozzles

3. Ingredients

- 1) PANAVIA F 2.0 Paste (A Paste/B Paste)

Principal ingredients

(1) A Paste

- 10-Methacryloyloxydecyl dihydrogen phosphate (MDP) (3-7%)
- Hydrophobic aromatic dimethacrylate (1-15%)
- Hydrophobic aliphatic dimethacrylate (1-10%)
- Hydrophilic aliphatic dimethacrylate (< 0.5%)
- Silanated silica filler (> 70%)
- Silanated colloidal silica (1-5%)
- dl-Camphorquinone (< 0.1%)
- Benzoyl peroxide (< 1%)
- Initiators (< 0.1%)

(2) B Paste

- Hydrophobic aromatic dimethacrylate (5-25%)
- Hydrophobic aliphatic dimethacrylate (1-10%)
- Hydrophilic aliphatic dimethacrylate (1-10%)
- Silanated barium glass filler (40-80%)
- Surface treated sodium fluoride (≤ 10%)
- Catalysts (< 1%)
- Accelerators (< 1%)
- Silanated titanium dioxide (< 10%)
- Pigments (< 0.1%)

The total amount of inorganic filler is approx. 59 vol%. The particle size of inorganic fillers ranges from 0.04 µm to 19 µm.

- 2) ED PRIMER II

Principal ingredients

(1) Liquid A

- 2-Hydroxyethyl methacrylate (HEMA) (30-50%)
- 10-Methacryloyloxydecyl dihydrogen phosphate (MDP) (20-40%)
- Water (5-45%)
- N-Methacryloyl-5-aminosalicylic acid (5-NMSA) (< 1.5%)
- Accelerators (1-10%)

(2) Liquid B

- N-Methacryloyl-5-aminosalicylic acid (5-NMSA) (< 1%)
- Water (> 80%)
- Catalysts (1-10%)
- Accelerators (1-10%)

- 3) ALLOY PRIMER

Principal ingredients

- Acetone (> 90%)
- 10-Methacryloyloxydecyl dihydrogen phosphate (MDP) (< 0.5%)
- 6-(4-Vinylbenzyl-N-propyl)amino-1,3,5-triazine-2,4-dithione (< 1%)

- 4) OXYGUARD II

Principal ingredients

- Glycerol
- Polyethyleneglycol
- Catalysts
- Accelerators
- Dyes

Units in parentheses are mass %.

VIII. RELATED PRODUCTS

The following products are necessary for specific procedures.

1) CLEARFIL PORCELAIN BOND ACTIVATOR

* This product contains a silane coupling agent. Mixing it with CLEARFIL SE BOND PRIMER or CLEARFIL LINER BOND 2V PRIMER or CLEARFIL NEW BOND or CLEARFIL PHOTO BOND improves the bond strength to porcelain or cured composite.

2) K ETCHANT GEL

* This phosphoric acid solution is used for pretreatment of uncut enamel and porcelain.

IX. CLINICAL PROCEDURES

1. Basic procedure (use of adhesive cement)

[NOTE]

Use the mixed paste as soon as possible after dispensing and mixing.

1) Dispensing the pastes

- Align marking on the nut with the reference line on the plunger and turn the syringe to dispense the necessary amount of A paste minimum turn of the syringe should be half a turn.
- Equal amounts of A paste and B paste should be dispensed.
- The amount of paste dispensed at last rotation of the syringe could be inaccurate. Therefore, discard the syringe before using the last portion.
- The necessary amount of paste for a typical applications is:

Number of rotations of syringe	Applications
Half a turn	Inlays and onlays
1 turn	Crowns

[NOTE]

- If the paste is dispensed by turning the plunger a quarter of a turn, the performance of the product could be impaired when the paste hardens.
- If not used immediately, it should be covered with a light-blocking plate.

2) Mixing A paste and B paste

Mix sufficient A paste and B paste on the mixing plate for 20 seconds. Be sure there is no water mist on the mixing plate or spatula before using them; the presence of water could shorten the working time of the mixed paste.

[CAUTION]

- The working time of mixed paste could vary if mixing is insufficient.
- The paste must be used within 3 minutes after mixing.

[REMARKS]

The working times of PANAVIA F 2.0 Paste from dispensation to the completion of cementation are :

Working time of PANAVIA F 2.0

	Working steps	Working time
1.	Dispensing the pastes (by turning the syringe the same amount of turn)	15 min.
2.	Mixing the pastes (for 20 sec.)	3 min.
3.	Placing the restoration under pressure	60 sec.
	---In the case of root canal	40 sec.
4.	Light curing	
	---Conventional halogen, LED ^{*1}	20 sec.
	---Plasma arc, fast halogen ^{*2}	5 sec.
	Applying OXYGUARD II	3 min.

Light intensity of approved curing lights (400~500nm)

*1) Conventional halogen (>250mW/cm²), LED (>160mW/cm²)

*2) Plasma arc (>2000mW/cm²), fast halogen (>550mW/cm²)

2. Standard procedure I (indications 1 to 4: for cementation)

[1] Surface treatment of restoration

1. Precious metals (crowns, bridges, inlays and onlays)

1) Sandblasting (as necessary)

Sandblast the restoration surface using 30-50 micron alumina particles at an air pressure of 4.2 - 7 kg/cm² (60-100 PSI); 2-3 seconds per cm² will remove the luster producing a matte finish.

2) Ultrasonic cleaning

Clean the restoration surface in an ultrasonic unit for 2 minutes.

3) Applying ALLOY PRIMER

Apply a thin coat of ALLOY PRIMER to the precious metal alloy.

[CAUTION]

If the adherend surface is contaminated with saliva or blood after ultrasonic cleaning, clean the adherend surface in the ultrasonic unit using a neutral detergent and then wash it for 1 minute with running water.

2. Non-precious metals

1) Sandblasting (as necessary).

Sandblast the restoration surface using 30-50 micron alumina powder at an air pressure of 4.2-7 kg/cm² (60-100 PSI); 2-3 seconds per cm² will remove the luster producing a matte finish.

2) Ultrasonic cleaning

Clean the restoration surface in an ultrasonic unit for 2 minutes.

[CAUTION]

If the adherend surface is contaminated with saliva or blood after ultrasonic cleaning, clean the adherend surface in an ultrasonic unit using a neutral detergent and wash it for one minute with running water.

3. Porcelain and cured composite restoration (inlays, onlays, crowns and veneers)

1) Sandblasting

Sandblast the restoration surface using 30-50 micron alumina particles at a low air pressure. (1-2 kg/cm² (14-28 PSI))

2) Etching with phosphoric acid

Etch the adherend surface with K ETCHANT GEL.

3) Rinse and dry

After etching with phosphoric acid, rinse the adherend surface with water and dry.

4) Silane coupling treatment

Silanate the adherend surface using the following:

Application of the mixture of CLEARFIL PORCELAIN BOND ACTIVATOR and CLEARFIL SE BOND PRIMER or CLEARFIL LINER BOND 2V PRIMER or CLEARFIL PHOTO BOND or CLEARFIL NEW BOND.

[NOTE]

After treating the restoration surface, proceed to cementation quickly.

[2] Adherend surface treatment

1. Cleaning cavity or abutment tooth surface

- Remove temporary sealing agent or temporary cementation agent from the adherend surface.
- When cementing to uncut enamel or using with adhesive bridge or porcelain laminate veneers, apply K ETCHANT GEL to the adherend surface for 10 seconds.

2. Adherend surface treatment

1) Mixing ED PRIMER II

Dispense one drop each of Liquid A and Liquid B into well of the mixing dish and mix.

2) Applying ALLOY PRIMER

If a precious metal abutment tooth is used, apply ALLOY PRIMER to its metal surface.

[CAUTION]

If the adherend surface is contaminated with saliva or blood after ALLOY PRIMER is applied, clean the adherend surface with a cotton pledget moistened with alcohol and apply ALLOY PRIMER again.

3) Applying ED PRIMER II

Apply ED PRIMER II to the entire tooth surface (enamel and dentin) of the adherend surface or metal or composite resin abutment tooth with a disposable brush tip or sponge and leave it in place for 30 seconds.

4) Drying

Using a sponge or paper point, remove excess primer to prevent the formation of a pool of the primer in the root canal or at the corners of the cavity. Dry the primer completely with gentle air flow. Remember that a pool of the primer will cause quick polymerization of the adhesive cement. Also do not rinse. To prevent the primer from splattering, it is good practice to dry while using a vacuum.

[CAUTION]

ED PRIMER II should be applied on the entire surface of the tooth structure. Do not apply it to the restoration.

[3] Preparing PANAVIA F 2.0 Paste

Prepare the adhesive cement according to the basic clinical procedure. Refer to IX.1."basic procedure".

[4] Cementation

1. Applying mixed paste to the restoration

Apply the mixed paste to the restoration.

[CAUTION]

DO NOT apply PANAVIA F 2.0 Paste to the tooth surface primed with ED PRIMER II as this will accelerate the set of PANAVIA F 2.0 Paste.

2. Cementing restoration

Cement the restoration to the cavity or the abutment tooth.

Cementation should be completed within 60 seconds.

[CAUTION]

When the adhesive cement comes in contact with ED PRIMER II, the polymerization of the adhesive cement is accelerated.

3. Removing excess paste

Any excess PANAVIA F 2.0 Paste remaining at the margin can be removed with an explorer or small scaler. The restoration can then be finished and polished with pumice and water.

4. Curing cement margin

Cure the mixed paste along the cement margin, using either of the following two methods.

① Light curing method

When it is possible to light cure adhesive cement along the cement margin, such as inlays and onlays, light cure each section of the cement margin for 20 seconds by conventional halogen curing lights (>250m W/cm²) or LED curing lights (>160mW/cm²). If plasma arc curing lights (>2000mW/cm²) or fast halogen curing lights (>550mW/cm²) are used, each section of the cement margin can be cured for 5 seconds.

[CAUTION]

The Opaque paste should not be light cured but allowed to cure by using OXYGUARD II. It has a low curing depth.

② OXYGUARD II

Use OXYGUARD II to cure the mixed paste as follows:

With a disposable brush tip, apply OXYGUARD II to the margin.

After 3 minutes remove OXYGUARD II with a cotton roll and water spray.

[5] Finishing

Remove excess cement adhered to tooth surface by polishing.

3. Standard procedure II (Indication 5: core build up)

[NOTE]

This procedure is for use with a preformed post and composite resin core building-up. For the cementation of metal cores, refer to standard procedure 1. And according to the instruction for use of the post and composite resin.

[1] Surface treatment of the post

1. Sandblasting

Sand blast the post as necessary.

[NOTE]

Some preformed posts do not require sandblasting. Refer to the Instruction for use of the specific post.

2. Applying ALLOY PRIMER

Apply ALLOY PRIMER to the post if it is precious metal alloy.

[2] Cleaning cavity and preparing root canal

Remove the temporary sealing agent from the cavity and filling material from the root canal. Using a Pizo reamer, prepare and clean the root canal opening.

[3] Tooth surface treatment

1. Mixing ED PRIMER II

Dispense one drop each of Liquid A and Liquid B on the mixing dish and mix.

2. Applying ED PRIMER II

Using a sponge or a broach cotton, apply the mixture to the root canal, on the root surface and the tooth structure. Leave it in place for 30 seconds.

3. Removing excess primer (the same step is also needed in the case of metal cores)

Using a sponge, broach cotton or paper point, remove excess primer to prevent the primer from pooling at the corners of the cavity and inside the root canal.

4. Drying

Dry the primer with gentle air flow. It is good practice to dry while using a vacuum to prevent the primer from splattering.

[CAUTION]

Dry the primer completely. A pool of the primer at the corners of the cavity or inside the root canal will cause quick polymerization of the mixed paste.

[4] Preparing PANAVIA F 2.0 Paste

Prepare the adhesive cement according to the basic clinical procedure. Refer to IX.1. "basic procedure".

[5] Seating the post

1. Apply the mixed paste to the post.

[REMARKS]

The mixed paste is applied to the metal post for cementation.

2. Seating the post into the root canal

After applying the mixed paste to the post, insert it into the root canal quickly. It is advisable to lightly vibrate the post while inserting it into the root canal to prevent the entrapment of air bubbles.

[CAUTION]

If multiple posts are to be fitted to a single tooth, use caution to prevent excess cement from flowing into other root canals.

[CAUTION]

Never use a lentulo spiral for loading the adhesive cement into the root canal. If the adhesive cement is loaded into the root canal using a composite resin syringe, the polymerization of the cement is accelerated. It is necessary to fit the post as quickly as possible.

3. Spreading excess cement

Using a small brush, spread excess cement over the remaining crown and post head.

4. Curing adhesive cement

Light cure the adhesive cement on the remaining crown and post head for 20 seconds by conventional halogen curing lights (>250mW/cm²) or LED curing lights (>160mW/cm²). If plasma arc curing lights (>2000mW/cm²) or fast halogen curing lights (>550mW/cm²) are used, each cement margin has to be cured for 5 seconds.

[NOTE]

If it is difficult to light cure if the opaque cement is used; use the core build-up composite resin.

[6] Build-up composite resin

Build-up the composite resin for preparation of the abutment tooth according to the instruction for use.

[7] Curing and finishing of composite resin

After curing the composite resin, prepare an abutment tooth.

4. Standard procedure III (indications 6: Amalgam Bonding)

[1] Cleaning of tooth structure

Cleaning cavity or abutment surface

Remove temporary sealing material or temporary cementation material from the adherend surface.

[2] Treatment of tooth surface

Treating with ED PRIMER II

1) Preparing ED PRIMER II

Dispense one drop each of Liquid A and Liquid B on the mixing plate and mix.

2) Applying ED PRIMER II

Apply ED PRIMER II to the entire adherend tooth surface (enamel and dentin), metal, or composite resin abutment with a small brush or sponge and leave it in place for 30 seconds.

3) Drying

Using a sponge or paper point, remove excess primer to prevent the primer from pooling in the corners of the cavity. Dry the primer completely by using gentle air flow. Remember that a pool of the primer will cause quick polymerization of the adhesive cement. To prevent the primer from splattering, it is good practice to dry while using a vacuum.

[3] Preparing of adhesive cement

Prepare the adhesive cement according to the basic clinical procedure. Refer to IX.1. "basic procedure".

[4] Placing the amalgam

1) Apply the adhesive cement to the cavity

Apply a thin, even layer of the mixed adhesive cement to the entire cavity surface primed with ED PRIMER II being careful to avoid entrapping air.

[CAUTION]

Because ED PRIMER II accelerates the set of the adhesive cement, the adhesive cement should be applied to the primed cavity quickly.

2) Amalgam filling

The triturated amalgam should be condensed on the unset adhesive cement. Occlusal carving can be accomplished in the normal manner.

[5] Removing excess cement

Any slight excess of PANAVIA F 2.0 Paste remaining at the margin can be removed with an explorer or small scaler.

[6] Curing adhesive cement

Cure the mixed paste along the cement margin, using either of the following two methods.

① Light curing method

When it is possible to light cure adhesive cement along the cement margin, such as inlays and onlays, light cure each section of the cement margin for 20 seconds by conventional halogen curing lights (>250mW/cm²) or LED curing lights (>160mW/cm²). If plasma arc curing lights (>2000mW/cm²) or fast halogen curing lights (>550mW/cm²) are used, each section of the cement margin has to be cured for 5 seconds.

[CAUTION]

The opaque paste should not be light cured but allowed to cure by using OXYGUARD II. It has a low curing depth.

② OXYGUARD II

Use OXYGUARD II to cure the mixed paste as follows:

With a disposable brush tip apply OXYGUARD II to the margin.

After 3 minutes remove OXYGUARD II with a cotton roll and water spray.

[7] Finishing

Remove excess cement adhered to tooth structure by polishing.

[WARRANTY]

Kuraray Noritake Dental Inc. will replace any product that is proved to be defective. Kuraray Noritake Dental Inc. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.


[NOTE]

If a serious incident attributable to this product occurs, report it to the manufacturer's authorized representative shown below and the regulatory authorities of the country in which the user/patient resides.

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"CLEARFIL", "PANAVIA" and "OXYGUARD" are registered trademarks or trademarks of KURARAY CO., LTD.
„CLEARFIL“, „PANAVIA“ ir „OXYGUARD“ yra „KURARAY CO., LTD“ registruotieji prekės ženklai arba prekės ženklai.