PANAVIA™ Veneer LC

EN	ENGLISH	English	DENTAL ADHESIVE RESIN CEMENT SYSTEM
FR	FRANÇAIS	French	SYSTÈME DE CIMENTATION DENTAIRE PAR RÉSINE ADHÉSIVE
ES	ESPAÑOL	Spanish	SISTEMA DE CEMENTO DE RESINA ADHESIVO DENTAL
IT	ITALIANO	Italian	SISTEMA CON CEMENTO RESINA ADESIVO DENTALE
NL	NEDERLANDS	Dutch	TANDHEELKUNDIG ADHESIEF COMPOSITCEMENTSYSTEEM
DE	DEUTSCH	German	DENTALES ADHÄSIVES KUNSTSTOFFZEMENTSYSTEM
SV	SVENSKA	Swedish	DENTAL VIDHÅFTANDE KOMPOSITCEMENT SYSTEM
NO	NORSK	Norwegian	DENTALT ADHESIVT RESINSEMENTSYSTEM
FI	SUOMI	Finnish	HAMMASLÄÄKETIETEELLINEN ADHESIIVINEN RESIINISEMENTTIJÄRJESTELMÄ
DA	DANSK	Danish	DENTALT ADHÆSIVT RESINCEMENTSYSTEM
PT	PORTUGUÊS	Portuguese	SISTEMA DE CIMENTO DE RESINA ADESIVO DENTÁRIO
EL	ΕΛΛΗΝΙΚΑ	Greek	ΟΔΟΝΤΙΑΤΡΙΚΗ ΣΥΓΚΟΛΛΗΤΙΚΗ ΚΟΝΙΑ ΠΗΤΙΝΗΣ
TR	TÜRKÇE	Turkish	DENTAL ADEZİV REZİN SİMAN SİSTEMİ
PL	POLSKI	Polish	STOMATOLOGICZNY SYSTEM ADHEZYJNEGO CEMENTU NA BAZIE ŻYWICY
RO	ROMÂNĂ	Romanian	SISTEM CIMENT DENTAR ADEZIV PE BAZĂ DE RĂȘINI
HR	HRVATSKI	Croatian	DENTALNI ADHEZIVNI SUSTAV NA BAZI CEMENTNE SMOLE
HU	MAGYAR	Hungarian	FOGÁSZATI ADHEZÍV MŰGYANTA CEMENTRENDSZER
SL	SLOVENŠČINA	Slovenian	DENTALNI LEPILNI SISTEM IZ KOMPOZITNEGA CEMENTA
CS	ČESKY	Czech	DENTÁLNÍ ADHEZIVNÍ SYSTÉM PLASTOVÉHO CEMENTU
SK	SLOVENSKY	Slovak	LEPIACI SYSTÉM Z UMELEJ HMOTY NA CEMENTOVANIE ZUBOV
BG	БЪЛГАРСКИ	Bulgarian	СТОМАТОЛОГИЧНА СИСТЕМА С АДХЕЗИВЕН КОМПОЗИТЕН ЦИМЕНТ
UK	УКРАЇНСЬКА	Ukrainian	АДГЕЗИВНА КОМПОЗИТНА ЦЕМЕНТНА СИСТЕМА
ET	EESTI	Estonian	DENTAALNE ADHESIIVNE VAIKTSEMENDISÜSTEEM
LV	LATVIEŠU	Latvian	ZOBĀRSTNIECĪBAS ADHEZĪVĀ PLASTMASAS CEMENTA SISTĒMA
LT	LIETUVIŠKAI	Lithuanian	ODONTOLOGINĖ ADHEZYVINĖ DERVINIO CEMENTO SISTEMA

Kuraray Noritake Dental Inc.

1621 Sakazu, Kurashiki, Okayama 710-0801, Japan

EMERGO EUROPE

Westervoortsedijk 60,
6827 AT Arnhem, The Netherlands

Kuraray Europe GmbH (EU Importer)


Philipp-Reis-Str. 4,
65795 Hattersheim am Main, Germany

Phone:+49 (0)69 305 35835 Fax:+49 (0)69 305 98 35835

URL:<https://www.kuraraynoritake.eu>

PANAVIA™ Veneer LC

I. INTRODUCTION

PANAVIA Veneer LC is a light-cure adhesive resin cement system. It consists of the PANAVIA Veneer LC Paste (a light-cure cement paste), CLEARFIL CERAMIC PRIMER PLUS (a silane, zirconia & metal primer), K-ETCHANT Syringe (phosphoric acid), and PANAVIA V5 Try-in Paste. PANAVIA V5 Tooth Primer or CLEARFIL Universal Bond Quick can be chosen for tooth treatment.

The cement paste is a light-cure, resin based material which provides color stability and has radiopacity equal to or greater than 1mm aluminum. It is supplied in Kuraray's ergonomic syringe and dispensed via an angled applicator tip (16G) into an inlay cavity or onto an onlay or a laminate veneer. It is available in 4 shades; Universal (A2), Clear, Brown (A4) and White. It is classified as a Type 2 and Class 2 (Group 1) material by ISO 4049. Shade adaptation may be checked with the PANAVIA V5 Try-in Pastes before cementation.

The general clinical benefit of PANAVIA Veneer LC is to restore tooth function for the following INDICATIONS FOR USE.

II. INDICATIONS FOR USE

Cementation of ceramic and composite inlays, onlays and laminate veneers with less than 2mm thickness.

III. CONTRAINDICATIONS

Patients with a history of hypersensitivity to methacrylate monomers. Patients known to be allergic to any of the ingredients contained in this product.

IV. POSSIBLE SIDE EFFECTS

- [1] The oral mucosal membrane may turn whitish when contacted by the product due to the coagulation of protein. This is usually a temporary phenomenon that will disappear in a few days. Instruct patients to avoid irritating the affected area while brushing.
- [2] K-ETCHANT Syringe (phosphoric acid) may cause inflammation or erosion due to its chemistry. Use caution to prevent the product from coming in contact with the skin or getting into the eye.

V. INCOMPATIBILITIES

- [1] Do not use eugenol-containing materials for pulp protection or temporary sealing, since the eugenol may cause discoloration and might retard the curing process.
- [2] Do not use hemostatic agents containing ferric compounds since these materials may impair adhesion and may cause discoloration of the tooth margin or surrounding gingiva due to remaining ferric ions.
- [3] Do not use a hydrogen peroxide solution for cleaning cavities since it may weaken the bonding to the tooth structure.
- [4] Do not use a sodium hypochlorite solution since it may weaken the bonding to the tooth structure.

VI. PRECAUTIONS

1. Safety precautions

1. This product contains substances that may cause allergic reactions. Avoid use of the product in patients with known allergies to methacrylate monomers or any other components.
2. If the patient demonstrates a hypersensitivity reaction, such as rash, eczema, features of inflammation, ulcer, swelling, itching or numbness, discontinue use of the product and seek medical attention.
3. Avoid direct contact with the skin and/ or soft tissue to prevent hypersensitivity. Wear gloves or take appropriate precautions when using the product.
4. Exercise caution to prevent the product from coming in contact with the skin or getting into the eye. Before using the product, cover the patient's eyes with a towel or safety glasses to protect them in the event of splashing material.
5. If the product comes in contact with human body tissues, take the following actions:
 - <If the product gets in the eye>
 - Immediately wash the eye with copious amounts of water and consult a physician.
 - <If the product comes in contact with the skin or the oral mucosa>
 - Immediately wipe the area with a cotton pellet or a gauze pad moistened with alcohol, and rinse with copious amounts of water.
6. Exercise caution to prevent the patient from accidentally swallowing the product.
7. Do not reuse the applicator brush, applicator tip (16G) and needle tip. To prevent cross-contamination, avoid using dispensed liquid/liquid and those brush and tip for different patients. The brush and tip are single-use and should be discarded after using.
8. Clean the cavity sufficiently to prevent poor bonding. If the adherent surface is contaminated with saliva or blood, wash it thoroughly and dry. We advise to use KATANA Cleaner for optimal performance of adhesion. Then re-apply the primer or bonding agent.
9. This product contains the Ytterbium Trifluoride (less than 15%). Using on children under the age of 6 may have a potential risk of fluorosis.
10. Dispose of this product as a medical waste to prevent infection. The tip must be disposed of after covering the tip of the needle to prevent injury.
11. The product should be used only for restorations that are less than 2.0 mm thickness and of sufficient translucency. Metal oxide ceramic (e.g. zirconia) should not be used because it may not have sufficient translucency. However, zirconia of our products with high translucency (KATANA Zirconia STML/UTML, KATANA Zirconia Block STML) can be used with a thickness of less than 1.2 mm. KATANA Zirconia YML can be used with a thickness of less than 1.2 mm only as laminate veneers.

2. Handling and manipulation precautions

[Common precautions]

1. The product must not be used for any purposes other than specified in [II.INDICATIONS FOR USE].
2. The use of this product is restricted to dental professionals.
3. Do not use the product as a provisional cement. This material is designed to use as a permanent cement.
4. Amalgam or other lining materials remaining in the cavity will prevent the passage of light and the polymerization of the product. Completely remove any lining material when preparing the cavity.
5. Use a pulp capping agent in a cavity close to the pulp or in the event of accidental pulp exposure.
6. Use a rubber dam to prevent contamination and to control moisture.
7. When using temporary materials containing tannin or magnesia, completely remove it to prevent the discoloration.
8. When using hemostatic agents containing aluminum chloride, minimize its quantity, and use caution to prevent contact with the adherent surface. Failure to do so might weaken the bond strength to the tooth structure.
9. Do not use hemostatic agents containing ferric sulphate. Ferric sulphate may cause the discoloration and weaken the bond strength to the tooth structure.
10. Do not mix the product with other materials. The mixed materials may cause a change in physical properties, including a possible decrease in the expected results.
11. Wear gloves or take other appropriate protective measures to prevent the occurrence of hypersensitivity that may result from contact with methacrylate monomers or any other components in the product.
12. The liquid bottles should be tightly capped immediately after use to reduce the evaporation of the volatile solvent.
13. If the containers and/ or instruments for this product are damaged, protect yourself from any danger and immediately discontinue their use.

[PANAVIA Veneer LC Paste]

1. Be careful to prevent unnecessary exposure to direct sunlight or operating lights. The paste contains a light-cure catalyst that is highly photo-reactive. During cementation, adjust the angle and/ or distance of the dental light to reduce the intensity of light entering the oral cavity to prevent premature polymerization of the paste.
2. Excess cement can be removed after tack light-curing it for 1 second. When removing the excess cement, hold the restoration in place to avoid the possibility of lifting the restoration, since there could be some insufficiently cured resin cement.
3. When attaching the applicator tip (16G), turn the tip clockwise and attach it securely, thereby preventing paste from leaking at the junction of the applicator tip (16G) and the syringe. In addition, make sure there is no residual paste at the junction of the syringe, which could cause the applicator tip (16G) to fall off.
4. Before wiping the residual paste adhering to the applicator tip (16G) or the junction of the syringe with an alcohol gauze pad, squeeze the pad to remove excess alcohol. Use of excessive alcohol in the gauze pad can cause penetration of alcohol into the tip and dilute the paste. In such cases, it may cause a decrease in physical properties from the expected results.
5. To expel air from the applicator tip (16G) nozzle and prevent to mix air bubbles to the paste, set the nozzle upward and push the plunger slowly (outside of the patient's mouth) until the paste reaches the nozzle.
6. After the paste has been dispensed, the syringe should be capped securely as soon as possible, to prevent the paste from hardening by ambient light, and to prevent foreign matter from entering the syringe.

[PANAVIA V5 Try-in Paste]

1. The use of the Try-in Paste should be limited to checking the shade matching with the PANAVIA V5 Paste or PANAVIA Veneer LC Paste.
2. The Try-in Paste does not set. Do not use it to cement restorations.
3. A transparent liquid may appear at the tip of the Try-in Paste syringe. If this transparent liquid is present, it should be dispensed and discarded, since the separated liquid may affect the shade matching.
4. Shade evaluation with Try-in Paste should be performed using the same approximate thickness of Try-in Paste as the hardened cement.
5. After use, thoroughly wash the Try-in Paste from the restoration and tooth surface with water to avoid poor adhesion.

[PANAVIA V5 Tooth Primer]

1. Use within 5 minutes after dispensing.
2. Only use with PANAVIA V5 Paste or PANAVIA Veneer LC Paste.
3. If the treated surface is contaminated, wash it with water and dry, or clean with alcohol, or KATANA Cleaner, and treat it again with this primer.

[CLEARFIL CERAMIC PRIMER PLUS]

1. Use at once after dispensing. It contains volatile ethanol. As the solvent evaporates, the viscosity increases and may make it difficult to apply.
2. Perform cementation soon after treating the restoration surface with this primer.
3. If the pre-treated/treated surface is contaminated with saliva or blood, wash it with water, dry, clean with K-ETCHANT Syringe or KATANA Cleaner, and re-treat it with this primer.

[K-ETCHANT Syringe (phosphoric acid)]

1. Be careful not to contaminate the treated surface with saliva or blood. If the treated surface is contaminated, re-treat.
2. If the product adheres to clothing, wash it off with water.
3. After each use, remove the needle tip from the syringe and recap the syringe immediately and tightly.
4. When using Tooth Primer, etching vital dentin may cause post-operative sensitivity.

[Dental light-curing unit]

1. Use it according to the Instructions for Use for the dental light-curing unit.
2. Do not look directly at the light source. Protective glasses are recommended.
3. Low light intensity causes poor adhesion. Check the lamp for service life and the dental curing light guide tip for contamination. It is advisable to check the dental curing light intensity using an appropriate light evaluating device at appropriate intervals.
4. The emitting tip of the dental curing unit should be held as near and vertical to the resin surface as possible. If a large resin surface has to be light-cured, it is advisable to divide the area into several sections and light-cure each section separately.

- Check the conditions required to cure the paste by referring to the light-curing times listed in these instructions for use before using the product.

3. Storage precautions

- The product must be used by the expiration date indicated on the package.
- The product, except for PANAVIA V5 Try-in Paste and K-ETCHANT Syringe, must be stored (2-8°C/ 36-46°F) when not in use. PANAVIA V5 Try-in Paste and K-ETCHANT Syringe should be stored at (2-25°C/ 36-77°F).
All components stored in the refrigerator should be brought to room temperature for 15 minutes before use in order to restore its normal viscosity and curing properties.
- The product should be kept away from extreme heat, direct sunlight or a flame.
- The product must be stored in a proper place where only dental practitioners can access.

VII. COMPONENTS

Please see the outside of the package for contents and quantity.

<Principal ingredients>

- PANAVIA Veneer LC Paste: Universal (A2), Clear, Brown (A4), White
Silanated spherical silica filler (40-70%), Urethane dimethacrylate (UDMA) (10-25%), Ytterbium trifluoride (10-20%), Silanated silica filler (5-15%), Triethyleneglycol dimethacrylate (TEGDMA) (5-10%), Hydrophilic aliphatic dimethacrylate (1-5%), Hydrophilic amide monomer (< 3%), Initiators (< 1%), Accelerators (< 0.1%), Phenyl bis(2,4,6-trimethylbenzoyl)-phosphine oxide (< 0.1%), dl-Camphorquinone (< 0.1%), Pigments (< 0.1%)

The total amount of inorganic filler is approximately 47 vol%.

The particle size of inorganic fillers range from 0.05 µm to 8 µm.

- PANAVIA V5 Try-in Paste: Universal (A2), Clear, Brown (A4), White
Glycerol, Silanated colloidal silica, Silanated silica, Colloidal silica, Pigments
- PANAVIA V5 Tooth Primer
10-Methacryloyloxydecyl dihydrogen phosphate (MDP) (10-20%), 2-Hydroxyethyl methacrylate (25-50%), Hydrophilic aliphatic dimethacrylate (10-20%), N,N'-Dimethylaminoethyl methacrylate (< 3%), Accelerators (< 1%), Water (10-50%)
- CLEARFIL CERAMIC PRIMER PLUS
3-Methacryloxypropyl trimethoxysilane (< 10%), 10-Methacryloyloxydecyl dihydrogen phosphate (MDP) (< 2%), Ethanol (> 80%)
- K-ETCHANT Syringe
Phosphoric acid, Water, Colloidal Silica, Pigment
- Accessories
Applicator tip (16G)(For Paste), Needle tip (E) (for K-ETCHANT Syringe), Applicator brush (fine <silver>), Mixing dish*
*Consumables

Units in parentheses are mass %.

VIII. CLINICAL PROCEDURES

1. Cleaning the prepared tooth (tooth, composite resin)

Remove the temporary materials from the prepared tooth in the usual manner; clean the cavity using moisture control.

[NOTE]

KATANA Cleaner is available to clean the prepared tooth, when contamination occurs. Apply it for 10 seconds with rubbing motion, then rinse with water until the color has been disappeared and dry.

2. Trial fitting and adjusting a restoration

- Trial fit the prosthetic restoration to check its fit on the prepared tooth, as necessary.
- If necessary, apply the selected shade of Try-in Paste onto the cementation surface of the restoration and trial fit the restoration onto the prepared tooth. Remove excess Try-in Paste from the margins with a brush. The shades of Try-in Paste correspond to those of the hardened cement paste. Check the shade for best color matching, and then remove the restoration. Using water, completely wash the Try-in Paste from the surface of the restoration and the prepared tooth surface.

3. Conditioning the prosthetic restoration surface

Follow the Instructions for Use of the restorative material.

In the absence of specific instructions, we recommend the following procedures and application of CLEARFIL CERAMIC PRIMER PLUS:

If the adherent surface is silica-based ceramic (e.g. conventional porcelain, lithium disilicate), hybrid ceramic or composite resin:
Based on the type of restoration, Acid Treatment or Blasting Treatment may be used:
Acid Treatment (e.g. conventional porcelain, lithium disilicate):

- Etch the adherent surface with a hydrofluoric acid solution in accordance with the Instructions for Use of the solution.
- Rinse the adherent surface with water and dry.

Blasting Treatment (e.g. hybrid ceramic, composite resin):

- Roughen the adherent surface by blasting with alumina powder (30-50 µm) using air pressure of 0.1-0.2 MPa (1-2 kgf/cm², 15-29 psi).
The air pressure and powder size should be properly adjusted to suit the material and/ or shape of the prosthetic restoration, using caution to prevent chipping.
- Clean the prosthetic restoration in an ultrasonic cleaning unit for 2 minutes.

[NOTE]

If the restorations are tried-in after acid treatment or blasting treatment, apply K-ETCHANT Syringe (phosphoric acid) for 5 seconds or KATANA Cleaner for 10 seconds with a rubbing motion then rinse with water until the color has been disappeared and dry the restorations to remove contaminants.

If the adherent surface is metal oxide ceramic (KATANA Zirconia STML/UTML/YML, KATANA Zirconia Block STML):

- Roughen the adherent surface by blasting with alumina powder (30-50 µm) using air pressure of 0.1-0.4 MPa (1-4 kgf/cm², 15-58 psi).
The air pressure and powder size should be properly adjusted to suit the material and/ or shape of the prosthetic restoration, using caution to prevent chipping.
- Clean the prosthetic restoration in an ultrasonic cleaning unit for 2 minutes followed by drying it with an air stream.

[NOTE]

If the restorations are tried-in after acid treatment or blasting treatment, apply KATANA Cleaner for 10 seconds with a rubbing motion then rinse with water until the color has been disappeared and dry the restorations to remove contaminants.

4. Application of CLEARFIL CERAMIC PRIMER PLUS to the prosthetic restoration

Apply the primer to the adherent surface of the restoration with an applicator brush. After application, dry the entire adherent surface sufficiently using mild oil free air flow.

5. Treatment of the prepared tooth

Choose either PANAVIA V5 Tooth Primer or the CLEARFIL Universal Bond Quick.

[NOTE]

KATANA Cleaner is available to clean the prepared tooth, when contamination occurs. Apply it for 10 seconds with rubbing motion, then rinse with water until the color has been disappeared and dry.

For PANAVIA V5 Tooth Primer

- If the adherent surface is uncut enamel or when cementing laminate veneers, apply phosphoric acid (K-ETCHANT Syringe) to the enamel surface only and leave it for 10 seconds; rinse and dry.
- Apply Tooth Primer to the entire prepared tooth surface, with an applicator brush and leave it in place for 20 seconds. Use caution not to allow saliva or exudates to contact the treated surfaces.
- Thoroughly dry the entire adherent surface sufficiently by blowing mild, oil-free air. Use a vacuum aspirator to prevent the Tooth Primer liquid from dispersing.

For CLEARFIL Universal Bond Quick

- Choose one of the three etching procedures before applying the bond. If the adherent surface includes uncut enamel or when cementing laminate veneers, application of phosphoric acid (K-ETCHANT Syringe) to the enamel surface is required< (1)-b or (1)-c >.
 - Self-etching procedure
Move to below (2) without etching with K-ETCHANT Syringe.
 - Selective enamel etching procedure
Apply K-ETCHANT Syringe to the uncut and/or cut enamel. Leave it in place for 10 seconds, then rinse and dry.
 - Total-etching procedure
Apply K-ETCHANT Syringe to the entire cavity (enamel and dentin), leave it in place for 10 seconds, then rinse and dry.
- Dispense the necessary amount of BOND into a well of the dispensing dish immediately before application.
- Apply BOND with a rubbing motion to the entire preparation with the applicator brush. No waiting time is required. Use caution not to allow saliva or exudate to contact the treated surfaces.
- Dry the entire cavity wall sufficiently by blowing mild air for more than 5 seconds until BOND does not move. Use a vacuum aspirator to prevent BOND from scattering.
- When cementing laminate veneers, there is no need to light-cure BOND prior to seating the laminate veneers. It should be cured after seating the laminate veneers with the PANAVIA Veneer LC Paste. When cementing inlays or onlays, BOND should be light-cured before seating the restoration. Please confirm the curing time by referencing the following table 1:

Table 1: Light-curing time of each light source

Type of light source (Light intensity)	Curing time
High-intensity BLUE LED* (More than 1500 mW/cm ²)	5 sec.
BLUE LED* (800-1400 mW/cm ²)	10 sec.
Halogen lamp (More than 400 mW/cm ²)	10 sec.

The effective wavelength range of each dental curing unit must be 400 - 515 nm.

* Peak of emission spectrum: 450 - 480 nm.

6. Syringe preparation of the cement paste

Select the suitable shade and remove the cap from the selected syringe and attach an applicator tip (16G) securely.

Cover the entire syringe with a disposable barrier (e.g. a poly bag) to prevent saliva and blood contamination.
Disinfect the syringe by wiping with an absorbent cotton with alcohol both before and after use.

7. Cementing the prosthetic restoration

- Apply the paste over the entire adherent surface of either the prosthetic restoration, or the prepared tooth.
- Place the prosthetic restoration on the prepared tooth.

8. Removing the excess cement

Remove any excess cement using either of the following two methods. Hold the restoration in place during excess removal.

Tack-Cure Clean-Up Technique:

Light-cure any excess cement in several spots for 1 second at a distance of 10-15 mm. Hold the prosthetic restoration in position, remove the tack-cured excess cement using a dental instrument. It is advisable to determine in advance the light-curing time of the excess cement by light-curing some paste on a mixing pad.

Wet Clean-Up Technique:

Any excess cement remaining at the margins can be removed with a small brush, foam pellet, dental floss or dental explorer.

[For curing margins of the restoration, it is beneficial to cover the margins with a protective gel (e.g. PANAVIA F 2.0 OXYGUARD II) to prevent the formation of an oxygen inhibition layer according to the Instructions for Use.]

9. Final curing

Light-cure the entire surface and margins of the prosthetic restoration that are less than 2.0 mm thick and have sufficient translucency. If the area you want to light-cure is wider than the light emitting tip diameter, extend the exposure process into a few applications.

Metal oxide ceramic (e.g. zirconia) should not be used because it may not have sufficient translucency. However, zirconia of our products with high translucency (KATANA Zirconia STML/UTML, KATANA Zirconia Block STML) can be used with a thickness of less than 1.2 mm. Please confirm the curing time by referencing the following table 2:

KATANA Zirconia YML can be used with a thickness of less than 1.2 mm only as laminate veneers.

Table 2: Light-curing time of each Restoration type

Restoration type and shade		Maximum restoration thickness	Type of curing light (Light intensity)		
			High-intensity BLUE LED (≥ 1500 mW/cm ²)	BLUE LED (1000~1400 mW/cm ²)	Halogen lamp (≥ 400 mW/cm ²)
Laminate veneer					
Silica-based ceramic (e.g. conventional porcelain, lithium disilicate)	BL shades, A1, A2, A3, A3.5, B1, B2	2.0 mm	Three times for 3 sec. or Twice for 5 sec.	10 sec.	20 sec.
	Other shades			20 sec.	
Hybrid ceramic, composite resin	All shades	1.2 mm			
KATANA Zirconia STML/UTML/YML, KATANA Zirconia Block STML					
Inlays, onlays					
Silica-based ceramic, hybrid ceramics, composite resin	All shades	2.0 mm	Three times for 3 sec. or Twice for 5 sec.	10 sec.	20 sec.
		1.2 mm			
KATANA Zirconia STML/UTML, KATANA Zirconia Block STML					
Shade and Depth of cure		Universal (A2), Clear, Brown (A4)			1.5 mm
		White			1.0 mm

10. Polishing the margins

Polish the margins using appropriate instruments for the polishing of composite resins.

[WARRANTY]

Kuraray Noritake Dental Inc. will replace any product that is proven to be defective. Kuraray Noritake Dental Inc. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

[NOTE]


If a serious incident attributable to this product occurs, report it to the manufacturer's authorized representative shown below and the regulatory authorities of the country in which the user/patient resides.

[NOTE]

"PANAVIA", "CLEARFIL" and "OXYGUARD" are registered trademarks or trademarks of KURARAY CO., LTD.

"KATANA" is a registered trademark or a trademark of NORITAKE. CO., LIMITED.

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