Noritake Super Porcelain EX-3

Rx Only

EN	ENGLISH	English	Dental Ceramic
FR	FRANÇAIS	French	Céramique Dentaire
ES	ESPAÑOL	Spanish	Cerámica Dental
IT	ITALIANO	Italian	Ceramica Dentale
NL	NEDERLANDS	Dutch	DENTAAL PORSELEIN
DE	DEUTSCH	German	Dentalkeramik
SV	SVENSKA	Swedish	Dentalkeramik
NO	NORSK	Norwegian	Dentalkeramikk
FI	SUOMI	Finnish	Hammaskeramiikka
DA	DANSK	Danish	Dental Keramik
PT	PORTUGUÊS	Portuguese	Cerâmica Dentária
EL	ΕΛΛΗΝΙΚΑ	Greek	Οδοντικό Κεραμικό Υλικό
TR	TÜRKÇE	Turkish	Dental Seramik
PL	POLSKI	Polish	Dentystyczny materiał ceramiczny
RO	ROMÂNĂ	Romanian	Ceramică dentară
HR	HRVATSKI	Croatian	Dentalna keramika
HU	MAGYAR	Hungarian	Fogászati kerámia
SL	SLOVENŠČINA	Slovenian	Dentalna keramika
CS	ČESKY	Czech	Dentální keramika
SK	SLOVENSKY	Slovak	Dentálna keramika
BG	БЪЛГАРСКИ	Bulgarian	Стоматологична керамика
UK	УКРАЇНСЬКА	Ukrainian	Стоматологічна кераміка
ET	EESTI	Estonian	Dentaalkeraamika
LV	LATVIEŠU	Latvian	Zobārstniecības keramika
LT	LIETUVIŠKAI	Lithuanian	Odontologinė keramika

UK REP	EN	UK Responsible Person	PL	Osoba odpowiedzialna w Wielkiej Brytanii
	FR	Personne responsable au Royaume-Uni	RO	Persoana responsabilă din Marea Britanie
	ES	Responsable en el Reino Unido	HR	Osoba odgovorna za Ujedinjeno Kraljevstvo (UK)
	IT	Responsabile per il Regno Unito	HU	Felelős személy az Egyesült Királyságban
	NL	Verantwoordelijke in het VK	SL	Odgovorna oseba v Združenem kraljestvu
	DE	Für UK verantwortliche Person	CS	Odpovědná osoba ve Spojeném království
	SV	Ansvarig person i Storbritannien	SK	Zodpovedná osoba v Spojenom kráľovstve
	NO	Ansvarlig person i Storbritannia	BG	Отговорно лице за Обединеното кралство
	FI	Vastuuhenkilö Yhdistyneessä kuningaskunnassa	UK	Відповідальна особа у Великій Британії
	DA	Ansvarlig person i Det Forenede Kongerige	ET	Suurbritannia vastutav isik
	PT	Responsável no Reino Unido	LV	Apvienotās Karalistes atbildīgā persona
	EL	Υπεύθυνο Πρόσωπο για το Ηνωμένο Βασίλειο	LT	Atsakingas asmuo Jungtinėje Karalystėje
	TR	Birleşik Krallık Sorumlusu		

Kuraray Noritake Dental Inc.

300 Higashiyama, Miyoshi-cho, Miyoshi, Aichi 470-0293, Japan
TEL +81-561-32-8953 FAX +81-561-32-8976

EMERGO EUROPE

Westervoortsedijk 60, 6827 AT Arnhem, The Netherlands

Kuraray Europe GmbH (EU Importer)

Philipp-Reis-Str. 4,
65795 Hattersheim am Main, Germany
Phone:+49 (0)69 305 35835 Fax:+49 (0)69 305 98 35835
URL:https://www.kuraraynoritake.eu

Noritake Super Porcelain EX-3

I. Introduction

This product is a dental porcelain, which has a potassium-aluminosilicate glass formulation as its primary chemical component. With Opagues available in both powder or paste, this product is suitable for a wide range of clinical applications, with variable shades. The general clinical benefit of this product is to restore tooth function for the following INDICATIONS FOR USE.

II. Indications for Use

For use in prosthetic dentistry to create a porcelain prosthesis on a frame of dental alloy

III. Composition

Paste type: Potassium-aluminosilicate glass (30-40 %), inorganic pigments (0-34 %), glycerol (10-30 %), butane-1,3-diol (1-8 %)

Powder type: Potassium-aluminosilicate glass (60-97 %), inorganic pigments (0-18 %) PASTE OPAQUE LIQUID: Glycerol (>99 %)

UP LIQUID: Glycerol (60-80 %), butane-1,3-diol (15-40 %), water (<1 %)

OPAQUE LIQUID: Water (70-90 %), polyethylene glycol (5-30 %), additives (<5 %)

FORMING LIQUID: Water (70-90 %), polyethylene glycol (5-30 %), additives (<5 %)

MEISTER LIQUID: Water (50-70 %), propane-1,2-diol (30-50 %), additives (<5 %)

IS LIQUID: 2-Phenoxyethanol (>99 %)

ES LIQUID: Propane-1,2-diol (>99 %)

ADDMATE FORMING LIQUID: Water (>99 %), oxidizing agents (<1 %)

Units in parentheses are mass %.

IV. Type and Class (ISO6872:2015)

Type:1 / Class:1

V. Physical Properties

See Table 2.

VI. Directions for Use

See Table 1 for the Baking schedule.

See the "TECHNICAL INSTRUCTIONS" for detailed information about the product's use, types and color variations.

① Preparation of metal framework

- Confirm the thickness of the metal framework and check the fit between the metal framework and the abutment by using the plaster model to evaluate the inner surfaces, margins, etc. Adjust the contours of the framework using a carbide bur, etc. so that the porcelain can be uniformly layered.
- After contouring, sandblast the facing surface of the framework according to the Instructions for Use supplied with the metal used. Clean the surface with an ultrasonic cleaner for approximately 10 minutes.
- Perform degassing according to the Instructions for Use that apply to the metal used.

② Opaque application and baking

[Paste Opaque]

- Apply a thin coat of a target shade of Paste Opaque, Paste Opaque Modifier or a mixture of them to the entire dry surface of the metal framework by rubbing. Apply the first layer of the paste thick enough so that approximately 70 percent of the metal coloration is masked, and bake it. If Universal Paste Opaque is used, apply to a minimum of 80 percent coverage and bake it according to the Baking schedule. For metal framework made of cobalt-chrome alloy, nickel-chrome alloy without Be and Noble alloys containing 75 - 85 percent palladium and copper, use NP Bonder of Paste Opaque. First, apply a thin coat of NP Bonder to the entire dry surface of the metal framework by rubbing. Then apply NP Bonder until approximately 70 percent of the metal coloration is masked and then bake it. (No color adjustment is available, in this case.)
- After baking, apply the second layer of the target shade of Paste Opaque, Paste Opaque Modifier, or a mixture. Apply until 100 percent of the metal coloration is masked and bake it again.

[Powder Opaque]

- Wet the surface of the metal framework slightly, using a brush moistened with water. Mix a target shade of Powder Opaque, Opaque Modifier or a mixture with OPAQUE LIQUID. Apply the first thin layer of the mixture over the entire surface of the metal framework, then bake it.
- After baking, apply the second layer (about 0.3 mm) of the mixture (prepared in the above step) until 100 percent of the metal coloration is masked, then bake it.

③ Internal Stain application and baking

- If the shade of a baked opaque needs to be adjusted, apply a mixture of Internal Stain and IS LIQUID, followed by additional baking. It is possible to use Internal Stain for the Cervical, Body, Enamel and Translucent building-up.

④ Cervical and Body building-up

- When fabricating the porcelain margin of the framework, first apply Margin on the marginal area.
- Cervical: Mix Cervical or the mixture of Body and Cervical with FORMING LIQUID or purified water. Apply the mixture on the cervical area of the framework.
- Crown contour: Mix the Body with FORMING LIQUID or purified water. Apply the mixture on the framework to shape a proper crown contour. Adjust the shade using Modifier and/or Mamelon, if necessary. If there is not enough space available to layer on the opaque, use Opacious Body.
- Cut the lingual, labial, and proximal surfaces back to make space for the building-up Enamel, then create the mamelon structure. Make sure that the Body porcelain remains at least 0.8 mm thick after it is cut back.

⑤ Enamel building-up

- Mix Enamel with FORMING LIQUID or purified water. Over-build the mixture one-third of the way from the incisal edge, so that the original contour can be restored. Be careful to avoid over-building the Enamel.
- Use Tissue to reproduce a gingival shade.

⑥ Translucent building-up and baking

- Mix Translucent and/or Luster (or Clear Cervical, if needed) with FORMING LIQUID or purified water. Over-build the mixture on the entire surface of the crown, approximately 10 percent thicker than the target contour of the crown. This will allow for the possibility of shrinkage due to baking. Then bake it.

⑦ Morphological Correction

- Using a diamond bur and a paper-abrasive material, adjust the crown and surface contours; clean the entire surface.

⑧ Glaze and Stain application, baking and finishing

- Reproduce the appropriate glaze using either method 1 or 2 below.
 - 1) Perform self-glaze baking.
 - 2) Mix External Stain Glaze with ES LIQUID, and apply a thin coat of the mixture. Then bake it.
- If the shade needs to be readjusted, apply a mixture of External Stain and ES LIQUID, and bake. After baking, polish and finish the restoration using silicon points, if needed.

⑨ Additional building-up

- Building-up the margins using Margin MRP, followed by re-baking, if there is not enough thickness. Building-up using Add-on, followed by re-baking if there is not enough contour of the crown.
- Addmate may be used for post-ceramic solder correction.

VII. Remarks on Handling

Contraindications:

- If the patient is hypersensitive to potassium-aluminosilicate glass or any other components, this product must not be used.

Caution:

1. If the patient or the dental professional demonstrates a hypersensitivity reaction, such as rash, dermatitis, etc., discontinue use of the product and seek medical attention immediately.
2. When mixing materials or grinding a restoration being fabricated, use an approved dust mask and vacuum with air filter to protect your lungs from inhaling the dust.
3. When mixing materials or grinding a restoration being fabricated, use safety glasses to prevent the dust from getting into your eyes. If the dust gets into your eyes, immediately rinse with copious amounts of water and consult a physician.
4. Do not use for any purposes except for fabricating dental restorations. This product is for dental application only.
5. Do not touch the items heated by the furnace with your bare hands.
6. PASTE OPAQUE, PASTE OPAQUE MODIFIER, PASTE OPAQUE LIQUID, UP LIQUID, IS LIQUID and ES LIQUID are flammable. Do not use them near fire or high-temperature objects.
7. IS LIQUID is a stimulant liquid. When using it, have good ventilation and if necessary wear an approved dust mask, safety glasses and gloves.
8. Stains and plaque can accumulate on a restoration, once it has been put in position in the patient's oral cavity, depending on the patient's eating habits or food choices. Instruct patients about the importance of cleaning their restorations.

Precautions in use:

1. For metal framework made of cobalt-chrome alloy, nickel-chrome alloy without Be and Noble alloys containing 75 - 85 percent palladium and copper, use NP Bonder of Paste Opaque for the first application.
2. When using a porcelain furnace for the first time, perform a baking test run in advance, since the working temperatures of porcelain furnaces may vary from one device to another.
3. Adjust the temperature and time according to the size of the restoration, referring to "Table 1. Baking schedule".
4. Use only the product with the LIQUIDS named in the Technical Instructions; otherwise, the restoration might change color when it is baked.
5. Do not use Internal Stain as an alternative to External Stain. Do not use External Stain as an alternative to Internal Stain.
6. Paste Opaque and Paste Opaque Modifier have liquid on the surface. Do not discard this liquid or stir the liquid into the paste. Tilt the container to reveal a portion of Paste Opaque where there is no liquid, and put this paste onto a pallet and blend to a smooth consistency.
7. Use caution when mixing two shades of Paste for color adjustment to avoid air bubbles from entering the mixture.
8. Paste Opaque will harden gradually over time after it has been dispensed onto the pallet. Use as soon as possible after dispensing. Do not use hardened Paste Opaque.
9. Do not mix Powder Opaque with Paste Opaque.
10. Do not mix this product with other porcelain materials or use it in unauthorized combinations.
11. Dispose of this product as a medical waste to prevent infection.

Storage:

1. Store in a cool and dry place. Keep away from direct sunlight.
2. The product should be stored at 1-30°C (33.8-86°F).
3. Place the cap securely back on the container after use.
4. The product must be stored in an appropriate place where only dental personnel have access.
5. The product must be used by the expiration date indicated on the package.
6. Do not store Paste Opaque, Paste Opaque Modifier, PASTE OPAQUE LIQUID, UP LIQUID, IS LIQUID and ES LIQUID near fire or high-temperature objects.

[WARRANTY]

Kuraray Noritake Dental Inc. will replace any product that is proven to be defective. Kuraray Noritake Dental Inc. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

[NOTE]

If a serious incident attributable to this product occurs, report it to the manufacturer's authorized representative shown below and the regulatory authorities of the country in which the user/patient resides. Refill does not include instructions for use (IFU). We provide the IFU in electronic format on the dedicated eIFU website.

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*1 96 kPa = 72 cmHg (29 inchesHg)

*2 This is the baking temperature for cases when Opacious Body, Body, Cervical, Clear Cervical, Enamel, Speed Enamel, Translucent, Luster, Mamelon, Modifier and Tissue are applied additionally.

*3 The baking conditions are established on the assumption that 1-3 units are baked.

In cases of 4-6 units, the dry-out time should be 10 min, release vacuum 925°C(1696°F), high temperature 935°C(1715°F).

In cases of over 7 units, the dry-out time should be 15 min, release vacuum 930°C(1706°F), high temperature 940°C(1724°F).

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