DUAL-CURED CORE BUILD-UP MATERIAL

CLEARFIL™ DC CORE PLUS KIT

DUAL-CURED CORE BUILD-UP MATERIAL

CLEARFIL™ DC CORE PLUS

ENGLISH

CLEARFIL™ DC CORE PLUS (PASTE)

[Precaution of the bonding system used in combination with CLEARFIL DC CORE PLUS]

Only CLEARFIL TRI-S BOND PLUS or dual-cured bonding agents can be used with CLEARFIL DC CORE PLUS. Do not use a light-cured bonding system except CLEARFIL TRI-S BOND PLUS. In the areas where the light cannot reach easily (e.g. adherent areas masked by the post), "CLEARFIL TRI-S BOND PLUS" can be self-cured in contact with CLEARFIL DC CORE PLUS.

[Precaution of ambient light]

Be careful to prevent unnecessary exposure to direct sunlight or surgery operating lights, otherwise the paste inside the tip may harden, leading to a shortened working time.

[Precaution of light curing]

Although CLEARFIL DC CORE PLUS is a dual-cured core build-up material, it must be light-cured in the following 3 procedures.

1. Application of BOND

Light-cure BOND with a dental curing unit.


Fig.1

2. Post-cementation using CLEARFIL DC CORE PLUS

Cure PASTE using the dental curing unit for the specified length of time, to fix the post in place. It is advisable to light-cure the margins well for strong bond.


Fig.2

3. Core build-up

After placing PASTE around the post, light-cure the paste from both the lingual and the labial (buccal) sides, using a dental curing unit.


Fig.3

CLEARFIL™ TRI-S BOND PLUS * (BOND)

- Use the light blocking plate to avoid exposing the material to an operating light or natural light, and use within 7 minutes after dispensing. As the volatile ethanol which is contained in the BOND evaporates, the viscosity increases, thereby making it difficult to apply.
- During the application of the BOND to the adherent surface and the duration of leaving it in place for 10 seconds after the application, move the lighting-spot out of the mouth or turn off the light to prevent the applied BOND from being exposed to the operating light because the BOND may gelate when left under such a light. In addition, after treatment for 10 seconds, dry sufficiently as soon as possible by blowing mild air.

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Kuraray Noritake Dental Inc.
Kuraray Europe GmbH