

WELCOME TO THE FUTURE OF CERAMIC LAYERING



DT Ioulianos Moustakis
Germany



MDT Andreas Chatzimpatzakis
Greece

Producing dental restorations that are not recognizable as such – this is probably the ultimate goal of every dental technician. For a long time, pursuing this goal was complicated by core materials whose optical properties were very different from those of natural teeth. The dark metal or opaque zirconia substructures had to be masked by applying multiple layers of intensively coloured ceramic powders, topped by more translucent porcelains imitating the enamel.

The rise of modern, tooth-coloured core materials such as lithium disilicate and zirconia has changed the game. With a core that is highly aesthetic, translucent and close to the final shade, it became much easier to produce a restoration that is virtually indistinguishable from the adjacent teeth. The thickness of the porcelain layer decreased as did the number of shades to be combined and necessary bakes to be conducted. The use of the existing porcelain systems for the new micro-layering techniques posed several new challenges: those systems originally developed for opaque zirconia were indicated for the more translucent zirconia core materials, but usually not for lithium disilicate. Moreover, the complexity of the systems made their use unnecessarily complicated for inexperienced users.

Consequently, Kuraray Noritake Dental Inc. developed a new porcelain system for micro-layering on zirconia and lithium disilicate core materials. The portfolio of CERABIEN™ MiLai, which refers to micro-layering and the Japanese word for future (mirai), consists of 15 internal stains (13 tooth colours including Bright to boost the translucent and Fluoro to boost the fluorescent effect, and two tissue colours) and 16 porcelains (12 tooth porcelains and four tissue porcelains). Hence, it enables dental technicians to implement a modernized version of the original Internal Live Stain Technique developed by Hitoshi Aoshima in the early 1990s in a porcelain layer of minimal thickness.

The following demo case is used to show how to achieve lifelike aesthetic restorations based on aesthetic zirconia and on lithium disilicate. Illustrating each step, the case allow users to anticipate how much time and effort can be saved compared to traditional layering techniques.

ADVANCED APPROACH WITH GUM AREAS ON ZIRCONIA

In this case, a highly complex ten-unit bridge with gum parts in the anterior region had to be produced. The selected framework material was KATANA™ Zirconia HTML Plus (Kuraray Noritake Dental Inc.), which offers a multi-layered colour structure, an optimized translucency and the high flexural strength required for long-span bridges. The restoration was milled in an anatomically reduced design and the surface texture was optimized with rotating instruments before sintering [Fig. 3a]. After the final sintering procedure, the restoration had a favourably high translucency in the incisal region and a natural shade structure [Figs. 3b and 3c]. In the first step of the micro-layering procedure, the application of the CERABIEN™ MiLai Internal Stains was planned

and carried out [Figs. 3d and 3e]. Subsequently, different layers of CERABIEN™ MiLai Porcelain were applied. The images 3f to 3h reveal which shades were combined and illustrate the procedure, while the outcome before and after the last bake is shown in Figures 3i to 3k. In the next step, the gum areas were completed using the CERABIEN™ MiLai tissue porcelains Tissue 4, 5 and 6 in the order and locations described in Figures 3l to 3o. In the final layer, Tissue 1 was mixed with ELT1 to imitate the labial frenulum and with LTx to create a smooth transition to the natural gingiva [Figs. 3p and 3q]. The final restoration is shown in Figure 3r.



Fig. 3a: Milled restoration after surface texturing.



Fig. 3b: Shade and translucency of the sintered zirconia restoration.




Fig. 3c: Highly translucent bridge on the model.




Fig. 3d: Chroma map for the application of CERABIEN™ MiLai Internal Stains.

Fig. 3e: Applying a mixture of Bright, Salmon Pink and Tissue Pink to the gum area.






Fig. 3f: Application of CERABIEN™ MiLai E2 to add translucency to the structure.




Fig. 3g: Application of Tx and a mixture of Tx and CCV-2 to individualize the cervical and incisal areas while boosting the translucency of the enamel in the middle and incisal third.




Fig. 3h: Adding a final layer of LT1 for additional translucency and opalescence.



Fig. 3i: Appearance of the ten-unit bridge before the bake – labial view.



Fig. 3j: Appearance of the ten-unit bridge before the bake – palatal view.

Fig. 3k: Appearance of the ten-unit bridge after the bake.




Fig. 3l: Application of small amounts of Tissue 5 ...




Fig. 3m: ... covered with a layer of Tissue 6 alternating with Tissue 5.



Fig. 3n: Following another bake, Tissue 5 is applied in the proximal areas.




Fig. 3o: How to combine Tissue 6 and Tissue 4 in the next layer.




Fig. 3p: How to complete the tissue layer with Tissue 1, locally mixed with ELT1 or LTx.



Fig. 3q: Restoration before the final bake.



Fig. 3r: Final ten-unit bridge ready for placement.



IOULIANOS MOUSTAKIS

Dental Technician/Photographer

1985 - 1987 Studied at the School of Dental Technology (SBIE) in Athens / Greece
 1997 - 1998 Master school in Berlin
 2007 - Education as Maxillofacial prosthetic technician (IASPE)
 2010 - Advanced education in Functional diagnosis temporomandibular joint
 2011 - 2012 Curriculum implant prosthetics for dental technicians (DGZI)
 2013 - 2014 Education as a graphic designer at the Media Design Hochschule (MDH) in Berlin
 2015 - 2017 Education as a photographer at the Photocentrum of the Gilberto Bosques VHS Friedrichshain-Kreuzberg
 2014 - 2016 - 2018 Further training at Noritake/Kuraray in Japan
 2019 - International Instructor Noritake/Kuraray
 2019 - Certified trainer of Teeth Morphology carving (Osaka Ceramic Training Center)
 1998 - Implant Dental Studio - Athens/Greece
 2010 - Zirkler & Moustakis Dental Technology - Falkensee/Germany
 2020 - Giuliano Dental design - Falkensee/Germany

Publications in Dental Journals

2014 - 5/2014 Dental Dialogue/Germany
 2015 - 10/2015 The International Journal of Dental Technology/Japan
 2018 - 1/2018 Cosmetic Dentistry/Germany
 2018 - 4/2018 Zahntechnik Zeitung/Germany
 2018 - 5/2018 Das Dental Labor/Germany
 2018 - 5/2018 Dental Dialogue/Italy
 2018 - 10/2018 Laborama/Greece
 2019 - 1/2019 LabLine/Hungary
 2019 - 3-4/2019 Dental Technologies/UK
 2020 - 4/2020 LabLine/Hungary
 2021 - 1+2 LabLine/Hungary
 2021 - 5/2021 + 12/2021 Quintessenz Zahntechnik/Germany
 2021 - 4/2021 QDRP France

Competitions

2013 – 6th place at the 8th KunstZahnWerk contest by Candulor
 2017 – 5th place at the 10th KunstZahnWerk contest by Candulor
 2017 – 1st place at the 10th KunstZahnWerk contest by Candulor as "Best Documentation"
 2020 – 1st place at the 4th Panthera Master Cup by Panthera Dental

Memberships

NGSC Noritake Greek Study Club
 DGZI German Society of Dental Implantology
 IASPE International Association for Surgical Prosthetics and Epithetics
 Key Opinion Leader (KOL) at company MPF Brush Company
 Key Opinion Leader (KOL) at company Candulor
 Key Opinion Leader (KOL) at company Kuraray/Noritake



ANDREAS CHATZIMPATZAKIS

Andreas graduated from the Dental Technology Institute (TEI) of Athens in 1999. During his studies he followed a program at the Helsinki Polytechnic Department of Dental Technique, where he trained on implant superstructures and all ceramic prosthetic restorations.

From the year 2000, he is running the ACH Dental Laboratory in Athens, Greece, specialized on refractory veneers, zirconia and long span implant prosthesis.

ACH Dental Laboratory is Co-operating lab with the National and Kapodistrian University of Athens for the MSc degree in Dental Laboratory Materials.

From 2016 he is key opinion leader for the MPF Brush.co.

On 2017 he visits Japan where he trained from Hitoshi Aoshima, Naoto Yuasa and Kazunabu Yamanda and becomes International Trainer for Kuraray – Noritake company.

In 2018 he became Editor-in-chief for the dental technician magazine "LABORAMA" published by OMNIPRESS co.

On 2019 he studies carving, morphology and all ceramic restorations at the Osaka Ceramic Training Center by Shigeo Kataoka.

On 2019 he establishes the Dental Technicians' Coaching Services and coaches dental technicians to improve their work.

Andreas has also conducted several lectures and hands on seminars in Greece and abroad and published articles in Greek and international magazines.

His lecture "An exciting journey ... to be a dental technician" is about inspiring dental technicians to improve their work by observing and emulate natural teeth using the internal live stain technique.