



SELF-ADHESIVE RESIN CEMENT

PANAVIA SA Cement Universal

- ✓ Suitable for prosthetic restorations, even LiDiSi cementation
- ✓ Easy, gum-friendly excess removal
- ✓ No need for separate primers
- Storage at room temperature

PANAVIA SA Cement Universal

PANAVIA™ SA Cement Universal offers a wide range of indications for adhesion and is suitable for crowns/bridges, inlays/onlays, posts, and even adhesion bridges. Thanks to its powerful singlestep bonding formula, cementation is easier than ever before.

CLEARFIL™ UNIVERSAL BOND QUICK

CLEARFIL $^{\mathbb{T}}$ Universal Bond Quick can be used in combination with PANAVIA $^{\mathbb{T}}$ SA Cement Universal to enhance bond strength to tooth structure.

A BREAKTHROUGH FROM THE INNOVATORS OF ADHESION AND SILANE TECHNOLOGY.



The unique LCSi monomer creates a strong chemical bond to Porcelain, Lithium Disilicate and Composite Resin.



The Original MDP monomer¹ creates a strong chemical bond to Enamel, Dentin, Metal Alloy and Zirconia.

PANAVIA[™] SA Cement Universal is a forerunner in self-adhesive resin cements containing the unique LCSi monomer that adheres to virtually every material including glass ceramics without the need for a separate primer. Using a unique silane coupling agent built into the paste - LCSi monomer - the cement delivers a strong, durable chemical bond to porcelain, lithium disilicate and composite resin. The original MDP monomer present in the paste too, allows for chemical reactiveness with zirconia, dentin and enamel.

ONE SINGLE PROCEDURE FOR ALL ADHESIVE CEMENTATIONS

Example: single crown

Regardless of the restorative material; no primer required. One single procedure regardless of the restorative material.



Conditioning the crown.*

*Clean and dry the tooth surface, and then trial fit the prosthetic restoration, followed by conditioning of the crown.



Dispense an equal amount of Paste A & B. Mix Paste A & B for 10 seconds.



Apply the cement to the crown.



Place the crown.



Light-cure for 2 to 5 seconds or let self-cure for 2 to 4 minutes, then remove the excess cement.



Light-cure the margin and maintain isolation for 5 minutes.*

*For a translucent restoration, fully light-cure.

SHADES



Universal (A2)



White



Translucent

Before using this product, be sure to read the Instructions for Use supplied with the product. The specifications and appearance of the product are subject to change without notice. Printed color can be slightly different from actual color. "PANAVIA" and "CLEARFIL" are trademarks of Kuraray Co., Ltd. "KATANA" is a trademark of NORITAKE CO., LIMITED.

