INDICATIONS

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Sealing of a prepared abutment tooth as a pretreatment for indirect restorations

Follow the standard procedure for isolation, moisture control, preparation of abutment tooth

- 1 Tooth pretreatment**
 Choose either etching procedure
 - a. Self-etching (Move to section 2)
 - b. Selective-etching*¹ **10sec.**
 - c. Total-etching*¹ **10sec.**
- 2 Apply BOND with a rubbing motion**
 No waiting time
- 3 Dry by blowing mild air**
 until BOND does not move*² **5sec. +**
- 4 Light-cure*³**
- 5** If necessary, place a thin coat of composite resin (e.g. CLEARFIL MAJESTY ES Flow) onto the tooth, light-cure according to the manufacturer's instructions.
- 6 Wipe the surface to remove the un-polymerized layer (oxygen inhibited layer)*⁴**

*¹ Refer to **Direct Restoration**
 *² Use a vacuum aspirator to prevent BOND from scattering.
 *³ Refer to Table 1 for light-curing time.
 *⁴ Use a cotton pellet or a gauze moistened with alcohol.

Intraoral repair of fractured restorations

- 1 Roughen, rinse and air dry**
- 2 Apply K-ETCHANT Syringe, then rinse and dry*¹** **5sec.**
- 3 Apply BOND with a rubbing motion**
 No waiting time
- 4 Dry by blowing mild air**
 until BOND does not move*² **5sec. +**
- 5 Light-cure*³**
- 6 Place composite resin*⁴, light-cure and finish**

*¹ This acid etching is not necessary for non-precious metal and metal oxide ceramic.
 *² Use a vacuum aspirator to prevent BOND from scattering.
 *³ Refer to Table 1 for light-curing time.
 *⁴ Use an opaque resin (e.g. CLEARFIL ST OPAQUER) to mask metal color.

Post cementation / Core build-ups with CLEARFIL DC CORE PLUS

Follow the standard procedures for isolation, moisture control and preparation of root canal and cavity

- 1 Post pretreatment**
 For glass fiber post
 [1] Apply K-ETCHANT Syringe, rinse and dry **5sec.**
 [2] Apply BOND, then dry by blowing mild air **5sec. Dry**
 For metal post
 [1] Blast with alumina powder, then ultrasonic clean and dry
- 2 Tooth pretreatment**
 Choose either etching procedure
 - a. Self-etching (Move to section 3)
 - b. Selective-etching*¹ **10sec.**
 - c. Total-etching*¹ **10sec.**
- 3 Apply BOND with a rubbing motion**
 No waiting time
- 4 Dry by blowing mild air and paper point until BOND does not move*²** **5sec. +**
- 5 Light-cure*³**
- 6 Post cementation and core build-up using CLEARFIL DC CORE PLUS according to the manufacturer's instructions**

*¹ Refer to **Direct Restoration**
 *² Use a vacuum aspirator to prevent BOND from scattering.
 *³ Refer to Table 1 for light-curing time.

Post cementation / Core build-ups with other core material (except for CLEARFIL DC CORE PLUS)

Follow the standard procedures for isolation, moisture control and preparation of root canal

1 Post pretreatment

For glass fiber post

[1] Apply K-ETCHANT Syringe, rinse and dry **5sec.**

[2] Apply the mixture of BOND and CLEARFIL DC Activator*1, then dry by blowing mild air **5sec. Dry**

[3] Light-cure*2 **Note**

For metal post

[1] Blast with alumina powder, then ultrasonic clean and dry **5sec.**

2 Tooth pretreatment

Choose either etching procedure

a. Self-etching (Move to section 3)

b. Selective-etching*3 **10sec.**

c. Total-etching*3 **10sec.**

3 Apply the mixture*1 with a rubbing motion

No waiting time

4 Dry by blowing mild air and paper point until the mixture does not move*4

5sec. +

5 Light-cure*2

Note

Note: Working time will be dramatically shortened when not light-curing
*4 Use a vacuum aspirator to prevent the mixture from scattering.

6 Post cementation and build-up by core material according to the manufacturer's instructions

*1 Dispense one drop each of BOND and CLEARFIL DC Activator and mix them.
*2 Refer to Table 1 for light-curing time.

*3 Refer to Direct Restoration

Cementation of indirect restorations with Kuraray's self-adhesive cements

Clean and dry the tooth surface, and then trial fit the prosthetic restoration

1 Surface preparation of prosthetic restorations

Silica-based glass ceramic (e.g. lithium disilicate)

Apply a hydrofluoric acid, then wash and dry

Metal oxide ceramics (e.g. zirconia), metals or composite resins

Blast with alumina powder (30-50µm/ 0.1-0.4MPa/ 14-58 PSI/ 1-4 bar), then ultrasonic clean and dry

2 Tooth Pretreatment

Choose either etching procedure

a. Self-etching (Move to section 3)

b. Selective-etching*1 **10sec.**

c. Total-etching*1 **10sec.**

3 Apply BOND with a rubbing motion

No waiting time

4 Dry by blowing mild air until BOND does not move*2

5sec. +

*2 Use a vacuum aspirator to prevent BOND from scattering.

5 Cementation using Kuraray's self-adhesive cements according to the manufacturer's instructions

Note: When using a partial light-curing (or "Tack-Cure") technique, the setting time of the excess cement will be shorter.

*1 Refer to Direct Restoration

Cementation of indirect restorations with self-adhesive resin cement without any specific instructions to pretreatment the adherent surface

Clean and dry the tooth surface, and then trial fit the prosthetic restoration

1 Surface preparation of prosthetic restorations

Silica-based glass ceramic (e.g. lithium disilicate)

Apply a hydrofluoric acid, then wash and dry

Metal oxide ceramics (e.g. zirconia), metals or composite resins

Blast with alumina powder (30-50µm/ 0.1-0.4MPa/ 14-58 PSI/ 1-4 bar), then ultrasonic clean and dry

2 Apply the mixture of BOND and CLEARFIL DC Activator*1, then dry until the mixture does not move*2

5sec. + Dry

Note

Note: Working time will be dramatically shortened when not light-curing

*1 Dispense one drop each of BOND and CLEARFIL DC Activator and mix them.
*2 Use a vacuum aspirator to prevent the mixture from scattering.
*3 Refer to Table 1 for light-curing time.

3 Light-cure*3

Note

4 Tooth Pretreatment

Choose either etching procedure

a. Self-etching (Move to section 5)

b. Selective-etching*4 **10sec.**

c. Total-etching*4 **10sec.**

5 Apply BOND with a rubbing motion

No waiting time

6 Dry by blowing mild air until BOND does not move*5

5sec. +

*5 Use a vacuum aspirator to prevent BOND from scattering.

7 Light-cure*3

Note

8 Cementation using resin cement according to the manufacturer's instructions

Note: When using a partial light-curing (or "Tack-Cure") technique, the setting time of the excess cement will be shorter.

*4 Refer to Direct Restoration