EN	ENGLISH	English	Dental CAD/CAM material
FR	FRANÇAIS	French	Consommables pour CFAO (CAD/CAM)
ES	ESPAÑOL	Spanish	CAD CAM CONSUMIBLE
IT	ITALIANO	Italian	Materiale dentale CAD/CAM
NL	NEDERLANDS	Dutch	Tandheelkundig CAD/CAM-materiaal
DE	DEUTSCH	German	CAD/CAM-Dentalmaterial
SV	SVENSKA	Swedish	Dentalt CAD/CAM-material or CAD/CAM block
NO	NORSK	Norwegian	CAD/CAM-dentalmateriale
FI	SUOMI	Finnish	Hammaslääketieteellinen CAD/CAM-materiaali
DA	DANSK	Danish	CAD/CAM-dentalmateriale
PT	PORTUGUÊS	Portuguese	CAD/CAM CONSUMIVEL
EL	ΕΛΛΗΝΙΚΑ	Greek	Οδοντιατρικό υλικό CAD/CAM
TR	TÜRKÇE	Turkish	Dental CAD/CAM materyali
PL	POLSKI	Polish	Materiał stomatologiczny do CAD/CAM
RO	ROMÂNĂ	Romanian	Material dentar CAD/CAM
HR	HRVATSKI	Croatian	Stomatološki CAD/CAM materijal
HU	MAGYAR	Hungarian	Fogászati CAD/CAM anyag
SL	SLOVENŠČINA	Slovenian	Dentalni CAD/CAM-material
CS	ČESKY	Czech	Dentální materiál CAD/CAM
SK	SLOVENSKY	Slovak	Dentálny materiál CAD/CAM
BG	БЪЛГАРСКИ	Bulgarian	Стоматологичен CAD/CAM материал
UK	УКРАЇНСЬКА	Ukrainian	Стоматологічний матеріал для CAD/CAM
ET	EESTI	Estonian	CAD/CAM-hambaravimaterjal
LV	LATVIEŠU	Latvian	CAD/CAM zobārstniecības materiāli
LT	LIETUVIŠKAI	Lithuanian	Odontologinė CAD/CAM medžiaga

Kuraray Noritake Dental Inc.

1621 Sakazu, Kurashiki, Okayama 710-0801, Japan


EMERGO EUROPE

Westervoortsedijk 60,
6827 AT Arnhem, The Netherlands

Kuraray Europe GmbH (EU Importer)

Philipp-Reis-Str. 4,
65795 Hattersheim am Main, Germany
Phone: +49 (0)69 305 35835 Fax: +49 (0)69 305 98 35835
URL: <https://www.kuraraynoritake.eu>

I. INTRODUCTION

KATANA AVENCIA Block 2 is CAD/CAM blocks of composite resin material that is milled to fabricate dental restorations, including inlays and crowns. It has 2 types: Mono Layered block and Multi Layered block. The general clinical benefit of this product is to restore tooth function in the following INDICATIONS FOR USE.

II. INDICATIONS FOR USE

KATANA AVENCIA Block 2 is indicated for the following uses:

Fabrication of inlays, onlays, veneers and full contour crown restorations using a dental CAD/CAM system.

III. CONTRAINDICATIONS

- [1] Patients with a history of hypersensitivity to methacrylate monomers
- [2] This product should NOT be used when malocclusion, clenching or bruxism conditions are present.

IV. PRECAUTIONS

1. Safety precautions

1. This product contains substances that may cause allergic reactions. Avoid use of the product in patients with known allergies to methacrylate monomers or any other components.
2. If the patient demonstrates a hypersensitivity reaction, such as rash, eczema, features of inflammation, ulcer, swelling, itching or numbness, discontinue use of the product and seek medical attention.
3. Avoid direct contact with the skin and/ or soft tissue to prevent hypersensitivity. Wear gloves or take appropriate precautions when using the product.
4. Single use only. Do not reuse the product.
5. Exercise caution to prevent the patient from accidentally swallowing the product.
6. If the instruments of this product are damaged, protect yourself from any danger and immediately discontinue their use.
7. Do not use the product if there is a crack noticed after removing it from the package. Do not use if you find a crack in the final prosthesis.
8. Dispose of this product as a medical waste to prevent infection.

2. Handling and manipulation precautions

[Precautions when using the product]

1. Multi Layered block consists of a Body layer, two transition layers and an Enamel layer. The layer on the same side as the notch in the metallic holder is the Body; the layer on the opposite side is the Enamel. Make sure that Multi Layered block can be designed and milled in the right direction (Body side and Enamel side) before using. It may be milled in the reverse direction depending on the CAD/CAM machine setting.
2. The product must not be used for any purposes other than specified in [II.INDICATIONS FOR USE].
3. The use of this product is restricted to dental professionals.
4. Use safety equipment, such as a local vacuum unit, protective shields, and an appropriately approved dustproof mask, to prevent the operator from inhaling dust when shaping, finishing or polishing the product.
5. For bonding, use a resin cement (e.g. PANAVIA V5 or PANAVIA SA Cement Universal). Before bonding, sandblast the adherent surface and then condition it with an etching agent (e.g. K-ETCHANT Syringe) and with a silane coupling agent (e.g. CLEARFIL CERAMIC PRIMER PLUS), according to the Instructions for Use.
6. Before applying additional resin cement, characterizing or glazing the surface, treat the adherent surface by etching it with an etching agent and then condition it with a silane coupling agent, according to the Instructions for Use.
7. Shape the outlines of inlays or onlays with appropriate smooth curves and locate the margins so that there is no contact with the opposing teeth.
8. The margins of crowns should be prepared with a deep chamfer. The margins of inlays/onlays should be prepared with a butt joint, while chamber margins should be prepared with a deep chamfer.
9. The unilateral bevel of the axial surface and the adjacent surface of the crowns should be set at 6°-10°.
10. When preparing the abutment of the crown, avoid the following:
 - Deep shoulder or "J" margin
 - Knife Edge
 - Parallel Axis
 - Groove
 - Rough Margin
 - Undercut
 - Sharp Incisal Top
11. When preparing the cavity of the inlay/onlay, avoid the following:
 - Bevel
 - Complex external shape
 - Narrow isthmus
 - Sharp line angle and point angle
 - Undercut
 - Parallel cavity sidewall
12. When using this product to fabricate a restoration, be sure to attain the thicknesses listed in the following table:

Restoration		Thickness
Anterior Crown	Incisal	>1.5 mm
	Circumferential wall	1.0-1.2 mm
	Margin	>0.8 mm
Posterior Crown	Fissure / Cusps / Circumferential wall	>1.5 mm
	Margin	>1.0 mm
Inlay	Fissure	>1.0 mm
	Isthmus	>1.5 mm
Onlay	Fissure	>1.0 mm
	Cusps / Isthmus	>1.5 mm
Veneer	Incisal	0.5-1.5 mm
	Labial wall	>0.6 mm
	Margin	>0.4 mm

13. When using this product to fabricate a restoration, the width of isthmus should be more than 1.5 mm.
14. This product is a CAD/CAM composite resin material. When milling the product, use the appropriate CAD/CAM milling machine and follow the Instructions for Use.
15. Do not bake this product in a porcelain furnace or similar item.

3. Storage precautions

1. The product must be used by the expiration date indicated on the package.
2. The product must be stored at 2-25°C/ 36-77°F when not in use. Store in a cool and dry place. Keep away from direct sunlight.
3. The product must be stored in a proper place where only dental practitioners can access.

V. COMPONENTS

1. Shades

KATANA AVENCIA Block 2 is available in the following 4 shades:

Mono Layered block: A1 HT, A2 HT, A3 HT, A3.5 HT
Multi Layered block: A1 ML, A2 ML, A3 ML, A3.5 ML

2. Sizes

KATANA AVENCIA Block 2 is available in the following 2 sizes:

Size 12 (10×12×15 mm)
Size 14L (14.5×14.5×18 mm)

3. Metallic holder

KATANA AVENCIA Block 2 is 2 kinds of metallic holder available: Universal and CEREC.

4. Components

See outer package for the items included in the package.

5. Ingredients

Principal ingredients:

- Mixed filler with colloidal silica and aluminum oxide (50-80%)
- Cured resins consisting of methacrylate monomer (Copolymer of Urethane dimethacrylate (10-30%) and other methacrylate monomers (5-20%))
- Pigments (< 1%)

Units in parentheses are mass %.

VI. CLINICAL PROCEDURES

1. Preparation

Prepare the restoration area in the usual manner.

2. Designing and processing the prosthesis

Select the product shade and size that suits the case. Design and process it using a CAD/CAM system. See the Instructions for Use.

3. Shaping and polishing

Shape the prosthesis in the usual manner and polish it with silicon points and other available methods.

4. Sandblasting

Roughen the adherend surface by blasting with alumina powder (30-50 μm) using air pressure of 0.1-0.2 MPa (1-2 kgf/cm², 14-29 psi).

5. Bonding

After a trial fit, treat the adherent surface with an etching agent and a silane coupling agent according to the Instructions for Use. Use a resin cement to bond the prosthesis.

[WARRANTY]


Kuraray Noritake Dental Inc. will replace any product that is proven to be defective. Kuraray Noritake Dental Inc. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

[NOTE]

- If a serious incident attributable to this product occurs, report it to the manufacturer through the EU Importer, and to the regulatory authorities of the country in which the user/patient resides.
- The Summary of Safety and Clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed>. Basic UDI-DI: 458211076BlockMH

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"KATANA" is a registered trademark or trademark of NORITAKE CO., LIMITED.
"AVENCIA", "PANAVIA" and "CLEARFIL" are registered trademarks or trademarks of KURARAY CO., LTD.
„Kuraray Noritake Dental Inc.“ pakeičia visus produktus su įrodymais trūkumais.