

# TEETHMATE F-1

## KIT COMPONENTS

### TEETHMATE F-1 Resin, 2 x 2,5 ml

Light cured pit and fissure sealant containing adhesion monomer (MDP) and a fluoride releasing polymer. Available in 3 different shades: natural, opaque and red.



### K-ETCHANT GEL, 6 ml

Thixotropic etching agent containing 40 % phosphoric acid



## ACCESSORIES

Brush tip handle, 50 disposable brush tips, 2 x 25 applicator nozzles, plate



## RELATED PRODUCTS

### CARIES DETECTOR

Caries Detector is a clinical aid in the diagnosis of infected carious dentine. Caries Detector penetrates the broken cross links in the collagen fibres, which have been destroyed by bacteria, staining this area in a red colour. Caries Detector is also a conservative aid, for it will not penetrate the inner layer – which is partly decalcified and will remineralize if left intact. Thus the doubt over where the carious dentine begins and ends, is no longer so difficult to establish when Caries Detector is used routinely in Conservative Dentistry.

### FEATURES AND BENEFITS

- Fast and easy to use
- Does not stain remineralizable carious-affected dentin
- Helps clinicians to preserve healthy tissue



## CONVENIENT FLOW CONTROL



TEETHMATE F-1 is a fluoride-releasing,

TEETHMATE F-1

light-curing, strongly adhering pit and

fissure sealant in a newly developed

application system for easy "flow-control".

The stable co-polymer releasing fluoride

is a unique invention by Kuraray.



### Literature:

- Charbeneau, C.T., Dennison, J.B., Ryge, G.: A filled pit and fissure sealant: 18-month results. J Am Dent Assoc Aug;55(2):299-306 (1977)
- Garcia-Godoy, F., Abarzua, I., De Goes, M.F., Chan, D.C.: Fluoride release from fissure sealants. J Clin Pediatr Dent 22(1):45-49 (1997)
- Iritoda, Y., Matsumura, Y., Kito, H., Nakano, T., Toyama, T., Nakagaki, H., Tsuchiya, T.: Effect of sealant viscosity on the penetration of resin into etched human enamel. Oper Dent Jul-Aug;25(4): 274-82 (2000)
- Kula, K., Thompson, V., Kula, T., Nelson, S., Selvaagi, R., Liao, R.: In vitro effect of topical fluorides on sealant materials. J Esthet Dent Jul-Aug 4(4):121-7 (1992)
- Mizuno, Y., Ohnori, I.: An in vitro study on the fluoride releasing resin sealant. Pediatric Dental Journal, 1(1): 89-93 (1991)
- Rock, W.P., Weatherill, S., Anderson, R.J.: Retention of three fissure sealant resins. The effects of etching curing method. Results over 3 years. Br Dent J Apr 21;168(8):323-5
- Shinji, H., Uchimura, N., Ishida, M., Motokawa, W., Miyazaki, K., Garcia-Godoy, F.: Enamel fluoride uptake from a fluoride releasing sealant. Am J Dent 11:58-60 (1998)
- Stach, D.J., Hatch, R.A., Tilliss, T.S., Cross-Poline, G.N.: Change in occlusal height resulting from placement of pit and fissure sealants. J Prosthet Dent Nov; 68(5):750-3 (1992)



Manufactured by KURARAY CO., LTD. 1-12-39 Umeda, Kita-ku, Osaka 530-8611, JAPAN  
Kuraray European Union Representative KURARAY EUROPE GmbH, Schiess Strasse 68, D-40549 Dusseldorf, eMail: dental@kuraray.de

© 2007 Kuraray Co., Ltd. All rights reserved. MDP, MONOMER, LONG TERM DURABILITY, HIGH BOND STRENGTH, and TEETHMATE F-1 are trademarks of Kuraray Co., Ltd. in Japan and other countries. REALITY RESEARCHING CO. is an independent research institute. Caries Detector is a registered trademark of Kuraray Co., Ltd. in Japan and other countries. K-ETCHANT GEL is a registered trademark of Kuraray Co., Ltd. in Japan and other countries. The logo and design of the product are the property of Kuraray Co., Ltd. in Japan and other countries. The logo and design of the product are the property of Kuraray Co., Ltd. in Japan and other countries. The logo and design of the product are the property of Kuraray Co., Ltd. in Japan and other countries.

# TEETHMATE F-1

