



## COLOR COMBINATION TABLE

### PASTE OPAQUE 6g / 2.4ml



### PASTE OPAQUE MODIFIER 3g / 1.2ml




### UNIVERSAL PASTE OPAQUE 6g




### UNIVERSAL PASTE OPAQUE MODIFIER 3g




### POWDER OPAQUE 10g / 50g / 200g




### POWDER OPAQUE MODIFIER 10g / 50g




### MARGIN 10g / 50g



### OPACIOUS BODY 10g / 50g / 200g



### CERVICAL 10g / 50g / 200g



## NORITAKE SUPER PORCELAIN EX-3 LINE-UPS

**Universal Paste Opaque:** Same color as Paste Opaque and Paste Opaque Modifier.


**Value Shade:** Value Shades are additional shades that expand the existing Noritake Super Porcelain EX-3 selection of body and Opacious Body porcelain powder to include the VITA® 3D-Master® shade.

**Speed Enamel (S1, S2, S3, S4):** Specially developed for the two layer build-up technique.


**n color Shade:** n color shade is formulated to achieve chroma intensified Body and Paste Opaque (Universal Paste Opaque, Powder Opaque) of the shade. With intensified chroma, it even can be used in a case where there is insufficient space for porcelain build-up.

## COLOR COMBINATION TABLE


### BODY 10g / 50g / 200g



### ENAMEL 10g / 50g / 200g



### SPEED ENAMEL 10g / 50g / 200g




### TRANSLUCENT 10g / 50g / 200g



### LUSTER 10g / 50g / 200g



### CLEAR CERVICAL 10g / 50g / 200g




### MAMELON 10g / 50g



### TISSUE 10g / 50g



### MODIFIER 10g / 50g



### ADD-ON 10g / 50g



### INTERNAL STAIN 3g



### EXTERNAL STAIN 3g (Glaze 10g / 30g)



### ADDMATE 10g

