

*Super
Porcelain* **EX-3**

Porcelain for metal ceramic dental restorations

Noritake Super Porcelain EX-3
For True-to-life result



MIMIC NATURE

A TECHNOLOGICAL MASTERPIECE THAT REFLECTS THE NATURALNESS OF DENTITION THROUGH YOUR HANDS.


Noritake Super Porcelain EX-3 is a high-quality porcelain system that was developed as a result of the extensive knowledge and ceramic technology of Kuraray Noritake Dental.

PRECISION OF THE NATURAL DENTITION

Kuraray Noritake Dental's shade concept creates an intra-oral harmony by diminishing the greenish color and adding a pinkish hue.

FLUORESCENCE CREATES THE NATURAL LOOK

Noritake Super Porcelain EX-3 features a clean fluorescent pale blue which creates the look of a natural tooth.



Noritake Super Porcelain EX-3 / Natural Tooth under UV Light

CLINICAL DATA

Survival Analysis of PFM (with Noritake Super Porcelain EX-3) posterior crowns

Group	Total # Crowns	Mean Survival (SE) (Days)	Probability of Survival in %
PFM	1080	2694 (5)	99.3

Fusun Ozer, DMD, PhD others, "A retrospective survey on long-term survival of posterior zirconia and porcelain-fused-to-metal crowns in private practice", Quintessence Int 2014;45:31-28

Coping system: Lava, 3M ESPE; Procera, Nobel Biocare; KATANA™, Kuraray Noritake Dental Fabricated by one dental laboratory (Cusp, Boston)

CASE 1

Porcelain Fused to Metal Restoration



Tooth preparation



Postoperative view

Fabrication by Naoki Hayashi, RDT

CASE 2

Porcelain Laminate Veneer Restoration



Preoperative view



Postoperative view


Fabrication by Nondas Vlachopoulos, RDT

CASE 3

PFM Incisor Crowns using Noritake Super Porcelain EX-3



Preoperative view




1 PFM. 3 Laminates on refractory



Postoperative view

Fabrication by Daniele Rondoni, RDT


Noritake Super Porcelain EX-3 chromatic map



by Daniele Rondoni, RDT


COLOR COMBINATION TABLE

PASTE OPAQUE 6g / 2.4ml




*When the different types of alloy is used, NP Bonder color will result in different color outcome.


PASTE OPAQUE MODIFIER 3g / 1.2ml




UNIVERSAL PASTE OPAQUE 6g




UNIVERSAL PASTE OPAQUE MODIFIER 3g




POWDER OPAQUE 10g / 50g / 200g




POWDER OPAQUE MODIFIER 10g / 50g




MARGIN 10g / 50g



OPACIOUS BODY 10g / 50g 200g



CERVICAL 10g / 50g / 200g



NORITAKE SUPER PORCELAIN EX-3 LINE-UPS

Universal Paste Opaque: Same color as Paste Opaque and Paste Opaque Modifier.


Value Shade: Value Shades are additional shades that expand the existing Noritake Super Porcelain EX-3 selection of body and Opacious Body porcelain powder to include the VITA® 3D-Master® shade.

Speed Enamel (S1, S2, S3, S4): Specially developed for the two layer build-up technique.


n color Shade: n color shade is formulated to achieve chroma intensified Body and Paste Opaque (Universal Paste Opaque, Powder Opaque) of the shade. With intensified chroma, it even can be used in a case where there is insufficient space for porcelain build-up.

COLOR COMBINATION TABLE


BODY 10g / 50g / 200g




ENAMEL 10g / 50g / 200g




SPEED ENAMEL 10g / 50g / 200g




TRANSLUCENT 10g / 50g / 200g




LUSTER 10g / 50g / 200g



CLEAR CERVICAL 10g / 50g / 200g




MAMELON 10g / 50g



TISSUE 10g / 50g



MODIFIER 10g / 50g



ADD-ON 10g / 50g




INTERNAL STAIN 3g



EXTERNAL STAIN 3g (Glaze 10g / 30g)



ADDMATE 10g



Your contact

**Kuraray Europe Benelux,
UK & Ireland**

Margadantstraat 20, 1976 DN IJmuiden, The Netherlands
PO-box 299, 1970 AG IJmuiden, The Netherlands
infoukireland@kuraray.com
Website www.kuraraynoritake.eu



EU Importeur

Kuraray Europe GmbH

Philipp-Reis-Strasse 4,
65795 Hattersheim am Main, Germany
Phone +49 (0)69 305 35 835
Fax +49 (0)69 305 98 35 835
www.kuraraynoritake.eu
centralmarketing@kuraray.com

- Before using this product, be sure to read the Instructions for Use supplied with the product.
- The specifications and appearance of the product are subject to change without notice.
- Printed color can be slightly different from actual color.

"KATANA" and "CERABIEN" are trademarks of NORITAKE CO., LIMITED.



Kuraray Noritake Dental Inc.

300 Higashiyama, Miyoshi-cho, Miyoshi, Aichi 470-0293, Japan
Website www.kuraraynoritake.com

