

# CLEARFIL™ ST OPAQUER



CLEARFIL™ ST OPAQUER is a light-curing opaque resin, especially indicated to mask exposed metal and dark tooth tissue.

## INDICATIONS

- ✓ Masking of the metal during intraoral repair of porcelain fused-to-metal crowns
- ✓ Masking of metal in restoration of post and core crowns
- ✓ Lightening of stained or discolored teeth in porcelain or composite veneer restoration
- ✓ Masking of pulp capping materials

## BENEFITS

- ✓ Highly opaque
- ✓ Immediate coverage of metal in one coat

## PRODUCT ASSORTMENT

ARTICLE NUMBER	PRODUCT	PRODUCT COMPOSITION
#636-EU	CLEARFIL™ ST OPAQUER L	1 syringe: 4.0 g
#637-EU	CLEARFIL™ ST OPAQUER US	1 syringe: 4.0 g

## ALSO RECOMMENDED



CLEARFIL™ CERAMIC PRIMER PLUS is a universal prosthetic primer for all cement and repair indications. It provides excellent bond strength to all silica-based ceramics, zirconia, composites and even metals, all at a high level of procedural simplicity.



ALLOY PRIMER is a metal adhesive primer used to enhance the adhesion of gold, (noble) alloys and other metals to resin composites

For more information visit [kuraraynoritake.eu](http://kuraraynoritake.eu)



EU Importer

**Kuraray Europe GmbH**

Philipp-Reis-Strasse 4,  
65795 Hattersheim am Main, Germany

Phone +49 (0)69 305 35 835

Fax +49 (0)69 305 98 35 835

[www.kuraraynoritake.eu](http://www.kuraraynoritake.eu)

[centralmarketing@kuraray.com](mailto:centralmarketing@kuraray.com)

Before using this product, be sure to read the  
Instructions for Use supplied with the product.



**Kuraray Noritake Dental Inc.**

1621 Sakazu, Kurashiki, Okayama 710-0801, Japan  
Website [www.kuraraynoritake.com](http://www.kuraraynoritake.com)



06/2024