















Restorative composites

Overview

	Class I	Class II	Class III	Class IV	Class V	Light-curing	Self-curing	Radiopaque	Contains color shades	Page number
CLEARFIL MAJESTY™ Esthetic 	 ✓	 ✓	 ✓	 ✓	 ✓	✓		✓	Standard shades: 9, Complementary shades: 8	38
CLEARFIL MAJESTY™ Posterior 	✓	✓	(✓)	(✓)	(✓)	✓		✓	6	40
CLEARFIL MAJESTY™ Flow 	✓	✓	✓		✓	✓		✓	9	42
CLEARFIL™ AP-X 	✓	✓	✓	✓	✓	✓		✓	14	44
CLEARFIL™ PHOTO BRIGHT 			✓	✓	✓	✓			18	46
CLEARFIL™ FIL 			✓	✓	✓		✓		1	47
CLEARFIL™ PHOTO POSTERIOR 	✓	✓				✓		✓	6	48
CLEARFIL™ POSTERIOR 3 	✓	✓					✓	✓	1	49
CLEARFIL™ ST OPAQUER 			✓	✓		✓			2	50

(): possible