





Resin cements

Overview

	Type	Mixing	Indications					Adhesion to						Curing systems			Page number	
			Crowns/ bridges	Inlays/ onlays	Cores/ posts	Adhesion bridges	Veneers	Enamel	Dentin	Metal- oxide ceramics	Silica- based ceramics	Precious metal	Non- precious metal	Curing type	Anaerobic curing	With OXY- GUARD		Fluoride release
CLEARFIL™ SA CEMENT 	self- adhesive resin cement	auto mixing	✓	✓	✓			✓	✓	✓	✓	✓					✓	28
CLEARFIL™ ESTHETIC CEMENT 	self- etching primer / resin cement	auto mixing	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓						30
PANAVIA™ F 2.0 	self- etching primer / resin cement	hand mixing	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	32
PANAVIA™ 21 	self- etching primer / resin cement	hand mixing	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		34

- 1) with ED PRIMER II
- 2) with ED PRIMER
- 3) with CLEARFIL™ CERAMIC PRIMER
- 4) with CLEARFIL™ PORCELAIN BOND ACTIVATOR and CLEARFIL™ SE BOND PRIMER
- 5) with ALLOY PRIMER

* recommendation