


252FP6213-02 06/2019

Dental Zirconia



ENGLISH Caution in conjunction with usage:

When milling the disc, use caution when approaching the milling bur with the plastic ring to prevent detachment of the disc. Mill the disc as if leaving a zirconia part contacting the internal side of the plastic ring at 2 mm.



 **Kuraray Noritake Dental Inc.**
300 Higashiyama, Miyoshi-cho, Miyoshi, Aichi 470-0293, Japan
TEL +81-561-32-8953 FAX +81-561-32-8976



 **Kuraray Europe GmbH**
Philipp-Reis-Str. 4,
65795 Hattersheim am Main, Germany

US Distributed by
Kuraray America, Inc.
32 Old Slip, 7th Floor, New York, NY 10005

210×75mm